

# Aisys

## User's Reference Manual - Part 1 of 2

*Software Revision 3.X*

System Controls, Operation, Checkout, Alarms



## User Responsibility

This Product will perform in conformity with the description thereof contained in this User's Reference manual and accompanying labels and/or inserts, when assembled, operated, maintained, and repaired in accordance with the instructions provided. This Product must be checked periodically. A defective Product should not be used. Parts that are broken, missing, plainly worn, distorted, or contaminated should be replaced immediately. Should repair or replacement become necessary, Datex-Ohmeda recommends that a telephonic or written request for service advice be made to the nearest Datex-Ohmeda Customer Service Center. This Product or any of its parts should not be repaired other than in accordance with written instructions provided by Datex-Ohmeda and by Datex-Ohmeda trained personnel. The Product must not be altered without the prior written approval of Datex-Ohmeda. The user of this Product shall have the sole responsibility for any malfunction which results from improper use, faulty maintenance, improper repair, damage, or alteration by anyone other than Datex-Ohmeda.

**CAUTION** U.S. Federal law restricts this device to sale by or on the order of a licensed medical practitioner. Outside the U.S.A., check local laws for any restriction that may apply.

Datex-Ohmeda products have unit serial numbers with coded logic which indicates a product group code, the year of manufacture, and a sequential unit number for identification. The serial number can be in one of two formats.

AAAX11111	AAAXX111111AA
The <b>X</b> represents an alpha character indicating the year the product was manufactured; <b>H</b> = 2004, <b>J</b> = 2005, etc. <b>I</b> and <b>O</b> are not used.	The <b>XX</b> represents a number indicating the year the product was manufactured; <b>04</b> = 2004, <b>05</b> = 2005, etc.

**Aisys, Carestation, Advanced Breathing System, Aladin, Aladin<sub>2</sub>, ComWheel, D-fend, Disposable Multi Absorber, Easy-Fil, Reusable Multi Absorber, PSVPro, and SmartVent** are registered trademarks of Datex-Ohmeda, Inc.

Other brand names or product names used in this manual are trademarks or registered trademarks of their respective holders.



# Table of Contents

## 1 Introduction

What is an Aisys Carestation? .....	1-2
Symbols used in the manual or on the equipment .....	1-4
Typeface conventions used .....	1-7
Abbreviations .....	1-7

## 2 System Controls and Menus

System overview .....	2-2
Aladin cassette controls .....	2-5
Advanced breathing system (ABS) components .....	2-6
Optional ABS components .....	2-8
Non-circle circuit .....	2-9
ACGO (optional) .....	2-9
Scavenging the ACGO sample flow .....	2-10
Scavenging from an auxiliary manual breathing circuit .....	2-10
Scavenging a gas monitor sample flow .....	2-10
Display controls .....	2-11
Anesthesia system display .....	2-12
Using menus .....	2-14

## 3 Operation

Turning on the system .....	3-2
Start gas flow (start case) .....	3-3
Using quick start .....	3-3
Using default settings .....	3-3
Presetting ventilator and alarm settings .....	3-3
End case .....	3-4
Ventilator setup .....	3-5
Using quick keys .....	3-5
Using Vent Setup menu .....	3-5

Changing ventilator modes and settings .....	3-6
Manual ventilation cardiac bypass .....	3-6
VCV cardiac bypass (optional) .....	3-6
Gas setup .....	3-7
Using quick keys .....	3-7
Using Gas Setup menu .....	3-7
Changing gas and settings .....	3-7
Changing circuit type .....	3-7
Spirometry setup .....	3-8
Setting patient and sensor type .....	3-8
Selecting a data source .....	3-9
Setting spirometry split screen .....	3-9
Setting loop type .....	3-9
Viewing spirometry loops .....	3-10
Screen configuration .....	3-11
Split screen .....	3-11
Setting gas controls and screen configuration .....	3-11
Alarm setup .....	3-12
Setting alarm limits .....	3-12
MV/TV Alarms .....	3-12
CO2 Audio .....	3-12
Auto MV Limit .....	3-12
Viewing trends .....	3-13
Alternate O <sub>2</sub> control .....	3-14
EZchange Canister (optional) .....	3-15
Condenser (optional) .....	3-16
Passive AGSS (optional) .....	3-17
Active AGSS (optional) .....	3-18
Connecting active AGSS with a flow indicator .....	3-19
Connecting active adjustable AGSS .....	3-19

## 4 Preoperative Checkout

Every day before your first patient .....	4-2
Before every patient .....	4-3

## 5 Preoperative Tests

Inspect the system .....	5-2
Leak < 250 ml .....	5-3
No .....	5-3
Yes .....	5-3
Machine Check .....	5-3
Machine Check - System .....	5-4
Machine Check - Circuit .....	5-4
Machine Check - Circuit O2 .....	5-4
Individual Checks .....	5-5
System .....	5-5
Circuit .....	5-5
Circuit O2 Cell .....	5-6
Low P Leak .....	5-6
Low P Leak (machines with ACGO) .....	5-6
Agent Delivery .....	5-7
Positive low-pressure leak test (ACGO systems only) .....	5-8
Aladin cassette installation .....	5-9
Flow sensor calibration .....	5-9

## 6 Airway Modules

Airway modules .....	6-2
Connection to a patient .....	6-4
Parameters setup .....	6-5
Data source .....	6-5
CO <sub>2</sub> setup .....	6-5
O <sub>2</sub> setup .....	6-5
Agent setup .....	6-5
Spirometry setup .....	6-5
Automatic agent identification .....	6-6
Calibration .....	6-7

## 7 Vaporizer

Vaporizer .....	7-2
Aladin <sub>2</sub> cassette variants .....	7-2
Aladin cassette variants .....	7-3
Installing cassettes .....	7-5
Cleaning .....	7-5
Service .....	7-5
Draining halothane cassettes .....	7-5
Filling Aladin <sub>2</sub> cassettes .....	7-6
Filling with Easy-Fil system .....	7-6
Filling with Quik-Fil system .....	7-8
Filling with Saf-T-Fil bottle .....	7-9
Filling Aladin cassettes .....	7-11
Filling with keyed filler system .....	7-11
Filling with Quik-Fil system .....	7-13
Filling with Saf-T-Fil bottle .....	7-14

## 8 Alarms and Troubleshooting

Alarms .....	8-2
Alarm priorities .....	8-2
Display changes during alarms .....	8-2
Battery indicator .....	8-3
Internal failure .....	8-3
List of alarms .....	8-3
Alarm ranges and default values .....	8-10
Alarm tests .....	8-11
Breathing system problems .....	8-12
Electrical problems .....	8-13
Pneumatic problems .....	8-13

## Warranty

## Index

# 1 Introduction

<b>In this section</b>	What is an Aisys Carestation? . . . . .	1-2
	Symbols used in the manual or on the equipment . . . . .	1-4
	Typeface conventions used . . . . .	1-7
	Abbreviations . . . . .	1-7



## What is an Aisys Carestation?

Aisys means **A**nesthesia **I**ntegrated **S**ystem.


The Aisys Carestation is scalable, flexible, and functionally integrated, featuring the most advanced design, ventilation, respiratory monitoring, and breathing system. In addition, module bays allow for the physical integration of legacy Datex-Ohmeda patient monitors and supports mounting of other GE Healthcare monitors. Optionally, the open architecture design supports mounting of non-Datex Ohmeda patient monitors, record keeping, and connections to the hospital information system. Aisys Carestation manages the complexity of anesthesia while supporting a broad range of patients from neonates and compromised patients to routine cases. The INview movable display arm helps keep the anesthetist's focus on the patient by offering control of all hemodynamic, gas delivery, anesthetic agent and ventilation parameters.

Aisys Carestation uses SmartVent ventilation technology offering Volume Control Ventilation with tidal volume compensation, Pressure Control Ventilation, and electronic PEEP. The proven SmartVent also features optional Pressure Support Ventilation with an Apnea Backup (PSVPro) that is used for spontaneously breathing patients, Synchronized Intermittent Mandatory Ventilation (SIMV), Pressure control ventilation-volume guarantee (PCV-VG), and VCV cardiac bypass. The Advanced Breathing System (ABS) has a small circuit volume, optimized for low and minimal flow applications. The ABS, designed to provide a larger work surface, is elegantly integrated and easy to remove and clean.

The Aisys Carestation is designed for expansion and upgrades, so it is easy to add new technologies and ventilation capabilities without investing in a new system.

The Aisys Carestation is not suitable for use in a MRI environment.

**Note** Configurations available for this product depend on local market and standards requirements. Illustrations in this manual may not represent all configurations of the product. This manual does not cover the operation of every accessory. Refer to the accessory documentation for further information.



AB.75.100

Figure 1-1 ▪ Aisys Carestation























## Symbols used in the manual or on the equipment

Symbols replace words on the equipment, on the display, or in Datex-Ohmeda manuals.

Warnings and Cautions tell you about dangerous conditions that can occur if you do not follow all instructions in this manual.

Warnings tell about a condition that can cause injury to the operator or the patient.

Cautions tell about a condition that can cause damage to the equipment. Read and follow all warnings and cautions.

	On (power)		Off (power)
	Standby		O <sub>2</sub> Flush button
	Type BF equipment		Type B equipment
	Dangerous voltage		Alternating current
	Caution, ISO 7000-0434		Attention, refer to product instructions, IEC 60601-1
	Electrical input		Electrical output
	Pneumatic inlet		Pneumatic outlet
	Serial number		Stock number
	Equipotential		Lamp, lighting, illumination
	Variability		Variability in steps
	Suction bottle outlet		Vacuum inlet

**MAX**

Maximum

**+**

Plus, positive polarity



Bag position/manual ventilation



Inspiratory flow



Movement in one direction



Lock



Isolation transformer

**134°C**

Autoclavable



This way up



No battery/battery failure

⚠ ≈ cmH<sub>2</sub>O


APL settings are approximate



Alarm silence

**Rx ONLY**

Caution: federal law prohibits dispensing without prescription.



Pipeline



Read to center of float

**VACUUM**

Vacuum

**-**

Minus, negative polarity



Mechanical ventilation



Expiratory flow



Movement in two directions



Unlock



Low pressure leak test



Not autoclavable

**O<sub>2</sub>%**O<sub>2</sub> cell connection

Battery in use. Bar indicates amount of battery power remaining.

**AGSS**

Anesthetic Gas Scavenging System



Submenu



Pinch hazard



Cylinder

EZchange Canister (CO<sub>2</sub> bypass)



Systems with this mark agree with the European Council Directive (93/42/EEC) for Medical Devices when they are used as specified in their User's Reference manuals. The xxxx is the certification number of the Notified Body used by Datex-Ohmeda's Quality Systems.



Authorized representative in the European Community



Manufacturer

## Typeface conventions used

Names of hard keys on the display and modules are written in bold typeface; for example, **Normal Screen**.

Menu items are written in bold italic typeface; for example, ***Vent Setup***.

Messages that are displayed on the screen are enclosed in single quotes; for example, 'Check sample gas out.'

When referring to different sections and other documents, the names are written in italic typeface and enclosed in double quotes; for example, "*System Controls and Menus*."

## Abbreviations

Abbreviation	Definition
<b>C</b>	
Compl	Compliance
<b>E</b>	
ET	End-tidal concentration
EtCO <sub>2</sub>	End-tidal carbon dioxide
EtO <sub>2</sub>	End-tidal oxygen
<b>F</b>	
F-V	Flow-Volume loop
FI	Fraction of inspired gas
FI-ET	Difference between inspiratory and expiratory concentrations
FiO <sub>2</sub>	Fraction of inspired oxygen
<b>I</b>	
I:E	Inspiratory-expiratory ratio
Insp Pause	Inspiratory pause time
<b>M</b>	
MV	Minute volume
MVexp	Expired minute volume
MVinsp	Inspired minute volume
<b>O</b>	
O <sub>2</sub>	Oxygen

Abbreviation	Definition
<b>P</b>	
P-F	Pressure-Flow loop
P-V	Pressure-Volume loop
Pair	Air supply pressure
Paux	Auxiliary pressure
Paw	Airway pressure
PCV	Pressure controlled ventilation
PCV-VG	Pressure controlled ventilation-volume guaranteed
PEEP	Positive end expiratory pressure
PEEPe	Extrinsic positive end expiratory pressure
Pexp	Expiratory pressure
Pinsp	Inspiratory pressure
Plimit	High pressure limit
Pmax	Maximum pressure
Pmean	Mean pressure
PO <sub>2</sub>	Oxygen supply pressure
Ppeak	Peak pressure
Pplat	Plateau pressure
Psupp	Support pressure
PSVPro	Pressure supported ventilation with apnea backup
<b>R</b>	
Rate	Respiratory rate
Raw	Airway resistance
RR	Respiratory rate
<b>S</b>	
SIMV-PC	Synchronized intermittent mandatory ventilation-pressure controlled
SIMV-VC	Synchronized intermittent mandatory ventilation-volume controlled
SIMV/PSV	Synchronized intermittent mandatory ventilation/pressure supported ventilation
<b>T</b>	
TV	Tidal volume
TVexp	Expired tidal volume

Abbreviation	Definition
TV <sub>insp</sub>	Inspired tidal volume
<b>V</b>	
VCO <sub>2</sub>	Carbon dioxide production
VCV	Volume controlled ventilation
Vol	Volume






# 2 System Controls and Menus

<b>In this section</b>	System overview .....	2-2
	Advanced breathing system (ABS) components .....	2-6
	Display controls .....	2-11
	Anesthesia system display .....	2-12
	Using menus .....	2-14

## System overview

**WARNING** Explosion hazard. Do not use this system with flammable anesthetic agents.










- |   |  |
|---|--|
| 1. Airway module (optional)                           | 9. O <sub>2</sub> flush button*                        |
| 2. Datex-Ohmeda patient monitoring modules (optional) | 10. System switch*                                     |
| 3. Dovetail/GCX rails                                 | 11. Integrated suction (optional)*                     |
| 4. Light switch                                       | 12. Advanced breathing system                          |
| 5. Alternate O <sub>2</sub> control*                  | 13. Auxiliary O <sub>2</sub> flow control (optional)*  |
| 6. Mains indicator*                                   | 14. Aladin cassette and active bay                     |
| 7. Brake*   | 15. Anesthesia display                                 |
| 8. Aladin cassette storage bay                        | 16. Datex-Ohmeda patient monitoring display (optional) |

Figure 2-1 ■ Front view

\*These items are described in further detail in the following table.

AB 75.097

Item, Figure 2-1	Description
5 Alternate O <sub>2</sub> control	<p>Alternate O<sub>2</sub> control activates automatically in the case of certain failures or errors. It can also be activated manually by pushing the Alternate O<sub>2</sub> control button if the display fails. Turn the knob counterclockwise to increase the flow. Turn the knob clockwise to decrease the flow. See <i>“Alternate O<sub>2</sub> Control”</i> in the Operation section of this manual.</p> 
6 Mains indicator	The mains indicator will come on when AC power is connected.
7 Brake	<p>Push down to lock. Lift to release.</p> 
9 O <sub>2</sub> flush button	<p>Push the O<sub>2</sub> flush button to supply high flows of O<sub>2</sub> to the breathing system.</p> 
10 System switch	<p>Set the switch to the On (I) position to permit gas flow and to turn on the system.</p> 
11 Integrated suction	<p>Turn the switch to Max for full vacuum. Turn the switch to Off for no vacuum. Turn the switch to On for adjustable vacuum. When in adjustable vacuum, turn the knob clockwise to increase the vacuum and counterclockwise to decrease the vacuum.</p> 
13 Auxiliary O <sub>2</sub> flow control	<p>Turn the knob counterclockwise to increase the flow. Turn the knob clockwise to decrease the flow.</p> 




AB: 75.046

- |  |                                |
|--|--------------------------------|
| 1. Serial port                             | 7. Mains inlet                 |
| 2. Collection bottle connection            | 8. System circuit breaker      |
| 3. Cylinder wrench (key) storage           | 9. Equipotential stud          |
| 4. Cylinder yoke                           | 10. Outlet circuit breaker     |
| 5. AGSS (Anesthesia Gas Scavenging System) | 11. Isolated electrical outlet |
| 6. Pipeline connections                    |                                |

Figure 2-2 ■ Rear view

## Aladin cassette controls

The electronically controlled vaporizer consists of the internal electronic control unit and the Aladin agent cassette. See “Vaporizer” section for more information.




AB 60.038

1. Handle with release trigger
2. Lock
3. Liquid level indicator
4. Agent filling port

Figure 2-3 ▪ Aladin<sub>2</sub> cassette

## Advanced breathing system (ABS) components





AB.75.098

- |                               |   |
|-------------------------------|---|
| 1. Airway module (optional)   | 8. Leak test plug                             |
| 2. Expiratory check valve     | 9. Manual bag port                            |
| 3. Inspiratory check valve    | 10. Breathing system release                  |
| 4. Inspiratory flow sensor*   | 11. Adjustable pressure-limiting (APL) valve* |
| 5. Expiratory flow sensor*    | 12. Bag/Vent switch*                          |
| 6. Absorber canister          | 13. Bellows assembly                          |
| 7. Absorber canister release* |   |


Figure 2-4 ■ Advanced breathing system

\*These items are described in further detail in the following table.

Item, Figure 2-4	Description
4, 5      Flow sensor	Flow sensors provide volume measurements for some monitoring functions and tidal volume delivery.
7      Absorber canister release	<p>Push to remove the canister. This causes the breathing system to vent to the room (unless the EZchange Canister option is installed). Be sure to hold the canister by the handle before releasing the canister.</p> 
11      APL valve	<p>Adjusts breathing system pressure limit during manual ventilation. The scale shows approximate pressures. Above 30 cmH<sub>2</sub>O, the knob will click as it turns.</p> 
12      Bag/Vent switch	Selects between manual ventilation (bag) or mechanical ventilation (ventilator).



Optional ABS  
components




AB.75.049  
AB.82.043

- 1. Bag support arm\*
- 2. Auxiliary Common Gas Outlet (ACGO) switch\*
- 3. ACGO port
- 4. EZchange Canister system (CO<sub>2</sub> bypass)
- 5. EZchange Canister release\*
- 6. Condenser drain button\*
- 7. Condenser

Figure 2-5 ■ Breathing system options

\*These items are described in further detail in the following table.

Item, Figure 2-5		Description
1	Bag support arm	Squeeze the button to raise or lower the arm. 
2	ACGO switch	Set switch to ACGO position for fresh gas to flow through the ACGO port. The ACGO may be used to provide fresh gas to an auxiliary manual breathing circuit.
5	EZchange Canister release	Push to drop the canister to EZchange position. This seals the breathing circuit, permitting continued ventilation and rebreathing of exhaled gases. Be sure to hold the canister by the handle before releasing the canister.
6	Condenser drain button	Push to drain water out of the condenser.

**Non-circle circuit**

Fresh gas flow is diverted around the inspiratory check valve and out through the inspiratory port when non-circle ventilation is selected. This fresh gas source may be used with circuits without CO<sub>2</sub> absorbent capability (e.g., Mapleson variants). Mechanical ventilation is not available when using the non-circle circuit. Tidal volume monitoring is not available.

O<sub>2</sub> monitoring of fresh gas is available automatically when the non-circle circuit is selected if the system has the airway module option or the O<sub>2</sub> cell monitoring option.

Fresh gas oxygen concentration is displayed on the screen. Set the alarm limits appropriately. Note that fresh gas oxygen concentration may not reflect FiO<sub>2</sub> when using these types of circuits. Use an external O<sub>2</sub> monitor when using a rebreathing circuit with the non-circle circuit.

Systems with both an airway module and an O<sub>2</sub> cell will display the circuit O<sub>2</sub> value obtained from the airway module.

Do not use an external ventilator when using the non-circle circuit. Do not use the non-circle circuit to drive external ventilators or for jet ventilation.

**WARNING**

The maximum pressure at the non-circle circuit can be up to 27 kPa (4 psi). Use a breathing circuit with pressure relief.

**ACGO (optional)**

Fresh gas flow is directed through the Auxiliary Common Gas Outlet (ACGO) on the front of the machine when the ACGO switch is in the ACGO position. Mechanical ventilation is not available when operating an auxiliary manual breathing circuit with fresh gas from the ACGO. The Bag/Vent switch, APL valve, and bag arm are not part of the external circuit. Volume and pressure monitoring are not available.

O<sub>2</sub> monitoring of fresh gas is available automatically when the ACGO is selected if the system has the airway module option or the O<sub>2</sub> cell monitoring option. A sample of the fresh gas is diverted to the O<sub>2</sub> cell in the breathing system. The sample flow to the O<sub>2</sub> cell is dependent on the pressure in the external circuit. The sample flow reduces the fresh gas flow rate to the auxiliary breathing circuit equal to the amount diverted to the O<sub>2</sub> cell.

Fresh gas oxygen concentration is displayed on the screen. Set the alarm limits appropriately. Note that fresh gas oxygen concentration may not reflect FiO<sub>2</sub> during spontaneous breathing or in rebreathing circuits. Use an external O<sub>2</sub> monitor if using a rebreathing circuit on ACGO.

Systems with both an airway module and an O<sub>2</sub> cell will display the circuit O<sub>2</sub> value obtained from the airway module.

Do not use an external ventilator on the ACGO. Do not use the ACGO to drive external ventilators or for jet ventilation.

**WARNING**

The maximum pressure at the ACGO can be up to 55 kPa (8 psi). Use a breathing circuit with pressure relief.

## **Scavenging the ACGO sample flow**

A sample of the fresh gas is diverted to the airway module or the O<sub>2</sub> cell in the breathing system. If an auxiliary manual breathing circuit is to be used with N<sub>2</sub>O or volatile anesthetics, this sample flow should be scavenged.

1. Attach a circle breathing circuit to the inspiratory and expiratory ports. Occlude the circle circuit by connecting the Y-piece to the plug located to the rear of the expiratory port.
2. Check for clinically correct settings. Set the Bag/Vent switch to mechanical ventilation mode. (Mechanical ventilation will not start when ACGO is selected.) Alternatively, set the Bag/Vent switch to the bag mode, set the APL valve to MIN, and attach a bag.
3. The bellows, or bag, will fill slowly with the fresh gas sample flow and then spill to the AGSS.

## **Scavenging from an auxiliary manual breathing circuit**

Scavenge the exhaust if an auxiliary manual breathing circuit is used with N<sub>2</sub>O or volatile anesthetics.

An auxiliary inlet is available for active and passive AGSS units. It provides a female connection with 30 mm - 30 mm male connector (or a 30 mm - 19 mm male connector) into the auxiliary port under the breathing system. Do not use these connectors as an outlet for exhaust flow.


The auxiliary inlet is a convenience inlet to the air brake of active AGSS units. There is a reservoir to capture exhaust flows higher than the extract flow.

A separate exhaust hose is needed from the auxiliary manual breathing circuit to the disposal point for all AGSS units.

## **Scavenging a gas monitor sample flow**

Sample gas from a gas monitor can be scavenged using the sample gas return port or the AGSS. To scavenge from a gas monitor using the sample gas return port, connect the tubing from the monitor to the sample gas return port. To scavenge from a gas monitor using the AGSS, connect tubing from the monitor to the male luer inlet on the bottom of the AGSS underneath the breathing system.

## Display controls




AB.75.099

- |                       |  |
|-----------------------|--|
| 1. Silence Alarms key | Push to silence any active, silenceable high and medium priority alarms or to suspend/acknowledge any non-active medium or high priority alarms. Alarm is silenced for 120 seconds or alarm is suspended for 90 seconds. |
| 2. Menu keys          | Push to show corresponding menu.   |
| 3. ComWheel           | Push to select a menu item or confirm a setting. Turn clockwise or counterclockwise to scroll menu items or change settings.   |
| 4. Normal Screen key  | Push to remove all menus from the screen.  |
| 5. Quick keys         | Push to change corresponding gas setting or ventilator setting. Turn the ComWheel to make a change. Push the ComWheel to activate the change.  |
| 6. Timer keys         | Push to start or stop the timer. Push to reset the timer back to zero.   |
| 7. MV/TV Alarms key   | Push to turn off the MV and TV alarms. Push again to turn the MV and TV alarms back on.  |

Figure 2-6 ▪ Display controls


# Anesthesia system display



1. Electronic gas flow indicators
2. Alarm silence countdown
3. Alarm message fields
4. Waveform fields
5. General message field or timer field
6. Clock
7. Battery indicator field
8. Measured values field
9. Pipeline and cylinder supply or respiratory data
10. Ventilator settings
11. Ventilation mode
12. Gas and agent settings

Figure 2-7 • Normal view

When a menu key is selected, the menu field overlays the gas flow indicators and the waveform fields start at the right edge of the menu.




AB.75.007

1. Menu  
2. Waveform fields  
Figure 2-8 ▪ Menu view

## Using menus

Push a menu key to display the corresponding menu. Use the ComWheel to navigate through the menu.




AB-91.007

1. Menu title
2. Present selection
3. Adjustment window
4. Indicates submenu
5. Short instructions
6. Menu selections


Figure 2-9 ■ Example menu

1. Push the menu key to display the corresponding menu.



AB.91.004


2. Turn the ComWheel counterclockwise to highlight the next menu item.  
Turn the ComWheel clockwise to highlight the previous menu item.



AB.91.005




3. Push the ComWheel to enter the adjustment window or a submenu.



AB.91.006

4. Turn the ComWheel clockwise or counterclockwise to highlight the desired selection.




AB.91.007


5. Push the ComWheel to confirm the selection.

Xxxx Xxxx

Xxxxxx




Xxxxxx Xxxxxx




Xxxxxxx Xxxxxxx

Xx


Xxxxxx Xxxxxx



Xxxxxxxx Xxxxx



Xxxxxxxxxxxx



Xxxxxx Xxxxxx


Xxxx xxx Xxxxxx Xxxxxx.

AB.9.1.008


6. Select **Normal Screen** or push the **Normal Screen** key to exit the menu and return to the normal monitoring display. (Select **Previous Menu** to return to the last displayed menu, if available.)

Xxxx Xxxx

Xxxxxx




Xxxxxx Xxxxxx




Xxxxxxx Xxxxxxx

Xx


Xxxxxx Xxxxxx



Xxxxxxxx Xxxxx



Xxxxxxxxxxxx



Xxxxxx Xxxxxx

Xxxx xxx Xxxxxx Xxxxxx.


AB.9.1.009



## 3 Operation

**WARNING** Ventilator alarms indicate potential hazard conditions. All alarms that occur should be investigated to help ensure adequate patient safety.


 Maintain sufficient fresh gas flow when using sevoflurane.

 Desiccated (dehydrated) absorbent material may produce dangerous chemical reactions when exposed to inhalation anesthetics. Adequate precautions should be taken to ensure that absorbent does not dry out. Turn off all gases when finished using the system.

<b>In this section</b>	Turning on the system . . . . .	3-2
	Start gas flow (start case) . . . . .	3-3
	End case . . . . .	3-4
	Ventilator setup . . . . .	3-5
	Gas setup . . . . .	3-7
	Spirometry setup . . . . .	3-8
	Viewing spirometry loops . . . . .	3-10
	Screen configuration . . . . .	3-11
	Alarm setup . . . . .	3-12
	Viewing trends . . . . .	3-13
	Alternate O <sub>2</sub> control . . . . .	3-14
	EZchange Canister (optional) . . . . .	3-15
	Condenser (optional) . . . . .	3-16
	Passive AGSS (optional) . . . . .	3-17
	Active AGSS (optional) . . . . .	3-18

## Turning on the system

1. Plug the power cord into a electrical outlet. Make sure the system circuit breaker is on.
  - The mains indicator will come on when AC power is connected. Battery is charging if not already fully charged.



1. System switch
2. Mains indicator

Figure 3-1 ▪ Mains indicator and system switch

2. Check that the breathing system is properly connected.

**CAUTION** Do not turn on the system with the right-hand (inspiratory) port plugged.

3. Turn the System switch to On.
  - The display shows the power-up screen.
  - The system does a series of automated self tests.
4. Perform a **Machine Check** before the first case of the day. Perform the **Individual Checks** as needed before each case. See sections “Preoperative Checkout” and “Preoperative Tests.”

## Start gas flow (start case)

### Using quick start

To quickly start a case using the defaults shown on the **Start Case** menu, press any quick key once.

### Using default settings

1. Set the Bag/Vent switch to Bag.
2. Push the **Start/End Case** key to go to the **Start Case** menu or select **Start Case** from the **Checkout** menu.
3. Select the ventilator and alarm settings from **Defaults**, **Last Case**, or **Patient Weight**.
4. Select **Other Gas**, and select a balance gas to use with O<sub>2</sub>.
5. Select **Circuit**, and select **Circle** circuit or **Non-Circle** circuit.
6. To start a case using the ventilator and alarm settings selected, select **Start Case Now**. Gas flow starts.

### Presetting ventilator and alarm settings

1. Set the Bag/Vent switch to Bag.
2. Push the **Start/End Case** key to go to the **Start Case** menu or select **Start Case** from the **Checkout** menu.
3. Select the ventilator and alarm settings from **Defaults**, **Last Case**, or **Patient Weight**.
4. Select **Other Gas**, and select a balance gas to use with O<sub>2</sub>.
5. Select **Circuit**, and select **Circle** circuit or **Non-Circle** circuit.
6. To change the alarm settings before starting a case:
  - Push the **Alarm Setup** key.
  - To change the alarm settings, select **Adjust Settings**.
    - Change the settings.
    - Select **Back**.
  - Make any other changes to the alarm setup.
  - To start a case, select **Start Case Now**. Gas flow starts.
  - To change the ventilator settings, select the **Vent Setup** key. Follow the instructions in the next steps to change the ventilator settings.

#### Note

If changes are only made to the alarm settings and no changes are made to the ventilator settings, the ventilator settings are set to the selection made on the **Start Case** menu: Defaults, Last Case, or Patient Weight.

7. To change ventilator settings for the selected mode before starting a case:
  - Push the **Vent Setup** key.
  - Select **Adjust Settings**.
  - Change the settings.
  - Select **Exit**.
  - Select **Start Case Now**. Gas flow starts.
8. To change the ventilation mode before starting a case:
  - Push the **Vent Setup** key.
  - Select the desired mode.
  - Select **Confirm** or adjust the settings and select **Confirm**.
  - Select **Start Case Now**. Gas flow starts.

**Note** If changes are made only to the ventilator settings and no changes are made to the alarm settings, the alarm settings are set to the selection made in the **Start Case** menu: Defaults, Last Case, or Patient Weight.

## End case

1. Set the Bag/Vent switch to Bag.
2. Push the **Start/End Case** key.
3. Select **End Case Now** to put the system in standby. (Stops gas flow and ventilator alarms. Stops alarms on battery-backup patient monitors.)

## Ventilator setup

The system has six modes of mechanical ventilation:

- volume control (VCV).
- pressure control (PCV) (optional).
- synchronized intermittent mandatory/pressure support (SIMV/PSV) (optional).
- pressure support (PSVPro) (optional).
- synchronized intermittent mandatory-pressure control (SIMV-PC) (optional).
- pressure control ventilation-volume guarantee (PCV-VG) (optional).

### WARNING

Most anesthetic agents will cause patients to have reduced ventilatory responses to carbon dioxide and to hypoxemia. Therefore, triggered modes of ventilation may not produce adequate ventilation.



The use of neuromuscular blocking agents will reduce the patient's breathing response, which will interfere with triggering.

### Note

Refer to the “*Specifications and Theory of Operation*” section in Part 2 of the User's Reference manual for more information on ventilation modes.

## Using quick keys

The three main ventilator settings for each mode can easily be changed using the ventilator quick keys.

1. Push a ventilator quick key to select the corresponding ventilator setting.
2. Turn the ComWheel to make a change.
3. Push the ComWheel to activate (confirm) the change.

## Using Vent Setup menu

1. Push the **Vent Setup** key.
  - An arrow to the right of the mode indicates the current mode.
  - **Adjust Settings** is selected by default.
2. Push the ComWheel to enter the adjustment window for the selected mode.
3. Use the ComWheel to navigate through the adjustment window and to change a value.
4. Push the ComWheel to activate the change.
5. Push the **Normal Screen** key or select **Exit**.



## Changing ventilator modes and settings

1. Push the **Vent Setup** key.
2. Select the mode and push the ComWheel to enter the adjustment window.
3. Use the ComWheel to navigate through the adjustment window and to change a value.
4. Confirm the settings to activate the mode.

## Manual ventilation cardiac bypass

Manual ventilation cardiac bypass suspends alarms for patients on cardiac bypass when the ventilator is not mechanically ventilating. The volume, apnea, PEEP, low agent, low CO<sub>2</sub>, and respiratory rate alarms are suspended. The alarms are enabled when cardiac bypass is turned off or mechanical ventilation is started.

1. Set the Bag/Vent switch to Bag.
2. Push the **Main Menu** key.
3. Set **Cardiac Bypass** to **On**.
  - The general message 'Cardiac Bypass' shows in the waveforms and in the general message field when passive cardiac bypass is active.

## VCV cardiac bypass (optional)

Systems equipped with VCV cardiac bypass and enabled in the **Super User** menu can use the cardiac bypass feature during mechanical ventilation while in VCV mode. Only VCV mode is available while using VCV cardiac bypass. The volume, apnea, PEEP, low agent, low CO<sub>2</sub>, and respiratory rate alarms are suspended. The alarms are enabled when VCV cardiac bypass is turned off or mechanical ventilation is stopped.

1. Start mechanical ventilation in VCV mode.
2. Push the **Main Menu** key.
3. Set **Cardiac Bypass** to **On**.
  - PEEP is set to 5 cmH<sub>2</sub>O and TV is set to 170 ml when the VCV setting for TV is greater than 170 ml.
  - The general message 'VCV Cardiac Bypass' shows in the waveforms and in the general message field when VCV cardiac bypass is active.

## Gas setup

### Using quick keys

The O<sub>2</sub>%, total flow, and agent concentration can easily be changed using the gas quick keys.

1. Push a gas quick key to select the corresponding gas setting.
2. Turn the ComWheel to make a change.
3. Push the ComWheel to activate (confirm) the change.

### Using Gas Setup menu

1. Push the **Gas Setup** key.
  - An arrow to the right of the balance gas indicates the gas currently being used with O<sub>2</sub>.
  - **Adjust Settings** is selected by default.
2. Push the ComWheel to enter the adjustment window.
3. Use the ComWheel to navigate through the adjustment window and change a value.
4. Push the ComWheel to activate the change.
5. Push the **Normal Screen** key or select **Exit**.

### Changing gas and settings

1. Push the **Gas Setup** key.
2. Select the balance gas to use with O<sub>2</sub> and push the ComWheel to enter the adjustment window.
3. Use the ComWheel to navigate through the adjustment window and change a value.
4. Push the ComWheel to activate the change.

### Changing circuit type

1. Push the **Gas Setup** key.
  - An arrow to the right of the circuit type indicates the circuit currently in use.
2. Use the ComWheel to select **Circle** circuit or **Non-Circle** circuit and enter the adjustment window.
3. Confirm or change and confirm the **O<sub>2</sub>** and **Total Flow** settings.
4. Push the **Normal Screen** key or select **Exit**.

## Spirometry setup

Push the **Spirometry** key and select **Spiro Setup** to access the **Spirometry Setup** menu. Use the **Spirometry Setup** menu to:

- Set the patient and sensor type.
- Select a data source.
- View, save, and erase spirometry loops; and to set the loop scaling.
- Change the volume shown on the spirometry split screen to MVexp or TVexp.
- Set split-screen view.
- For quick access to the Paw and MVexp alarm limits.

**Note** The **Spirometry Setup** menu is also accessible through **Main Menu - Parameter Setup - Spirometry Setup**.

### Setting patient and sensor type

Patient and sensor type refer to the style of airway adapter used with the airway module. If spirometry data is obtained from the airway module, ensure that the sensor type matches the airway adapter used. If an airway module is not installed, **Patient and sensor type** is not selectable.

If the sensor type is not set correctly, the information displayed may not be accurate.

1. Push the **Spirometry** key.
2. Select **Spiro Setup - Patient and Sensor Type**.
3. Select **Adult** or **Pedi** depending on the sensor used.
  - Adult refers to the D-lite sensor.
  - Pedi refers to the Pedi-lite sensor.

## Selecting a data source

Several monitoring parameters can be obtained from the ventilator or the airway module. Information that is retrieved from the airway module is identified with the module data indicator. Refer to the “*Airway Modules*” for more information.




Figure 3-2 ▪ Airway module data indicator

1. Push the **Spirometry** key.
2. Select **Spiro Setup - Data Source**.
3. Select **Patient** or **Vent** as the primary source for information.
  - If **Patient** is selected, the airway module will be the first source for information if a D-lite or Pedi-lite sensor is properly connected. If the sensor is not properly connected, but the airway module is installed the waveform data is not displayed. If information is not available through the airway module, information will come from the internal ventilator sensors.
  - If **Vent** is selected, the internal sensors of the ventilator will be the first source for information.

## Setting spirometry split screen

Spirometry loops can be viewed alongside the waveforms on the normal screen. To set up the spirometry split screen:

1. Push the **Spirometry** key.
2. Select **Spiro Setup**.
3. Select **Split Screen - Spiro**.
4. Push the **Normal Screen** key.

## Setting loop type

Spirometry loop type shown in the spirometry split screen and in the spirometry split screen view can be set in the **Spirometry Setup** menu.

1. Push the **Spirometry** key.
2. Select **Spiro Setup - Loop Type** and choose the desired view.


**Note** The loop type can also be set in the **Spirometry Menu**. The **Spirometry Setup** menu can also be accessed through **Main Menu - Parameter Setup - Spirometry Setup**.

## Viewing spirometry loops

There are three types of spirometry loops: Pressure-Volume (P-V), Flow-Volume (F-V) and Pressure-Flow (P-F). View, save, and erase spirometry loops and set the loop scaling in the **Spirometry** menu. Push the **Spirometry** key to go to the **Spirometry** menu.

- To set the loop type, select **Loop Type** and the desired view.
- To store a reference loop to memory, select **Save Loop**.
- To view a saved loop, select **Reference Loop** and the time at which the loop was saved.
- To erase a saved loop, select **Erase Loop** and the time at which the loop was saved.
- To adjust the loop scaling, select **Scaling**.

**Note** The loop type can also be set in the **Spirometry Setup** menu. Go to **Main Menu - Parameter Setup - Spirometry Setup**. Select **Loop Type** and choose the desired view.



1. Volume axis
2. Pressure axis
3. Real-time loop
4. Reference loop (appears on display in white)

Figure 3-3 ▪ Example of a P-V loop

## Screen configuration

Waveform, fresh gas control style, split screen, sweep speed, time and date, and brightness are adjusted in the **Screen Setup** menu. The top waveform is fixed at Paw. The center waveform can be adjusted to show agent or flow, or it can be turned off. The bottom waveform can be set to show CO<sub>2</sub>, or it can be turned off. To change the display:

1. Push the **Main Menu** key.
2. Select **Screen Setup**.
3. Select the field to change, and make the change.
4. Push the **Normal Screen** key or select **Previous Menu**.

### Split screen

Gas and agent delivery, trends, spirometry loops, a metabolics chart, or a paw gauge can be shown as a split screen alongside the waveforms on the normal screen. Split screen settings cannot be saved as defaults. To change the normal screen to show a split screen:

1. Push the **Main Menu** key.
2. Select **Screen Setup - Split Screen**.
3. Select the desired view.
4. Push the **Normal Screen** key or select **Previous Menu**.

### Setting gas controls and screen configuration

There are two ways to configure gas controls: O2%/total flow and individual gas flow. The Super User selects the configuration or can allow the user to select the configuration.

Selecting **O2** shows O2% as the first quick key and total flow l/min as the second quick key. The gas control/screen configuration adjusts the balance automatically once the user changes either the O2% or the total flow is changed.

Selecting **Flow** shows balance gas l/min as the first quick key and O2 l/min as the second quick key. This gas control/screen configuration allows individual control of the gasses.

1. Push the **Main Menu** key.
2. Select **Screen Setup - Fresh Gas Controls**.
3. Select the desired view.
4. Push the **Normal Screen** key or select **Previous Menu**.

## Alarm setup

Alarm limits, alarm volume, and other alarm settings are adjusted in the Alarm Setup menu. Alarm history is also accessed through this menu. Selecting Default Limits loads the default settings as set by the Super User or the factory defaults if no Super User settings have been entered. Setting Leak Audio to Off silences audio alarms for small leaks.

### Setting alarm limits

1. Push the **Alarms Setup** key.
2. Select **Adjust Settings**.
3. Scroll to the desired alarm.
4. Select alarm limit and change the value.
5. Push the **Normal Screen** key or select **Back**.

### MV/TV Alarms

The MV and TV alarms can be turned off. When the volume alarms are disabled, 'MV/TV alarms Off' appears in the general message field. If the MV and TV alarms are set to **Off** during manual ventilation, the alarms remain off until the Bag/Vent switch is set to ventilator or the **MV/TV Alarms** is set to **On**.

1. Push the **Alarm Setup** key.
2. Select **MV/TV Alarms** and set to **Off**.

### CO2 Audio

Set **CO2 Audio** to **Off** to silence the audio alarms for 'Apnea' CO2 and 'ETCO2 low' when the Bag/Vent switch is set to bag or when the respiratory rate setting is greater than 60 during mechanical ventilation. The visual alarms for 'Apnea' and for 'ETCO2 low' are not affected. The audio alarms remain off until the **CO2 Audio** is set to **On** or the respiratory is less than or equal to 60 during mechanical ventilation.

1. Push the **Alarm Setup** key.
2. Select **CO2 Audio** and set to **Off**.

### Auto MV Limit

MV alarm limits can be calculated automatically for mechanical ventilation when in VCV mode. Set **Auto MV Limit** to **On** for automatic calculations of the MV alarm limits. The MV alarm limits continue to be automatically calculated until **Auto MV Limit** is set to **Off** or until an MV alarm limit is manually adjusted during mechanical ventilation.

If the automatic calculation of the low or high minute volume alarm limit exceeds the allowable limit, the minimum or maximum alarm limit is used.

1. Push the **Alarm Setup** key.
2. Select **Auto MV Limit** and set to **On**.

## Viewing trends

There are three views for patient trends: measured (numerical), settings, and graphical. Trend information is saved every one minute for the most recent 24 hours.

1. Push **Trends**.
2. Select the desired view.
  - An arrow indicates the current trend view.
3. Select **Cursor** to scroll through the current trend view.
4. Push the ComWheel to return the highlight to **Cursor**.
5. Select **Next Page** to view additional parameters.
6. Push the **Normal Screen** key or select **Previous Menu**.



## Alternate O<sub>2</sub> control

**WARNING** The Alternate O<sub>2</sub> control is not an auxiliary source of O<sub>2</sub>.


Alternate O<sub>2</sub> control delivers agent and O<sub>2</sub> through an independent pneumatic path to the selected patient circuit and is connected to the system O<sub>2</sub> supply. Alternate O<sub>2</sub> control activates automatically in the case of certain failures or errors. It can also be activated manually. O<sub>2</sub> flow is 0.5 l/min to 10 l/min as indicated on the flow tube.

**WARNING** When Alternate O<sub>2</sub> control is enabled, flow from the electronic mixer is stopped and the agent concentration is set to off. O<sub>2</sub> is flowing through the Alternate O<sub>2</sub> control to the breathing system. To activate anesthetic agent flow to the breathing system, set the agent to the desired concentration.

The agent delivery cannot be activated in the case of certain electronic or agent delivery failures. The Agent quick key will be blank if agent delivery is unavailable.

The Alternate O<sub>2</sub> control is available approximately 20 seconds after the system is turned on. The minimum flow can be preset to 0.5 to 10 l/min.

1. Push the Alternate O<sub>2</sub> ON/OFF button.
2. The O<sub>2</sub> flow is indicated on the flow tube.
3. Use the flow control to adjust the O<sub>2</sub> flow.
4. Set the agent to the desired concentration.



1. ON/OFF button
2. Flow tube
3. Flow control

Figure 3-4 ▪ Alternate O<sub>2</sub> control

## EZchange Canister (optional)

Push the absorber canister release to activate the EZchange Canister mode. The canister will swing down to the EZchange position. The EZchange Canister mode seals the breathing circuit when the canister holder is down. This permits continued ventilation and rebreathing of exhaled gases while easily replacing the absorber canister.

Systems with EZchange Canister have the following label on the canister holder. When the system is in EZchange position, the message 'CO<sub>2</sub> Absorber Out of Circuit' will be shown in the waveform area on the anesthesia display.




To return to absorber mode, reinsert the canister into the holder and push the canister back up and snap it into absorber position. When the canister is in the absorber position, the exhaled gas flows through the absorber, removing  $\text{CO}_2$ .

**Note** Check the absorber canister to ensure it has side rails. If the canister does not have side rails, it will not work on the EZchange Canister holder.



## Condenser (optional)

Visually check the condenser reservoir daily. If needed, drain the reservoir. Place a container under the reservoir. Push the drain button to empty any water in the condenser reservoir.



1. Drain button
2. Reservoir
3. Condenser

*Figure 3-5 ■ Condenser*

## Passive AGSS (optional)

**WARNING** Always verify the proper operation of any gas scavenging system; ensure the scavenging system is not occluded.

The passive AGSS (Anesthesia Gas Scavenging System) contains both positive and negative pressure relief valves to protect the breathing system and the patient. The outlet is a 30-mm connector on the bottom of the receiver.

There is also a connector that may be used for scavenging the sample from a gas monitor. The male luer inlet connection is located near the 30-mm connector.

Passive AGSS is intended primarily for use in operating room environments which have no active gas extraction system for waste gas disposal. The disposal system generally consists of large diameter tubing directly linking the passive AGSS with the building exterior. The tubing should be as large in diameter and as short as possible for the particular application.

Passive AGSS may also be used with a non-recirculating facility ventilation system for waste gas disposal. The tubing connection from passive AGSS to the non-recirculating facility ventilation system should be an open connection, essentially at atmospheric pressure, such as to an exhaust grill.

## Active AGSS (optional)

**WARNING** Always verify the proper operation of any gas scavenging system; ensure the scavenging system is not occluded.

There are several versions of the optional active AGSS (Anesthesia Gas Scavenging System) available depending on the hospital's type of waste gas disposal system.

Each version has a two-liter reservoir to capture peak exhaust flows that briefly exceed the extract flow. The disposal system normally entrains room air through an air brake (located in a receiver underneath the breathing system), but will spill from this port during extended periods of high exhaust flow. Its effectiveness is limited by the extract flow of the particular active AGSS device.

- The active low flow system is for use with high vacuum disposal systems. It requires a vacuum system capable of a continuous nominal flow of 36 l/min and 300 mmHg (12 inHg) or greater vacuum pressure. A flow indicator on the system indicates when the unit is in operation.
- The active high flow system is for use with low vacuum (blower type) disposal systems. This requires a system capable of providing a continuous nominal flow of 50 l/min. A flow indicator on the system indicates when the unit is in operation.
- Another version is the active adjustable flow. It provides the capability to adjust the flow with a needle valve (located in a receiver underneath the breathing system) and a visual indicator bag which should be properly inflated. It requires a vacuum system capable of a continuous nominal flow of 36 l/min and 300 mmHg (12 inHg) or greater vacuum pressure. A flow indicator on the system indicates when the unit is in operation.
- The active low flow system with a 12.7 mm hose barb connector is for use with low vacuum disposal system. It requires an external Venturi system with flowmeter and 36 l/min extract flow.
- The active low flow system with a 25 mm barb connector is for use with low vacuum disposal systems. It requires an external venturi/ejector system with 36 l/min extract flow. A flow indicator on the system indicates when the unit is in operation.
- The active low flow system with 30 mm ISO taper is for use with low vacuum disposal systems. It requires an external venturi/ejector system with 36 l/min extract flow. A flow indicator on the system indicates when the unit is in operation.

## Connecting active AGSS with a flow indicator

To use the optional active AGSS on a system that has a flow indicator, connect it as follows.

1. Connect the proper hose to the AGSS outlet connector on the bottom of the AGSS underneath the breathing system. Attach the other end to the hospital disposal system.
2. With the AGSS operating, verify that the flow indicator ball on the flow indicator rises to the green zone, indicating adequate flow.

### Note

The ball in the upper red zone indicates excessively high extraction flow. The ball in the lower red zone indicates extraction flow rate is too low or a blocked filter.


3. Complete the tests in the “*Preoperative Tests*” section of this manual.

## Connecting active adjustable AGSS

The active adjustable AGSS option flow rate is limited to 30 l/min with this option.

To use the optional active AGSS installed on the system which uses the three-liter bag as a visual indicator, connect it as follows.

1. Connect a disposal hose to the DISS connector on the needle valve on the bottom of the AGSS (underneath). The hose should be flexible and reinforced to help prevent kinking and crushing.
2. Attach the other end of the hose to the hospital disposal system.
3. Attach the three-liter bag to the 30-mm auxiliary 1 port on the bottom of the AGSS.
4. Use the needle valve to adjust the flow rate to match the amount of gas being scavenged. Use the visual indicator bag when adjusting the flow rate. The bag should remain partially inflated.



5. Complete the tests in the “*Preoperative Tests*” section of this manual.



# 4 Preoperative Checkout

**WARNING** Read each component’s User’s Reference manual and understand the following before using this system:

- All system connections.
- All warnings and cautions.
- How to use each system component.
- How to test each system component.

Before using the system:

- Complete all of the tests in this section.
- Test all other system components.

If a test fails, do not use the equipment. Have a Datex-Ohmeda trained service representative repair the equipment.

<b>In this section</b>	
Every day before your first patient .....	4-2
Before every patient .....	4-3



## Every day before your first patient

- ☐ Check that necessary emergency equipment is available and in good condition.
- ☐ Check that the equipment is not damaged and that components are correctly attached.
- ☐ Check that pipeline gas supplies are connected and cylinders are installed.
- ☐ Check that the amount of liquid agent in the cassette is adequate, and install the agent cassette into the active bay. The cassette is properly inserted when the agent is identified on the display.
- ☐ Check that the breathing circuit is correctly connected, not damaged, and the breathing system contains sufficient absorbent.
- ☐ Turn the System switch to On.
- ☐ Connect scavenging and verify operation.
- ☐ Calibrate the flow sensors.
- ☐ Do a **Machine Check** in the **Checkout** menu.
- ☐ To check additional Aladin cassettes, do an **Agent Delivery** check in the **Checkout** menu.
- ☐ Ensure adequate reserve O<sub>2</sub> supply is available.
- ☐ Check that the ventilator functions correctly:
  - Connect a test lung to the patient Y.
  - Set the ventilator to VCV mode and the settings to TV at 400 ml, RR at 12, I:E at 1:2, Tpause at Off, PEEP at Off, and Pmax at 40.
  - Set the gas flow to the minimum settings.
  - Start a case.
  - Set the Bag/Vent switch to ventilator.
  - Fill the bellows using O<sub>2</sub> flush.
  - Check that mechanical ventilation starts. Check that the bellows inflate and deflate. Check that the display shows the correct ventilator data. Check that there are no inappropriate alarms.
- ☐ Set the appropriate controls and alarms limits for the case.

## Before every patient

**Note** This check does not need to be done before the first case of the day if the “Every day before your first patient check” was done.

- ☐ Check that necessary emergency equipment is available and in good condition.
- ☐ If the cassette was changed since the last case:
  - Check that the amount of liquid agent in the cassette is adequate.
  - Check that the cassette is properly inserted. The cassette is properly inserted when the agent is identified on the display.
- ☐ Check that the breathing circuit is correctly connected, not damaged, and the breathing system contains sufficient absorbent.
- ☐ Leak check the breathing system:
  - Set the Bag/Vent switch to Bag, close the APL valve (set to 70), and occlude the patient Y.
  - Pressurize the breathing system to about 30 cmH<sub>2</sub>O using the O<sub>2</sub> flush button.
  - Ensure that the pressure remains fixed for at least 10 seconds.
- ☐ Check that the ventilator functions correctly:
  - Connect a test lung to the patient Y.
  - Set the ventilator to VCV mode and the settings to TV at 400 ml, RR at 12, I:E at 1:2, Tpause at Off, PEEP at Off, and Pmax at 40.
  - Set the gas flow to the minimum settings.
  - Start a case.
  - Set the Bag/Vent switch to ventilator.
  - Fill the bellows using O<sub>2</sub> flush.
  - Check that mechanical ventilation starts. Check that the bellows inflate and deflate. Check that the display shows the correct ventilator data. Check that there are no inappropriate alarms.
- ☐ Set the appropriate controls and alarms limits for the case.





# 5 Preoperative Tests

## In this section

Inspect the system .....	5-2
Leak < 250 ml .....	5-3
Machine Check .....	5-3
Individual Checks .....	5-5
Positive low-pressure leak test (ACGO systems only) .....	5-8
Aladin cassette installation .....	5-9
Flow sensor calibration .....	5-9

## Inspect the system

**WARNING** The top shelf weight limit is 45 kg (100 lb).

-  Make sure that the breathing circuit is correctly connected and not damaged. Replace the breathing circuit if it is damaged.
-  Do not leave gas cylinder valves open if the pipeline supply is in use. Cylinder supplies could be depleted, leaving an insufficient reserve supply in case of pipeline failure.

Before using the system, ensure that:

- The equipment is not damaged.
- Components are correctly attached.
- The breathing circuit is correctly connected, not damaged, and the breathing system contains sufficient absorbent.
- The Aladin cassette is locked in position and contains sufficient agent.
- Pipeline gas supplies are connected and the pressures are correct.
- Cylinder valves are closed.
- Models with cylinder supplies have a cylinder wrench (key) attached to the system.
- Models with cylinder supplies have a reserve supply of O<sub>2</sub> connected to the machine during system checkout.
- The necessary emergency equipment is available and in good condition.
- Equipment for airway maintenance, manual ventilation, tracheal intubation, and IV administration is available and in good condition.
- Applicable anesthetic and emergency drugs are available.
- If an optional O<sub>2</sub> flowmeter is present, ensure there is adequate flow.
- If an optional suction regulator is present, ensure there is adequate suction.
- The casters are not loose and the brakes are set and prevent movement.
- The power cord is connected to an electrical outlet. The mains indicator comes on when AC power is connected. If the indicator is not on, the system does not have mains (electrical) power. Use a different outlet, close the circuit breaker, or replace or connect the power cable.

## Leak < 250 ml

The **Leak < 250 ml** setting is used during the circuit leak check portion of the checkout procedures. This check tests for leaks in the machine, breathing circuit, patient circuit, and manual bag. The default setting is **No**.

**Note** Extraction of gas by external gas monitors may cause failure of the leak checks during tests.

**No** When **No** is selected, the leak test will pass for leaks below 250 ml at 3 kPa (30 cmH<sub>2</sub>O) pressure with no user interaction required. For leaks between 250 ml and 750 ml, the user can fix the leak and rerun the test or accept the leak and continue. For leaks above 750 ml, the test will fail and the user must fix the leak and rerun the test.

**Yes** Set to **Yes** to measure small leaks above 100 ml during the checkout procedures. Selecting **Yes** will display the measured leak at 3 kPa (30 cmH<sub>2</sub>O) pressure and result in the test taking somewhat longer.

## Machine Check

Perform **Machine Check** at the start of each day. The machine check runs automatically and beeps to indicate when it is finished or if interaction is required.

The **Machine Check** does a **Machine Check - System** check, **Machine Check - Circuit** check, and a **Machine Check - Circuit O<sub>2</sub>** cell check (if circuit O<sub>2</sub> cell is present). When one of the checks is completed, the next check begins.

1. Turn the System switch to On.
2. Select **Machine Check** and follow the instructions.
3. If a check fails, follow the instructions to perform a recheck or accept the results.
4. When the **Machine Check** is completed, start a case.

**Note** In case of a patient emergency, the **Machine Check** may be bypassed by selecting **Start Case** from the **Checkout** menu. The general message 'Please Do Checkout' is displayed if a **Machine Check** is not completed with passing results within 24 hours.

## Machine Check - System

The **Machine Check-System** checks the Bag/Vent switch, proper gas supply pressures, ventilator operation and leak, battery and electrical power, circuit compliance, flow control operation, and vaporizer operation. This is a two-step check.

1. Set the Bag/Vent switch to Vent.
2. Open the patient Y.
3. (ACGO option only.) Set the ACGO switch to Circle.
4. Select **Start**. The display shows the checks being run.
  - The system beeps when this portion of the check is done.
  - The results are shown on the display.
5. Make sure the bellows is fully collapsed.
6. Occlude the patient Y.
7. Select **Continue**. The display shows the checks being run.
8. When the check passes, the next check starts.

## Machine Check - Circuit

The **Machine Check-Circuit** checks the Bag/Vent switch, proper gas supply pressures, airway pressure measurement transducer, APL valve, and manual circuit leak.

1. Occlude the patient Y.
2. Set Bag/Vent switch to Bag.
3. Set the APL valve halfway between 30 and 70.
4. (ACGO option only.) Set the ACGO switch to Circle.
5. Select **Start**. The display shows the checks being run.
  - The system beeps when the check is done.
  - The results are shown on the display.
6. When the check passes, the next check starts.

## Machine Check - Circuit 02

The **Machine Check-Circuit 02** check measures the O<sub>2</sub>%.

1. Open the patient Y.
2. Set the Bag/Vent switch to Vent.
3. (ACGO option only.) Set the ACGO switch to Circle.
4. The display will show the O<sub>2</sub>%. Do not select **Done** when 21 is first displayed. Allow the reading to stabilize, then select **Done**. Calibrate the O<sub>2</sub> cell if necessary.

## Individual Checks

Individual checks allow the user to perform any combination of single checks. These checks are helpful if there is a specific problem/alarm and the user wishes to test only that portion of the system.

The checks do not automatically move on to the next check. After completing a check, do another check or start a case. All of the checks must be performed at least once within every 24-hour period. If a check fails, follow the instructions to perform a recheck or accept the results.

### System

The **System** check checks the Bag/Vent switch, proper gas supply pressures, ventilator operation and leak, battery and electrical power, circuit compliance, flow control operation, and vaporizer operation. This is a two-step check.

1. Set the Bag/Vent switch to Vent.
2. Open the patient Y.
3. (ACGO option only.) Set the ACGO switch to Circle.
4. Select **Start**. The display shows the checks being run.
  - The system beeps when this portion of the check is done.
  - The results are shown on the display.
5. Make sure the bellows is fully collapsed.
6. Occlude the patient Y.
7. Select **Continue**. The display shows the checks being run.
8. When the check passes, select **Back**.
9. Select another check or select **Start Case** to go to the **Start Case** menu.

### Circuit

The **Circuit** check checks the Bag/Vent switch, proper gas supply pressures, airway pressure measurement transducer, APL valve, and manual circuit leak.

1. Occlude the patient Y.
2. Set Bag/Vent switch to Bag.
3. Set the APL valve halfway between 30 and 70.
4. (ACGO option only.) Set the ACGO switch to Circle.
5. Select **Start**. The display shows the checks being run.
6. When the check passes, select **Back**.
7. Select another check or select **Start Case** to go to the **Start Case** menu.



## Circuit O<sub>2</sub> Cell

The **Circuit O<sub>2</sub> cell** check measures the O<sub>2</sub>%.

1. Open the patient Y.
2. Set the Bag/Vent switch to Vent.
3. (ACGO option only.) Set the ACGO switch to Circle.
4. The display will show the O<sub>2</sub>%. Do not select **Done** when 21 is first displayed. Allow the reading to stabilize, then select **Done**. Calibrate the O<sub>2</sub> cell if necessary.
5. Select another check or select **Start Case** to go to the **Start Case** menu.

## Low P Leak

The positive pressure **Low P Leak** check measures machine leaks before the breathing system and between the gas mixer and the common gas outlet. It measures low pressure pneumatic leaks with a pass/fail limit of 50 ml.

1. Occlude the inspiratory (right-hand) port.
2. Select **Start**.
3. The display shows the checks being run. The system beeps when the check is done.
4. Open the inspiratory port and reconnect the breathing circuit.
5. Select another check or select **Start Case** to go to the **Start Case** menu.

## Low P Leak (machines with ACGO)

The negative low P leak check measures machine leaks before the breathing system and between the gas mixer and the common gas outlet. It measures low pressure pneumatic leaks with a pass/fail limit of 50 ml.

1. Make sure the ACGO switch is set to ACGO.
2. Attach the squeeze bulb to the ACGO outlet.
3. Squeeze (collapse) the bulb.
4. If the bulb inflates in < 30 seconds, select **Fail**.
5. If the bulb remains collapsed, select **Pass**.
6. Remove the squeeze bulb from the ACGO outlet.

## Agent Delivery

The agent delivery check checks the agent delivery system and the cassette. This check will release agent to the circuit.

1. Insert a cassette, connect a patient circuit, and connect scavenging.
2. Set the Bag/Vent switch to Vent.
3. (ACGO option only.) Set the ACGO switch to Circle.
4. Occlude the patient Y.
5. Select **Start**.
6. The display shows the checks being run. The system beeps when the check is done.
7. Repeat the check once for each cassette
8. When the check passes, select **Back**.
9. Select another check or select **Start Case** to go to the **Start Case** menu.

## Positive low-pressure leak test (ACGO systems only)

**Note** For ACGO machines, perform either a negative **Low P Leak** check in the **Checkout** menu or positive low-pressure leak check depending on local requirements.

**CAUTION** Do a positive-pressure leak test at the ACGO port only.

1. Connect the leak test device to the ACGO port with the positive-pressure leak test adapter. Push the adapter into the ACGO port throughout the test to get a good seal.
2. Fully open the needle valve on the test device. Keep the test device flow tube vertical for accurate results.

**CAUTION** If the needle valve is not fully open, this test can damage the pressure gauge on the test device.

3. Turn the ACGO switch to the ACGO position.
4. Set Air as the balance gas.
5. Set O<sub>2</sub> to 100%.
6. Set the O<sub>2</sub> flow so that total flow through the flowmeter on the test device is 0.5 l/min.
7. Make sure that the pressure gauge on the test device reads zero.
8. Close the needle valve on the test device until the test gauge reads 20 kPa (3 psi) (BSI) or 3 kPa (0.4 psi) (ISO).
9. If the flow through the test device is less than 0.45 l/min (ISO) or 0.4 l/min (BSI), there is a low pressure leak in the anesthesia machine. See section 7, "Alarms and troubleshooting."

**WARNING** Agent mixtures from the low-pressure leak test stay in the system. Always flush the system with O<sub>2</sub> after the low-pressure leak test (1 l/min for one minute).

10. Remove the adapter and leak test device.
11. Set the O<sub>2</sub> flow to 1 l/min and continue flow for one minute.

## Aladin cassette installation

1. Check that the cassette is filled to the appropriate level.
2. Insert cassette into the active bay until a click is heard, ensuring the cassette is in position.
3. Turn the lock on the handle to the horizontal position (Aladin<sub>2</sub> cassettes only).
4. The cassette is properly inserted when the agent is identified on the display. Make sure that the displayed agent matches the cassette.

**Note** Store the cassettes in the cassette bay when they are not in use.

## Flow sensor calibration

**Important** Recalibrate the flow sensors if the room temperature changes by more than 5 °C.

1. Set the Bag/Vent switch to Bag.
2. Remove the flow sensor module.
3. Wait for 'No insp flow sensor' and 'No exp flow sensor' alarms to occur.
4. Reinsert the flow sensor module. Wait for alarms to clear.
5. Start mechanical ventilation when ready.



# 6 Airway Modules

<b>In this section</b>	Airway modules . . . . .	6-2
	Parameters setup . . . . .	6-5
	Automatic agent identification . . . . .	6-6
	Calibration . . . . .	6-7


## Airway modules

The optional compact airway modules measure and monitor gases delivered to the patient and exhaled through the breathing circuit. The modules consist of an infrared sensor for measuring CO<sub>2</sub>, N<sub>2</sub>O, and anesthetic agents; a paramagnetic O<sub>2</sub> sensor; and a gas sampling system with the D-fend water separation system.


Systems with both an airway module and an O<sub>2</sub> cell will display the patient inspired O<sub>2</sub> value obtained from the airway module.

Respiratory rate is the frequency of peak (end tidal) CO<sub>2</sub> measurements per minute. A breath is defined as a change in the CO<sub>2</sub> signal that exceeds 1% (8 mmHg). All concentrations are measured and displayed breath by breath.

**WARNING** Remove the airway sampling line from the patient's airway and seal the sample port while nebulized medications are being delivered. Nebulized medications interfere with accurate gas reading.

 If the **Data Source** is set to **Patient** in the spirometry **Setup** menu, no Paw, Flow, Agent, or CO<sub>2</sub> waveforms or numeric information is displayed during the airway gas module warm-up period (approximately 2 minutes). Set the **Data Source** to **Vent** to display the waveforms and numeric information during the airway gas module warm-up period.


**CAUTION** Use only cables and accessories approved by Datex-Ohmeda. Other cables and accessories may damage the system or interfere with measurement. Single-use accessories are not designed to be reused.

 Strong scavenging suction on the monitor exhaust port may change the operating pressure of the monitor and cause inaccurate readings or internal damage.

Use only airway modules that have anesthetic agent monitoring and O<sub>2</sub> monitoring on this system. The following modules can be used on this system: E-CAiO, E-CAiOV, E-CAiOVX, M-CAiO, M-CAiOV, and M-CAiOVX. (E series modules must be software version 4.5 and above. M series modules must be software revision 3.2 and above.)

Letters in the name of the airway modules stand for:

- E = plug-in gas module
- M = plug-in gas module
- C = CO<sub>2</sub> and N<sub>2</sub>O
- A = anesthetic agents
- i = agent identification
- O = patient O<sub>2</sub>
- V = patient spirometry
- X = Gas exchange



1. D-fend water trap
2. Sampling line connector
3. Water trap latch
4. Reference gas inlet
5. Sample gas outlet
6. Cooling fan

Figure 6-1 ▪ Compact airway module




## Connection to a patient

1. Check that the airway gas module is installed.
2. Check that the airway adapter connections are tight and that the adapter is correctly installed.
3. Check that the water trap container is empty and properly attached.

### WARNING

Before connecting the exhaust line to the sample gas outlet on the compact airway module, ensure the other end is connected to the sample gas return port on the anesthesia machine. Incorrect connections may cause patient injury.

4. Connect the exhaust line from the sample gas return port to the sample gas outlet if N<sub>2</sub>O or volatile agents are used. (Refer to the “Parts” section in Part 2 of the User’s Reference manual.)
5. Attach the gas sampling line to the sampling line connector on the water trap.
6. Turn the system on.
7. Connect the sampling line to the airway adapter. Take the gas sample as close to the patient’s airway as possible. Position the adapter’s sampling port upwards to prevent condensed water from entering the sampling line.



1. Airway module
2. Gas sampling line
3. Airway adapter with sampling line connector

Figure 6-2 ▪ Airway gases setup with compact airway module

## Parameters setup

Use the **Parameters Setup** menu to change the monitoring settings of the data source, CO<sub>2</sub>, O<sub>2</sub>, agent, and spirometry.

### Data source

Several monitoring parameters can be obtained from the ventilator or the airway module. Set to **Patient** or **Vent** to select the primary source for information. Information that is retrieved from the airway module is identified with the module data indicator.



Figure 6-3 ▪ Airway module data indicator

- If **Patient** is selected, the airway module will be the first source for information if a D-lite or Pedi-lite sensor is properly connected. If the sensor is not properly connected, but the airway module is installed the waveform data is not displayed. If information is not available through the airway module, information will come from the internal ventilator sensors.
- If **Vent** is selected, the internal sensors of the ventilator will be the first source for information.

### CO<sub>2</sub> setup

Change the size of the CO<sub>2</sub> waveform by changing the scale height. The scale units (% , kPa, and mmHg) are set by an authorized user. Select **CO<sub>2</sub> Alarm** or **Resp Rate Alarm** to access and change the alarm limits for the corresponding alarm.

### O<sub>2</sub> setup

Select **O<sub>2</sub> Alarm** to access and change the O<sub>2</sub> alarm limits.

### Agent setup

Change the size of the agent waveform by changing the scale height. Select **Agent Alarm** to access and change the agent alarm limits.

### Spirometry setup

Select **Spirometry Setup** to access the **Spirometry Setup** menu. Change the size of the Paw and Flow waveforms by changing the corresponding scale heights. Select **Paw Alarm** or **MV Alarm** to access and change the alarm limits for the corresponding alarm.

## **Automatic agent identification**

Airway modules with agent identification will automatically identify and select Halothane, Enflurane, Isoflurane, Sevoflurane, and Desflurane. The inspiratory and expiratory concentrations of the agent appear in the number field or the agent waveform field if selected.

Minimum concentration for the identification is 0.15% volume. The agent selection remains active even if the concentration decreases below 0.15% volume during the case.

Automatic agent identification operates after the normal warm up of the gas module (approximately five minutes).

# Calibration

Calibrate airway modules once every six months or whenever there are indications of errors in the gas readings. Use a Datex-Ohmeda calibration gas and regulator to calibrate the modules. See the “Parts” section in Part 2 of the User’s Reference manual for stock numbers of the calibration gas and regulator.

**WARNING** Only use Datex-Ohmeda calibration gas. Do not use any other calibration gases or the calibration will not succeed.

During gas calibration, % units are used for CO<sub>2</sub> regardless of selected measuring units.

1. Turn on the power. Let the module warm up for 30 minutes before starting calibration.
2. Attach the regulator to the calibration gas cylinder.
3. Attach a new sampling line to the water trap. Connect the loose end of the sampling line to the regulator on the calibration gas cylinder.
4. Push the **Main Menu** key.
5. Select **Calibration**.
6. Select **Airway Gas**.
7. Wait until ‘Feed Gas’ appears after each gas.
8. Open the regulator until the gauge reads between 5 to 7 psi and feed calibration gas until the message ‘OK’ or ‘Adjust’ appears.
  - If an error occurs during calibration or if no gas is fed, ‘Calib Error’ appears after each gas. Push the ComWheel to perform a new calibration.
9. If adjustments are needed:
  - Do not close the regulator until all the adjustments have been made.
  - Select the gas to be adjusted and press the ComWheel.
  - Use the ComWheel to change the value until it matches the calibration gas cylinder value. Push the ComWheel to confirm the change.
  - Repeat for each gas requiring adjustment.



# 7 Vaporizer

<b>In this section</b>	Vaporizer . . . . .	7-2
	Installing cassettes . . . . .	7-5
	Cleaning . . . . .	7-5
	Service . . . . .	7-5
	Draining halothane cassettes . . . . .	7-5
	Filling Aladin <sub>2</sub> cassettes . . . . .	7-6
	Filling Aladin cassettes . . . . .	7-11

## Vaporizer

The electronically controlled vaporizer consists of the internal electronic control unit and the Aladin agent cassette. The agent cassettes are color coded, have indexed filling ports, and are magnetically coded for each agent. The electronic control unit governs the flow through the agent cassette and the agent concentration in the fresh gas flow.


Both the Aladin<sub>2</sub> and the Aladin cassettes can be used on this system. Aladin<sub>2</sub> cassettes can be easily recognized by the lock on the front handle and the liquid level indicator.

Remove the cassette from the active bay when not administering anesthetic agent. Before removing a cassette, ensure that agent delivery is off. Store cassettes in the cassette bay when they are not in use. Store the cassettes and agent at the same temperature as the system.

### Aladin<sub>2</sub> cassette variants

There are three types of Aladin<sub>2</sub> cassette filler systems. Halothane, enflurane, and isoflurane use a color-coded, Easy-Fil mechanism. Sevoflurane cassettes are available with a color-coded, Easy-Fil or Quik-Fil mechanism. The desflurane cassettes have a filling mechanism that is compatible with Saf-T-Fil desflurane bottles.

Anesthetic agent	Filling system	Color code
Halothane	Easy-Fil	Red
Enflurane	Easy-Fil	Orange
Isoflurane	Easy-Fil	Purple
Sevoflurane	Easy-Fil or Quik-Fil	Yellow
Desflurane	Compatible with Saf-T-Fil	Blue



1. Handle with release trigger
2. Lock
3. Liquid level indicator
4. Agent filling port


Figure 7-1 ▪ Aladin<sub>2</sub> cassette

## Aladin cassette variants

There are three types of Aladin cassette filler systems. Halothane, enflurane, and isoflurane use color-coded, keyed fillers. Sevoflurane cassettes are available with a color-coded, keyed filler or a Quik-Fil mechanism. The desflurane cassettes have a filling mechanism that is compatible with Saf-T-Fil desflurane bottles.

**Note** Sevoflurane Quik-Fil Aladin cassettes are not available in all countries.

Anesthetic agent	Filling system	Color code
Halothane	Keyed	Red
Enflurane	Keyed	Orange
Isoflurane	Keyed	Purple
Sevoflurane	Keyed or Quik-Fil	Yellow
Desflurane	Compatible with Saf-T-Fil	Blue




AB.60.004

1. Handle with release trigger
2. Liquid level indicator
3. Agent filling port
4. Lock and fill wheel

*Figure 7-2 ▪ Aladin cassette for halothane, enflurane, isoflurane, and sevoflurane with keyed filler system*






1. Handle with release trigger
2. Liquid level indicator
3. Agent filling port

*Figure 7-3 ■ Aladin cassette for sevoflurane with Quik-Fil system*

**Note** Sevoflurane Quik-Fil Aladin cassettes are not available in all countries.



1. Handle with release trigger
2. Liquid level indicator
3. Agent filling port

*Figure 7-4 ■ Aladin cassette for desflurane with Saf-T-Fil compatible filler system*

## Installing cassettes

1. Check that the cassette is filled to the appropriate level.
2. Insert cassette into the active bay until a click is heard, ensuring the cassette is in position.
3. Turn the lock on the handle to the horizontal position (Aladin<sub>2</sub> cassettes only).
4. The cassette is properly inserted when the agent is identified on the display. Make sure that the displayed agent matches the cassette.

**Note** Store the cassettes in the cassette bay when they are not in use.

## Cleaning

Remove the cassette from the machine. Clean the cassette surface with a cloth moistened in mild soap solution.

**CAUTION** Do not wipe Aladin cassettes with alcohol-based detergents. This may damage the surface of the cassette.

## Service

All types of Aladin cassettes must be emptied before shipping. Contact a Datex-Ohmeda trained service representative for draining of cassettes. Package the cassettes in suitable wrapping when shipping.

## Draining halothane cassettes

The decomposition of halothane causes the release of halides, which may corrode metal components particularly in the presence of moisture. A preservative added to halothane by its manufacturers, to impede decomposition, can leave a residue which may cause cassette components to stick.

- Have the halothane cassettes drained every two weeks.
- If halothane is used infrequently, have the halothane cassette drained after use.
- Do not flow gas through the cassette to dry out the cassette.

**Note** Contact a Datex-Ohmeda trained service representative for draining of cassettes.

## Filling Aladin<sub>2</sub> cassettes

Unlock and remove the cassette from the machine before filling the cassette. Ensure the cassette remains in a horizontal position during filling. If the cassette is tilted, the flow of liquid is shut off to prevent overfilling.

Ensure that the valve pins on the back of the cassette do not come in contact with the machine or any other object that could depress the valve pins and let gas out of the cassette. Observe the liquid level indicator of the cassette during filling. When the liquid level reaches the full mark, stop filling.

**WARNING** Always remove the Aladin cassette from the machine and put it on a horizontal surface before filling. Never try to fill a cassette while it is in the machine.

### Filling with Easy-Fil system

The Easy-Fil system consists of three elements:

- The bottle collar.
- The bottle adapter.
- The agent filling port.




Figure 7-5 ▪ Aladin<sub>2</sub> cassette with Easy-Fil system

**WARNING** Do not open or press the cassette filling port or the gas connection valves with fingers or any kind of instrument. Anesthetic agent liquid or gas may squirt into the air.

1. Remove the cassette from the machine, and place it on a horizontal surface. Ensure that the valve pins on the back of the cassette do not come in contact with the machine or any other object.
2. Align the notches on the bottle adapter to the agent bottle collar, and tighten the adapter onto the bottle.

**WARNING**

To avoid exposure to anesthetic agent during filling, ensure that the bottle adapter is properly connected to the bottle.

3. Remove the filling port cap from the cassette by turning the cap counterclockwise.
4. Align the bottle adapter keys with the index slots in the filling port.
5. Push the agent bottle firmly into the filling port. Let the agent flow into the cassette.
6. Watch the liquid level indicator. When the cassette is full, remove the bottle from the filling port.
7. Put the cap back on the filling port.
8. Put the cap back on the agent bottle.

## Filling with Quik-Fil system




Figure 7-6 ▪ Aladin<sub>2</sub> cassette with Quik-Fil system

**WARNING** Do not open or press the cassette filling port or the gas connection valves with fingers or any kind of instrument. Anesthetic agent liquid or gas may squirt into the air.

1. Remove the cassette from the machine, and place it on a horizontal surface. Ensure that the valve pins on the back of the cassette do not come in contact with the machine or any other object.
2. Remove the yellow protective cap from the anesthetic agent bottle. Check that the filling mechanism on the bottle is not damaged.
3. Remove the filling port cap from the cassette by turning the cap counterclockwise.
4. Align the bottle nozzle keys with the index slots in the filling port.
5. Push the agent bottle firmly into the filling port. Let the agent flow into the cassette.
6. Watch the liquid level indicator. When the cassette is full, remove the bottle from the filling port.
7. Put the cap back on the filling port.
8. Put the cap back on the agent bottle.

## Filling with Saf-T-Fil bottle




Figure 7-7 ▪ Aladin<sub>2</sub> cassette with Saf-T-Fil mechanism

**CAUTION** Do not attempt to fill the cassette with desflurane that is warmer than 26 °C/78 °F. The overfill prevention system may lock and prevent cassette filling.

**WARNING** Do not store filled or partially filled desflurane cassettes above the normal working temperature of 35 °C/95 °F. Storage at high temperature may cause the overpressure valve to vent desflurane vapor into the ambient air.

⚠ Do not open or press the cassette filling port or the gas connection valves with fingers or any kind of instrument. Anesthetic agent liquid or gas may squirt into the air.

1. Remove the cassette from the machine, and place it on a horizontal surface. Ensure that the valve pins on the back of the cassette do not come in contact with the machine or any other object.
2. Remove the cap from the desflurane bottle. Check that the o-ring is correctly fitted on the bottle nozzle.


3. Insert the nozzle into the filling port, and push the bottle firmly against the spring pressure until it stops.
4. Turn the bottle upwards while keeping it firmly inserted.
5. Watch the liquid level indicator. Do not leave the bottle unattended while it is attached to the cassette.
6. When the liquid level reaches the full mark, lower the bottle to stop the flow of liquid.
7. Remove the bottle from the filling port.
8. Put the cap back on the agent bottle.

## Filling Aladin cassettes

Unlock and remove the cassette from the machine before filling the cassette. Ensure the cassette remains in a horizontal position during filling. If the cassette is tilted, the flow of liquid is shut off to prevent overfilling.

Ensure that the valve pins on the back of the cassette do not come in contact with the machine or any other object that could depress the valve pins and let gas out of the cassette. Observe the liquid level indicator of the cassette during filling. When the liquid level reaches the full mark, stop filling.

**WARNING** Always remove the Aladin cassette from the machine and put it on a horizontal surface before filling. Never try to fill a cassette while it is in the machine.

 Inaccurate anesthetic agent dosages may temporarily occur after filling if the temperature of the liquid is considerably different from the normal operating temperature.

### Filling with keyed filler system




Figure 7-8 ▪ Filling Aladin cassette with keyed filler system

**WARNING** Do not open or press the cassette filling port or the gas connection valves with fingers or any kind of instrument. Anesthetic agent liquid or gas may squirt into the air.

1. Remove the cassette from the machine, and place it on a horizontal surface. Ensure that the valve pins on the back of the cassette do not come in contact with the machine or any other object.



2. Align the notches on the bottle adapter to the agent bottle collar, and tighten the adapter onto the bottle.

**WARNING** To avoid exposure to anesthetic agent during filling, ensure that the bottle adapter is properly connected to the bottle.

3. Insert the keyed filler into the agent filling port on the cassette. Make sure that the square end of the filler is flush with the cassette body as shown.

**CAUTION** Ensure that the keyed filler is properly inserted into the filling port on the cassette. If air leaks into the system, the overfill protection mechanism of the cassette will not function properly. Overfilling the cassette may result in vaporizer delivery failure.

4. Lock the keyed filler by turning the lock and fill wheel fully clockwise.
5. Turn the bottle upside down. Check that the liquid flows evenly into the cassette. If it does not, tighten the wheel until the liquid starts to flow.
6. Watch the liquid level indicator. When the cassette is full, lower the bottle to stop the flow of liquid.
7. Turn the lock and fill wheel a half turn counterclockwise to close the filling port. Wait for the residual anesthetic agent to trickle back into the bottle.
8. Turn the lock and fill wheel fully counterclockwise, and remove the keyed filler from the filling port.
9. Remove the bottle adapter from the anesthetic agent bottle. Put the cap back on the agent bottle.

## Filling with Quik-Fil system




Figure 7-9 ▪ Filling Aladin cassette with Quik-Fil system

### WARNING

Do not open or press the cassette filling port or the gas connection valves with fingers or any kind of instrument. Anesthetic agent liquid or gas may squirt into the air.

1. Remove the cassette from the machine, and place it on a horizontal surface. Ensure that the valve pins on the back of the cassette do not come in contact with the machine or any other object.
2. Remove the yellow protection cap from the anesthetic agent bottle. Check that the filling mechanism of the bottle is not damaged.
3. Remove the filling port cap from the cassette by turning the cap counterclockwise.
4. Insert the nozzle of the bottle into the filling port.
5. Push the agent bottle firmly into the filling port. Let the agent flow into the cassette.
6. Watch the liquid level indicator. When the cassette is full, remove the bottle from the filling port.
7. Put the cap back on the filling port.
8. Put the cap back on the agent bottle.

## Filling with Saf-T-Fil bottle




Figure 7-10 ▪ Filling Aladin cassette with Saf-T-Fil bottle

**CAUTION** Do not attempt to fill the cassette with desflurane that is warmer than 26 °C/78 °F. The overfill prevention system may lock and prevent cassette filling.

**WARNING** Do not store filled or partially filled desflurane cassettes above the normal working temperature of 35 °C/95 °F. Storage at high temperature may cause the overpressure valve to vent desflurane vapor into the ambient air.

⚠ Do not open or press the cassette filling port or the gas connection valves with fingers or any kind of instrument. Anesthetic agent liquid or gas may squirt into the air.

1. Remove the cassette from the machine, and place it on a horizontal surface. Ensure that the valve pins on the back of the cassette do not come in contact with the machine or any other object.
2. Remove the cap from the desflurane bottle. Check that the o-ring is correctly fitted on the bottle nozzle.
3. Insert the nozzle into the filling port, and push the bottle firmly against the spring pressure until it stops.
4. Turn the bottle upwards while keeping it firmly inserted.
5. Watch the liquid level indicator. Do not leave the bottle unattended while it is attached to the cassette.

6. When the liquid level reaches the full mark, lower the bottle to stop the flow of liquid.
7. Remove the bottle from the filling port.
8. Put the cap back on the agent bottle.



# 8 Alarms and Troubleshooting

**CAUTION** No repair should ever be attempted by anyone not having experience in the repair of devices of this nature. Refer to “*Repair policy*” in the “*User Maintenance*” section in Part 2 of the User’s Reference manual.

**WARNING** If an alarm occurs, safeguard the patient first before performing troubleshooting or doing repair procedures.

<b>In this section</b>	Alarms .....	8-2
	List of alarms .....	8-3
	Alarm ranges and default values .....	8-10
	Alarm tests .....	8-11
	Breathing system problems .....	8-12
	Electrical problems .....	8-13
	Pneumatic problems .....	8-13

## Alarms

Alarms are divided into technical alarms and parameter alarms. These alarms may be high priority, medium priority, or informational. When an alarm occurs during a case, an alarm tone sounds and the alarm message is displayed in the alarm message field.

Technical alarms result from a technical problem and occur whether or not a patient is connected to the system. Parameter alarms are calculated limits and limits that are set by the user on the **Alarm Setup** menu. Parameter alarms occur only during a case.

Silencing an alarm stops the audible tone for 120 seconds. Pushing the **Silence Alarms** key when no medium or high priority alarms are active suspends audible alarm tones for 90 seconds.

### Alarm priorities

Alarm priority is indicated by the color of the alarm message and the alarm LED located next to the alarm silence button.

High-priority alarm messages appear in white text on a red background. During a high-priority alarm, the red LED flashes.

Medium-priority alarm messages appear in yellow text on a gray background. During a medium-priority alarm, the yellow LED flashes.

Informational alarms appear in white text on a gray background. During an informational alarm, the yellow LED is on solid.

When a high-priority alarm is active at the same time as a medium-priority alarm or an informational alarm, the red and yellow LEDs flash. When a high-priority alarm is active at the same time as an informational alarm, the red LED flashes and the yellow LED is on solid. When a medium-priority alarm is active at the same time as an informational alarm, the yellow LED flashes. Pushing the **Silence Alarms** key shuts off the LED until the end of the alarm silence.

### Display changes during alarms

Messages may appear in the waveform field during some alarms. If more than one alarm has a message, the message for the highest priority alarm is displayed. The message is removed when the alarm is resolved.

Messages for high-priority alarms use red text. Messages for medium-priority alarms use yellow text. Informational messages use white text.

Some device related alarm, such as 'Reverse flow' and 'No insp flow sensor,' will de-escalate priority when the alarm is acknowledged by pushing the **Silence Alarms** key.

When the O<sub>2</sub> pipeline supply pressure drops to less than 252 kPa, the lower-right corner of the display toggles between the O<sub>2</sub> supply information and the respiratory data.

Some patient parameter alarms, such as 'Ppeak high and FiO2 Low, will latch when the alarm condition is corrected. When an alarm is latched, it is displayed in white text on a black background. The parameter box will stop flashing. The alarm will remain in this condition until it is acknowledged by pushing the **Silence Alarms** key or until the alarm re-occurs. When the alarm is acknowledged, it is removed from the screen. If an alarm has latched and the alarm re-occurs before it is acknowledged, the alarm will revert to an active state.

## Battery indicator

The color and fill amount of the battery in use symbol indicates the amount of battery power remaining. Green indicates greater than 10 minutes of battery power remaining. Yellow indicates less than 5 minutes battery power remaining.

## Internal failure

'Internal problem prevents normal operation' will show on the display during a software or hardware failure that requires service. If this message occurs, contact a Datex-Ohmeda trained service representative.

## List of alarms

If the corrective action does not resolve the alarm message, contact a Datex-Ohmeda trained service representative

Circuit pressures and volumetric flows are measured by the ventilator and airway gas module. If the **Data Source** is set to **Patient**, the displayed waveforms and numeric information are measured by the airway gas module. Although not displayed, the ventilator measurements continue and if a measured value violates an alarm setting, the appropriate alarm occurs. The value highlighted in the parameter numeric box may not appear to have violated the alarm setting. Changing the **Data Source** to **Vent** will display the circuit pressures and volumetric flows measured by the ventilator.

Message	Priority	Cause	Action
AA, CO2 monitoring not connected	Medium	External Gas Monitor setting is set to No on the Install/Service menu.	Set External gas monitor setting to Yes if the system uses a stand-alone monitor for O2, AA, and CO2.
Adjust low MV limit	Medium	Low MV limit is off in SIMV/PSV, SIMV-PC, or PSVPro modes.	Increase Low MV alarm limit to improve patient disconnection detection.
Agent output not accurate. Schedule service.	Medium	Vaporizer detected a flow meter temperature sensor failure and is using a default temperature value.	Contact a Datex-Ohmeda trained service representative.
Air pressure low. Increase O2%.	High	Air pipeline pressure is less than 252 kPa (36 psi) and the air cylinder pressure dropped below 2633 kPa (381 psi) for one second.	Ensure the air pipeline and cylinder are properly connected.



Message	Priority	Cause	Action
Air supply pressure low	Medium	Air pipeline pressure is less than 252 kPa (36 psi) and the air cylinder pressure dropped below 2633 kPa (381 psi) for one second.	Ensure the air pipeline and cylinder are properly connected.
Apnea	Medium	Apnea time delay (10-30 seconds) has passed without a measured breath.	Check for leaks in the patient circuit. Check for patient disconnection.
Apnea >120 s	High	Apnea time exceeds 120 seconds.	Check for leaks or blockages in the breathing circuit. Ensure the Bag/Vent switch is in the Vent position. Check the patient.
Backup Mode active	Informational	No spontaneous breaths in set period of time (Backup Time) and 30 seconds have passed since starting PSVPro mode.	Select a new ventilation mode.
Breathing System loose	Informational	The breathing system is not latched.	Push the breathing system onto the frame and ensure it latches.
Calibrate flow sensors	Informational	Flow calibration failure.	Calibrate the flow sensors.
Calibrate O <sub>2</sub> sensor	Informational	Calibration failure or O <sub>2</sub> > 110%.	Calibrate the O <sub>2</sub> cell. Replace the O <sub>2</sub> cell if necessary.
Calibrate, dry, or replace flow sensors	Informational	Patient volume mismatch occurred.	Calibrate, dry or replace the flow sensors. Start a new case.
Cannot deliver agent setting at set flow	Informational	Output flow ≥ 6.0 l/min for > 10 seconds.	Reduce the flow.
Cannot identify cassette	Medium	The identification coding of the cassette cannot be interpreted.	Exchange the cassette. Contact a Datex-Ohmeda trained service representative to have the faulty cassette serviced.
Cannot monitor Air pipeline	Medium	Air Pipeline pressure is invalid.	Check pipeline supply pressure.
Cannot monitor gas supplies	Informational	Hardware failure.	Contact a Datex-Ohmeda trained service representative.
Cannot monitor O <sub>2</sub> pipeline	Medium	O <sub>2</sub> Pipeline pressure is invalid.	Check pipeline supply pressure.
Cannot read gas supply pressures	Medium	Transducer failure.	Contact a Datex-Ohmeda trained service representative.
Cassette overfilled, replace cassette	Medium	Agent level sensor indicates the cassette is overfilled.	Exchange the cassette with one filled to a proper level. Contact a Datex-Ohmeda trained service representative to drain the cassette.
Check agent level	Informational	Cassette reporting a value of 1/4 or empty.	Fill the cassette.
Check cassette. Set agent.	Medium	Cassette pressure is out of range. Agent flow interrupted or control failure. System or cassette temperature is out of range. Agent setting is Off.	Try another cassette. Remove and reinsert the cassette. Ensure the room temperature is within specified range. Set the agent to the desired concentration.
Check circuit connections	Medium	Breaths detected in circle circuit while non-circle circuit is selected.	Check circle circuit connections and settings.

Message	Priority	Cause	Action
Check D-Fend	Medium	Water trap not attached.	Check that the water trap is properly attached to the airway module.
Check flow sensors	Medium	System has detected an improper flow pattern in the breathing circuit.	Ensure the internal flow sensors are connected correctly.
Check sample gas out	Medium	Possible blockage in airway module sample gas outlet.	Check for blockage in the airway module sample gas outlet. Remove blockage.
Circuit leak	Medium	Vent TVexp < half of vent TVinsp for at least 30 seconds.	Check for leaks in the patient circuit. Calibrate flow sensors. If problem persists, replace flow sensors.
Circuit leak silenced	Informational	Setting on Alarm Setup menu. Vent TVexp < 50% of vent TVinsp for at least 30 seconds.	Message indicates that the Leak Audio alarm is turned off.
Circuitry >75C shutdown possible	Medium	Power supply temperature exceeds 75C.	Shut down system as soon as possible. Then, check cooling fans and filters.
Cooling fan needs service. System OK.	Medium	Fan reporting error.	Shut down system as soon as possible. Then, check cooling fans and filters.
Cooling fans failed. May overheat.	Medium	Fan reporting error.	Shut down system as soon as possible. Contact a Datex-Ohmeda trained service representative.
Display panel controls failure	Medium	Communication lost between panel and key pad.	Turn the system off and back on.
EtCO2 high	<b>High</b>	EtCO2 > high alarm limit.	Check the patient and EtCO2 settings. Check if absorbent needs to be changed.
EtCO2 low	<b>High</b>	EtCO2 < alarm limit.	Ensure the patient is properly intubated. Check for leaks or blockages in the patient circuit.
EtDES high	Medium	EtDES > alarm limit.	Set the alarm limits appropriately. Decrease the agent concentration.
EtDES low	Informational	EtDES < alarm limit.	Check the fill level on the vaporizers. Set alarm limit appropriately. Increase the agent concentration.
EtENF high	Medium	EtENF > alarm limit.	Set the alarm limits appropriately. Decrease the agent concentration.
EtENF low	Informational	EtENF < alarm limit.	Check the fill level on the vaporizers. Set alarm limit appropriately. Increase the agent concentration.
EtHAL high	Medium	EtHAL > alarm limit.	Set the alarm limits appropriately. Decrease the agent concentration.
EtHAL low	Informational	EtHAL < alarm limit.	Check the fill level on the vaporizers. Set alarm limit appropriately. Increase the agent concentration.
EtISO high	Medium	EtISO > alarm limit.	Set the alarm limits appropriately. Decrease the agent concentration.
EtISO low	Informational	EtISO < alarm limit.	Check the fill level on the vaporizers. Set alarm limit appropriately. Increase the agent concentration.
EtO2 high	Medium	EtO2 > high alarm limit.	Set the alarm limits appropriately. Decrease the O <sub>2</sub> concentration.

Message	Priority	Cause	Action
EtO2 low	Medium	EtO2 < low alarm limit.	Set the alarm limits appropriately. Increase the O <sub>2</sub> concentration.
EtSEV high	Medium	EtSEV > alarm limit.	Set the alarm limits appropriately. Decrease the agent concentration.
EtSEV low	Informational	EtSEV < alarm limit.	Check the fill level on the vaporizers. Set alarm limit appropriately. Increase the agent concentration.
FiCO2 high. Absorbent OK?	<b>High</b>	FiCO2 > alarm limit.	Check if absorbent needs to be changed. Check the patient.
FiDES high	Medium	FiDES > alarm limit.	Set the alarm limits appropriately. Decrease the agent concentration.
FiDES low	Informational	FiDES < alarm limit.	Check the fill level on the vaporizers. Set alarm limit appropriately. Increase the agent concentration.
FiENF high	Medium	FiENF > alarm limit.	Set the alarm limits appropriately. Decrease the agent concentration.
FiENF low	Informational	FiENF < alarm limit.	Check the fill level on the vaporizers. Set alarm limit appropriately. Increase the agent concentration.
FiHAL high	Medium	FiHAL > alarm limit.	Set the alarm limits appropriately. Decrease the agent concentration.
FiHAL low	Informational	FiHAL < alarm limit.	Check the fill level on the vaporizers. Set alarm limit appropriately. Increase the agent concentration.
FiISO high	Medium	FiISO > alarm limit.	Set the alarm limits appropriately. Decrease the agent concentration.
FiISO low	Informational	FiISO < alarm limit.	Check the fill level on the vaporizers. Set alarm limit appropriately. Increase the agent concentration.
FiO2 high	Medium	FiO2 > high alarm limit.	Check the O <sub>2</sub> setting. Recalibrate the O <sub>2</sub> cell and the airway module.
FiO2 low	<b>High</b>	FiO2 < low alarm limit.	Check O <sub>2</sub> setting. Check for leaks or blockages in the patient circuit.
FiSEV high	Medium	FiSEV > alarm limit.	Set the alarm limits appropriately. Decrease the agent concentration.
FiSEV low	Informational	FiSEV < alarm limit.	Check the fill level on the vaporizers. Set alarm limit appropriately. Increase the agent concentration.
For mech vent, set Bag/Vent switch	Informational	Bag/Vent switch is set to Bag and the ACGO switch is set to Circle.	Move switch to the Vent position to start mechanical ventilation.
Gas monitoring not available	Medium	Airway module hardware failure.	Replace airway module.
Increase low MV limit	Medium	Low MVexp limit < half of measured MVexp.	Increase Low MV alarm limit to improve patient disconnection detection.
Insert cassette	Informational	Cassette removal detected during active delivery.	Reinsert the cassette.
Inspiration stopped	Medium	High airway pressure.	Check system for blockages.
Internal failure. System may shut down.	<b>High</b>	Power controller software failure.	Contact a Datex-Ohmeda trained service representative.

Message	Priority	Cause	Action
Internal failure. System may shut down.	Medium	Power controller software failure.	Contact a Datex-Ohmeda trained service representative.
Memory (EEPROM) failure	Informational	Software error.	Contact a Datex-Ohmeda trained service representative.
Module fail. No CO <sub>2</sub> , AA, O <sub>2</sub> data.	Medium	Airway module hardware failure.	Replace module.
Module not compatible	Informational	The monitoring module detected is not compatible with system software.	Remove the incompatible module.
Move Bag/Vent Switch to Bag	Medium	Bag/Vent switch is in the wrong position to start case.	Move switch to the Bag position.
MVexp high	Medium	MVexp > MVexp high alarm limit (for nine breaths or 1 minute).	Change TV, RR, I:E, or PEEP to reset the minute volume to below the MVexp high alarm limit.
MVexp low	Medium	MVexp < MVexp low alarm limit (for nine breaths or 1 minute).	Change TV, RR, I:E, or PEEP to reset the minute volume to above the MVexp low alarm limit.
N <sub>2</sub> O supply pressure low	Informational	N <sub>2</sub> O pipeline pressure is less than 252 kPa (36 psi) and the N <sub>2</sub> O cylinder pressure is less than 2633 kPa (381 psi).	Ensure the N <sub>2</sub> O pipeline and cylinder are properly connected.
Negative airway pressure	<b>High</b>	Paw < -10 cmH <sub>2</sub> O	Check for blockages in the patient circuit.
No battery backup	Medium	Battery or charging failure.	Between cases turn the system circuit breaker off, then back on after 15 seconds to reset the system.
No battery backup for monitor	Informational	Monitor power cable is unplugged.	Make sure the monitor power cable is plugged in.
No exp flow sensor	Medium	Electrical signals show the flow sensor is not connected.	Connect the flow sensor. Replace the flow sensor if necessary.
No fresh gas flow!	<b>High</b>	Possible patient detected while system is in standby.	Disconnect the patient or start a case.
No fresh gas flow?	<b>High</b>	Possible fresh gas flow occlusion or loss of gas pressure.	Switch to circle circuit or bag the patient. Check pipeline supply connection.
No insp flow sensor	Medium	Electrical signals show the flow sensor is not connected.	Connect the flow sensor. Replace the flow sensor if necessary.
O <sub>2</sub> flush stuck on?	Informational	Switch is detected "on" continuously > 30 seconds.	Check flush valve. Ensure flush valve is not sticking.
O <sub>2</sub> monitoring not connected	Medium	O <sub>2</sub> cell not connected.	Install airway gas module or connect the O <sub>2</sub> cell.
O <sub>2</sub> supply pressure low	<b>High</b>	O <sub>2</sub> pipeline pressure is less than 252 kPa (36 psi) and the O <sub>2</sub> cylinder pressure dropped below 2633 kPa (381 psi) for one second.	Ensure the O <sub>2</sub> pipeline and cylinder are properly connected.
PEEP high. Blockage?	<b>High</b>	Paw ≥ sustained limit for 15 seconds. <sup>1</sup>	Check for blockages in the patient circuit.
Plug in power cable. On battery.	Medium	The mains supply is not connected or has failed and the system is using battery power.	Ventilate manually to save power. Make sure the power cable is plugged in and system circuit breaker is on.
Ppeak high	<b>High</b>	Paw > Pmax alarm limit	Check for blockages in the patient circuit.

Message	Priority	Cause	Action
Ppeak low. Leak?	Medium	Peak airway pressure < low P min + 4 cmH <sub>2</sub> O for 20 consecutive seconds if the set respiratory rate is four or higher and 30 seconds if the set respiratory rate is less than four breaths/min.	Check for leaks in the patient circuit.
Replace D-Fend	Medium	Buildup in airway module sample line.	Replace D-Fend.
Replace exp flow sensor	Informational	EEPROM calibration data read failure.	Replace the expiratory flow sensor.
Replace insp flow sensor	Informational	EEPROM calibration data read failure.	Replace the inspiratory flow sensor.
Replace O <sub>2</sub> sensor	Informational	O <sub>2</sub> < 5%.	Calibrate the O <sub>2</sub> cell. Replace the O <sub>2</sub> cell if necessary.
Reverse exp flow. Check valves OK?	Medium	Flow toward the patient seen in the expiratory flow sensor during inspiration for six breaths in a row.	Check the flow sensor condition. Replace the expiratory check valve between cases. Zero flow transducers between cases.
Reverse insp flow. Check valves OK?	Medium	Flow through inspiratory sensor during expiration for six breaths in a row.	Check the flow sensor condition. Replace the expiratory check valve between cases. Zero flow transducers between cases.
RR high	Medium	RR > high alarm limit.	Set the alarm limits appropriately or adjust the RR setting.
RR low	Medium	RR < low alarm limit.	Set the alarm limits appropriately or adjust the RR setting.
Sample line blocked	Medium	Airway module sample line is blocked.	Replace airway module sample line.
Service calibration advised	Informational	Calibration data corrupt.	Contact a Datex-Ohmeda trained service representative.
Set Alt O <sub>2</sub> flow. Check agent setting.	Medium	Multiple possible causes for the failure.	Contact a Datex-Ohmeda trained service representative.
System leak?	Informational	Leak detected between ventilator and patient circuit.	Check for leaks in the breathing system.
System shutdown in <5 min	<b>High</b>	Remaining battery power is between one and five minutes.	Plug in the power cable. Check that the system circuit breaker is on.
Try another cassette. Schedule service.	Medium	Vaporizer detected a cassette temperature failure. Delivery cannot continue.	Insert a different cassette. Contact a Datex-Ohmeda trained service representative.
Turn power Off and On for self tests	Informational	System has been operating for longer than 12 hours without a power up self test.	Turn power off and back on between cases to perform a self test.
Turn switch on to continue use	<b>High</b>	System is in therapy when System switch is turned to Standby.	Turn the System switch to On to continue therapy. The system will return to normal function. If System switch is not turned to On within 8 seconds, the system will shut down.
TV not achieved	Informational	Measured tidal volume is < set tidal volume.	Check for leaks in the patient circuit. Check for leaks in the breathing system.
TVexp high	Medium	TVexp > TVexp high alarm limit (for nine breaths).	Change TV, RR, I:E, or PEEP or change the tidal volume to reset TVexp high alarm limit.

Message	Priority	Cause	Action
TVexp low	Medium	TVexp < TVexp low alarm limit (for nine breaths).	Change TV, RR, I:E, or PEEP or change the tidal volume to reset TVexp low alarm limit.
Unable to drive bellows	Informational	Bellows is collapsed.	Check the drive gas. Increase fresh gas flow (or push the O <sub>2</sub> flush button) to fill the bellows.
Using battery. Power Controller fail.	Medium	Mains supply is OK, but the system is running on the battery.	Shut down the system as soon as possible, and contact a Datex-Ohmeda trained service representative.
Vaporizer failure	Medium	Internal issue with agent delivery hardware.	Change anesthesia method or use an alternate machine. Shutdown the system as soon as possible and contact a Datex-Ohmeda trained service representative.
Ventilate Manually	<b>High</b>	Software or hardware failure prevents mechanical ventilation.	Use a manual bag to ventilate the patient or use an alternate machine. Shutdown the system as soon as possible and contact a Datex-Ohmeda trained service representative.
Ventilate Manually	Medium	Ventilator monitoring only. Software or hardware failure.	Contact a Datex-Ohmeda trained service representative.
Ventilator has no drive gas	<b>High</b>	Drive gas supply is not sufficient to mechanically ventilate.	Check drive gas supply. Ventilate manually until drive gas supply is restored.
Vol and Apnea monitoring off	Informational	Non Circle ACGO selected. Non Circle SCGO selected.	Message will clear when circle circuit is selected.
Vol vent only. No PEEP or PSV.	Medium	Manifold pressure error. Pressure control unavailable. Medium priority alarm when Bag/Vent switch is in Vent and running PCV, PSVPro, or SIMV-PC mode.	Use volume control ventilation mode. Shut down system as soon as possible. Contact a Datex-Ohmeda trained service representative.
Vol vent only. No PEEP or PSV.	Informational	Manifold pressure error. Pressure control unavailable. Informational priority alarm when Bag/Vent switch is in Vent and not running PCV, PSVPro, or SIMV-PC mode; Bag/Vent switch is in Bag; or Non-circle circuit or ACGO is selected.	Continue to use volume control ventilation mode or ventilate manually. Shut down system as soon as possible. Contact a Datex-Ohmeda trained service representative.
Volume sensors disagree	Informational	TVexp > larger of (TVinsp + 3* 'set TV') or (TVinsp + 100)) for six breaths.	Calibrate the flow sensors. Replace the flow sensors if the message does not clear.

<sup>1</sup>The sustained pressure threshold is calculated from the pressure limit setting. The sustained limit is calculated as follows:

Mechanical Ventilation On - Volume Mode	For Plimit ≤ 30 cmH <sub>2</sub> O; the sustained pressure limit is 6 cmH <sub>2</sub> O plus Set PEEP. For Plimit between 30 and 60 cmH <sub>2</sub> O; the sustained pressure limit is 20% of the Plimit plus Set PEEP. For Plimit ≥ 60 cmH <sub>2</sub> O; the sustained pressure limit is 12 cmH <sub>2</sub> O plus Set PEEP.
Mechanical Ventilation On - Pressure Mode	Sustained pressure limit is 50% of set P <sub>insp</sub> or 4 cmH <sub>2</sub> O, whichever is greater, plus Set PEEP and P <sub>max</sub> - P <sub>min</sub> must be less than 50% of set P <sub>insp</sub> or 4 cmH <sub>2</sub> O, whichever is greater.
Mechanical Ventilation Off	For Plimit ≤ 60 cmH <sub>2</sub> O; the sustained pressure limit is 50% of the Plimit setting. For Plimit > 60 cmH <sub>2</sub> O; the sustained pressure limit is 30 cmH <sub>2</sub> O.

## Alarm ranges and default values

Alarm	Range	Increment	Default
EtCO2 High	0.1 – 15, Off % 0.1 – 15, Off kPa 0-115, Off mmHg	0.1% 0.1 kPa 1 mmHg	6.5% 6.5 kPa 50 mmHg
EtCO2 Low	Off, 0.0 – 14.9% Off, 0.0 – 14.9 kPa Off, 0-114 mmHg	0.1% 0.1 kPa 1 mmHg	Off Off Off
EtDES High	0.1 – 20, Off %	0.1%	Off
EtDES Low	Off, 0.1 – 19.9%	0.1%	Off
EtENF High	0.1 – 7.0, Off %	0.1%	Off
EtENF Low	Off, 0.1 – 6.9%	0.1%	Off
EtHAL High	0.1 – 7.0, Off %	0.1%	5%
EtHAL Low	Off, 0.1 – 6.9%	0.1%	Off
EtISO High	0.1 – 7.0, Off %	0.1%	Off
EtISO Low	Off, 0.1 – 6.9%	0.1%	Off
EtO2 High	18-100, Off %	1%	Off
EtO2 Low	Off, 0 – 99%	1%	Off
EtSEV High	0.1 – 10.0, Off %	0.1%	Off
EtSEV Low	Off, 0.1 – 9.9%	0.1%	Off
FiCO2 High	0.1 – 15, Off % 0.1 – 15, Off kPa 0-115, Off mmHg	0.1% 0.1 kPa 1 mmHg	Off Off Off
FiDES High	0.1 – 20.0, Off %	0.1%	15%
FiDES Low	Off, 0.1 – 19.9%	0.1%	Off
FiENF High	0.1-7.0, Off %	0.1%	5%
FiENF Low	Off, 0.1 – 6.9%	0.1%	Off
FiHAL High	0.1 – 7.0, Off %	0.1%	5%
FiHAL Low	Off, 0.1 – 6.9%	0.1%	Off
FiISO High	0.1 – 7.0, Off %	0.1%	5%
FiISO Low	Off, 0.1 – 6.9%	0.1%	Off
FiSEV High	0.1 – 10.0, Off %	0.1%	8%
FiSEV Low	Off, 0.1 – 9.9%	0.1%	Off
High FiO2	19 – 100, Off %	1%	Off
High MV	0.5 – 30.0, Off l/min	0.5 l/min	10 l/min
High RR	2 – 100, Off bpm	1 bpm	Off
High TV	20 – 1600, Off ml	20 ml	1000 ml
Low FiO2	18 – 99%	1%	21%
Low MV	Off, 0.1 – 10.0 l/min	0.1 l/min	2 l/min
Low RR	Off, 0 – 99 bpm	1 bpm	Off
Low TV	Off, 5 – 20 ml 20 – 1500 ml	5 ml 20 ml	Off
Pmax (only high)	12 -100 cmH <sub>2</sub> O, hPa 1.2 – 9.8 kPa 12 – 98 mbar 9 -74 mmHg	1 cmH <sub>2</sub> O, hPa 0.1 kPa 1 mbar 1 mmHg	40 cmH <sub>2</sub> O, hPa 4 kPa 4 mbar 30 mmHg

## Alarm tests

1. Connect a test lung to the patient connection.
2. Start a case.
3. Set the Bag/Vent switch to Vent.
4. Set the O<sub>2</sub> concentration to 30%, and allow the O<sub>2</sub> reading to stabilize.
5. Test the O<sub>2</sub> alarms:
  - Set the **FiO2 low** alarm limit to 50%. Make sure an **FiO2 low** alarm occurs.
  - Set the **FiO2 low** alarm limit back to 21% and make sure that the **FiO2 low** alarm cancels.
  - Set the **FiO2 high** alarm limit to 50%.
  - Push the O<sub>2</sub> flush button.
  - Make sure the **FiO2 high** alarm occurs.
  - Set the **FiO2 high** alarm limit back to 100%. Make sure that the **FiO2 high** alarm cancels.
6. Test the **MVexp low** alarm:
  - Go to the **Alarm Setup** menu.
  - Set the **MV low** alarm limit to greater than the measured minute volume.
  - Make sure that a **MVexp low** alarm occurs.
  - Set the **MV low** alarm limit to off.
7. Test the **Ppeak high** alarm:
  - Set the **Pmax** to less than the peak airway pressure.
  - Make sure that the **Ppeak high** alarm occurs.
  - Set the **Pmax** to the desired level.
8. Test the **PEEP high. Blockage?** alarm:
  - Close the APL valve.
  - Set the Bag/Vent switch to Bag. Mechanical ventilation stops.
  - Block the patient connection and push the O<sub>2</sub> flush button.
  - Make sure that the **PEEP high. Blockage?** alarm occurs after approximately 15 seconds.
9. Test the **Ppeak low. Leak?** alarm:
  - Unblock the patient connection.
  - Set the Bag/Vent switch to Vent.
  - Set the tidal volume and total flow to minimum.
  - Other alarms such as **MVexp low** can occur.
  - Make sure that the **Ppeak low. Leak?** alarm occurs.
10. Set all alarm limits to approved clinical values.



## Breathing system problems

Symptom	Problem	Solution
Gas scavenging flow is too low or too high.	Scavenging extract flow problem.	Use a different scavenging extraction system. Verify flow is within specification.
	Filter blockage. Active systems have a flow indicator.	Replace the filter. Refer to <i>Remove the AGSS receiver filter</i> in the “ <i>Cleaning and Sterilization</i> ” section in Part 2 of the User’s Reference manual.
The bellows fills when the Bag/Vent switch is set to Bag or the bag fills when the switch is set to Vent.	Leak through Bag/Vent switch.	Contact a Datex-Ohmeda trained service representative to repair the system.
The ventilator does not read the position of the Bag/Vent switch.	Ventilator or absorber malfunction.	Ventilate manually. Contact a Datex-Ohmeda trained service representative to repair the system.
APL valve does not operate correctly.	APL valve problem.	Replace APL Valve seal and diaphragm.
Large breathing system leak not quickly located (in Bag mode).	Bag hose not connected properly.	Ensure that the Bag Hose is connected to the Bag port (below the APL Valve).
	Absorber canister not installed correctly.	Reinstall the absorber canister, ensure both pins are engaged.
Bellows falls below top of indicator during “ <i>Bellows assembly test.</i> ”	Leak in the breathing system.	Check, clean, or reposition the pressure relief valve. If the problem persists, replace the pressure relief valve, bellows base, or bellows assembly.

## Electrical problems

**WARNING** If a circuit breaker opens frequently, do not use the system. Have a Datex-Ohmeda trained service representative repair the system.

Symptom	Problem	Solution
Mains indicator is not on.	The electrical power cable is not connected.	Connect the power cable.
	The inlet circuit breaker (switch) is off.	Turn the circuit breaker on.
	The power cable is damaged.	Replace the power cable.
	The electrical socket the power cable connects to has no power.	Use a different electrical socket.
	An internal fuse is open.	Have a Datex-Ohmeda trained service representative repair the system.
One electrical outlet does not have power.	The outlet circuit breaker is off.	Turn the circuit breaker on.
A circuit breaker opens frequently.	Equipment connected to the outlet(s) uses more current than the circuit breaker rating.	Use a different power supply for some of the equipment.
	The equipment connected to the outlet has a short.	Have a Datex-Ohmeda trained service representative repair the system.
The real time clock on the anesthesia display does not maintain accurate time and date.	The battery in the anesthesia display needs to be replaced.	Have a Datex-Ohmeda trained service representative replace the battery.

## Pneumatic problems

Symptom	Problem	Solution
High-pressure leak test fails.	Controls are not set correctly.	Ensure no gas is flowing, turn off the auxiliary flowmeter, and repeat the test.
	Incorrect cylinder connection.	Make sure that there is only one cylinder gasket, the gasket is in good condition, and the connection is tight.
Low-pressure leak with a vaporizer OFF.	Anesthesia machine problem.	Contact a Datex-Ohmeda trained service representative.



## Warranty

This Product is sold by Datex-Ohmeda under the warranties set forth in the following paragraphs. Such warranties are extended only with respect to the purchase of this Product directly from Datex-Ohmeda or Datex-Ohmeda's Authorized Dealers as new merchandise and are extended to the Buyer thereof, other than for the purpose of resale.

For a period of twelve (12) months from the date of original delivery to Buyer or to Buyer's order, but in no event for a period of more than two years from the date of original delivery by Datex-Ohmeda to a Datex-Ohmeda Authorized Dealer, this Product, other than its expendable parts, is warranted against functional defects in materials and workmanship and to conform to the description of the Product contained in this User's Reference manual and accompanying labels and/or inserts, provided that the same is properly operated under the conditions of normal use, that regular periodic maintenance and service is performed and that replacements and repairs are made in accordance with the instructions provided. This same warranty is made for a period of thirty (30) days with respect to expendable parts. The foregoing warranties shall not apply if the Product has been repaired other than by Datex-Ohmeda or in accordance with written instructions provided by Datex-Ohmeda, or altered by anyone other than Datex-Ohmeda, or if the Product has been subject to abuse, misuse, negligence, or accident.

Datex-Ohmeda's sole and exclusive obligation and Buyer's sole and exclusive remedy under the above warranties is limited to repairing or replacing, free of charge, at Datex-Ohmeda's option, a Product, which is telephonically reported to the nearest Datex-Ohmeda Customer Service Center and which, if so advised by Datex-Ohmeda, is thereafter returned with a statement of the observed deficiency, not later than seven (7) days after the expiration date of the applicable warranty, to the Datex-Ohmeda Customer Service and Distribution Center during normal business hours, transportation charges prepaid, and which, upon Datex-Ohmeda's examination, is found not to conform with above warranties. Datex-Ohmeda *shall not be otherwise liable for any damages including but not limited to incidental damages, consequential damages, or special damages.*

*There are no express or implied warranties which extend beyond the warranties hereinabove set forth. Datex-Ohmeda makes no warranty of merchantability or fitness for a particular purpose with respect to the product or parts thereof.*



# Index

## A

- Absorber canister release 2-7
- ACGO 2-8, 2-9
  - positive low-pressure leak test 5-8
  - scavenging the ACGO sample flow 2-10
- Advanced breathing system 2-6
- Agent
  - set measured 3-7
  - settings 2-12
- AGSS 2-4
  - active 3-18
  - connecting active with a flow indicator 3-19
  - connecting active without a flow indicator 3-19
  - passive 3-17
- Air
  - adjust 3-7
  - select 3-3
- Airway modules 6-2
  - agent identification 6-6
  - connection to a patient 6-4
- Aladin cassettes
  - Aladin 7-3
  - Aladin2 7-2
  - cleaning and disinfection 7-5
  - filling Aladin cassettes 7-11
  - filling Aladin2 cassettes 7-6
  - installing cassettes 7-5
- Alarm 8-2
  - priorities 8-2
  - ranges and defaults 8-10
  - silence countdown 2-12
  - silence key 2-11
  - tests 8-11
- Alternate O2 control 2-2, 2-3, 3-14
- Auxiliary O2 2-2, 2-3

## B

- Bag/Vent switch 2-6, 2-7
- Battery 1-5
  - indicator 8-3
  - indicator field 2-12
- Breathing system 2-2, 8-12

## C

- Calibration
  - airway module 6-7
  - flow sensor 5-9
- Cassette 2-2
  - controls 2-5
  - insert 5-9
  - test delivery 5-7
- Checkout 4-1
- Circuit
  - circle 3-7
  - non-circle 2-9, 3-7
- Condenser
  - drain button 2-8
  - operation 3-16
- Controls
  - alternate O2 3-14
  - display 2-11

## D

- Display
  - anesthesia system 2-12
  - controls 2-11

## E

- End case 3-4
- EZchange Canister 2-8
  - operation 3-15

## F

- Filling Aladin cassettes
  - filling with keyed filler system 7-11
  - filling with Quik-Fil system 7-13
  - filling with Saf-T-Fil bottle 7-14
- Filling Aladin2 cassettes
  - filling with Easy-Fil system 7-6
  - filling with Quik-Fil system 7-8
  - filling with Safe-T-Fil bottle 7-9

## G

- Gas setup 3-7

## **L**

Low Pressure Leak test 5-6

## **M**

Manual ventilation 2-7

Menu

- remove 2-17
- select item 2-15
- using 2-14

## **N**

Negative low pressure leak test 5-6

## **O**

O2 flow

- adjust 3-7
- start 3-3

## **P**

Preoperative checkout 4-2

Problems

- breathing system 8-12
- electrical 8-13
- pneumatic 8-13

## **Q**

Quick keys 2-11

- changing gas settings using 3-7
- changing ventilator settings using 3-5

## **S**

Scavenging

- a gas monitor sample flow 2-10
- ACGO sample flow 2-10
- from an auxiliary manual breathing circuit 2-10

Setup

- changing circuit type 3-7
- changing gas and settings 3-7
- gas 3-7
- vent 3-5

Spirometry

- menu functions 3-8
- reference loop 3-10
- split screen 3-9
- viewing loops 3-10

Split screen 3-11

Start case 3-3

Suction 2-3

System

- turn on 3-2

## **T**

Timer 2-11

- start/stop 2-11

Troubleshooting

- breathing system 8-12
- electrical 8-13
- Pneumatic 8-13

## **V**

Vaporizer 7-1

- controls 2-5

Vaporizer cassettes

- filling Aladin cassettes 7-11
- filling Aladin2 cassettes 7-6
- installing cassettes 7-5

Ventilator

- changing modes and settings 3-6
- modes 3-5
- settings 2-12
- setup 3-5

Aisys  
User's Reference Manual  
English  
M1076285  
05 06 18 05 02  
Printed in USA  
©Datex-Ohmeda, Inc.  
All rights reserved

