
S/5 Aespire Anesthesia Machine

Technical Reference Manual



Datex-Ohmeda products have unit serial numbers with coded logic which indicates a product group code, the year of manufacture and a sequential unit number for identification.





Aespire, ProTIVA, SmartVent, and Link-25 are registered trademarks of Datex-Ohmeda Inc.

Other brand names or product names used in this manual are trademarks or registered trademarks of their respective holders.

Covers the following:

S/5 Aespire anesthesia machine

S/5 Aespire 100 anesthesia machine

S/5 ProTIVA anesthesia machine

This document is not to be reproduced in any manner, nor are the contents to be disclosed to anyone, without the express authorization of the product service department, Datex-Ohmeda, Ohmeda Drive, PO Box 7550, Madison, Wisconsin, 53707.

1009-0356-000 07/04 i

^{© 2004} Datex-Ohmeda Inc.

Important

The information contained in this service manual pertains only to those models of products which are marketed by Datex-Ohmeda as of the effective date of this manual or the latest revision thereof. This service manual was prepared for exclusive use by Datex-Ohmeda service personnel in light of their training and experience as well as the availability to them of parts, proper tools and test equipment. Consequently, Datex-Ohmeda provides this service manual to its customers purely as a business convenience and for the customer's general information only without warranty of the results with respect to any application of such information. Furthermore, because of the wide variety of circumstances under which maintenance and repair activities may be performed and the unique nature of each individual's own experience, capacity, and qualifications, the fact that customer has received such information from Datex-Ohmeda does not imply in anyway that Datex-Ohmeda deems said individual to be qualified to perform any such maintenance or repair service. Moreover, it should not be assumed that every acceptable test and safety procedure or method, precaution, tool, equipment or device is referred to within, or that abnormal or unusual circumstances, may not warrant or suggest different or additional procedures or requirements.

This manual is subject to periodic review, update and revision. Customers are cautioned to obtain and consult the latest revision before undertaking any service of the equipment. Comments and suggestions on this manual are invited from our customers. Send your comments and suggestions to the Manager of Technical Communications, Datex-Ohmeda, Ohmeda Drive, PO Box 7550, Madison, Wisconsin 53707.

▲ CAUTION

Servicing of this product in accordance with this service manual should never be undertaken in the absence of proper tools, test equipment and the most recent revision to this service manual which is clearly and thoroughly understood.

Technical Competence

The procedures described in this service manual should be performed by trained and authorized personnel only. Maintenance should only be undertaken by competent individuals who have a general knowledge of and experience with devices of this nature. No repairs should ever be undertaken or attempted by anyone not having such qualifications.

Datex-Ohmeda strongly recommends using only genuine replacement parts, manufactured or sold by Datex-Ohmeda for all repair parts replacements.

Read completely through each step in every procedure before starting the procedure; any exceptions may result in a failure to properly and safely complete the attempted procedure.

ii 07/04 1009-0356-000

Table of Contents

| | Important | ii |
|-----------------------|--|--------|
| | Technical Competence | ii |
| 1 Introduction | | |
| | 1.1 What this manual includes | 1-2 |
| | 1.2 Standard service procedures | 1-2 |
| | 1.2.1 User's Reference Manuals | |
| | 1.2.2 Technical Reference Manuals | |
| | 1.3 What is an S/5 Aespire | |
| | 1.3.1 Aespire 100 | 1-3 |
| | 1.3.2 S/5 ProTIVA | 1-3 |
| | 1.4 Configuration options | 1-4 |
| | 1.4.1 Standard configuration | 1-4 |
| | 1.4.2 Options | 1-4 |
| | 1.5 Components | 1-4 |
| | 1.6 Symbols used in the manual or on the equipment | 1-9 |
| 2 Theory of Operation | | |
| | 2.1 Theory overview | 2-2 |
| | 2.2 Gas flow through the anesthesia machine | 2-2 |
| | 2.2.1 Overview | 2-2 |
| | 2.2.2 Physical connections | 2-6 |
| | 2.2.3 Suction regulators | 2-7 |
| | 2.2.4 System switch | 2-8 |
| | 2.2.5 Flow control | 2-9 |
| | 2.3 Flow through the breathing system | . 2-11 |
| | 2.3.1 Overview of flow paths | . 2-11 |
| | 2.3.2 Manual ventilation. | |
| | 2.3.3 Mechanical ventilation | |
| | 2.3.4 Fresh gas and O_2 flush flow | . 2-18 |

1009-0356-000 07/04 iii

4

3 Checkout Procedure

| | 3.1 Inspect the system | 3-2 |
|-------------------|--|------|
| | 3.2 Pipeline and cylinder tests | 3-3 |
| | 3.3 Flow control, pressure relief, O ₂ supply alarm, and flush flow tests | 3-4 |
| | 3.3.1 With O ₂ monitoring | 3-4 |
| | 3.3.2 Without O ₂ monitoring | 3-6 |
| | 3.3.3 Pressure relief tests | 3-8 |
| | 3.3.4 0 ₂ supply alarm test | 3-8 |
| | 3.3.5 Flush Flow Test | 3-8 |
| | 3.4 Vaporizer back pressure test | 3-10 |
| | 3.5 Low-pressure leak test | 3-11 |
| | 3.5.1 Negative low-pressure leak test | 3-11 |
| | 3.5.2 ISO or BSI standard low-pressure leak test | 3-12 |
| | 3.6 Alarm tests | 3-14 |
| | 3.7 Breathing system tests | 3-16 |
| | 3.8 Auxiliary 02 flowmeter tests | 3-18 |
| | 3.9 Integrated Suction Regulator tests | 3-18 |
| | 3.10 Power failure test | 3-19 |
| | 3.11 Electrical safety tests | 3-19 |
| Repair Procedures | 3 | |
| | 4.1 Servicing the ventilator | 4-3 |
| | 4.2 How to bleed gas pressure from the machine | 4-4 |
| | 4.3 How to remove the rear panels | 4-4 |
| | 4.3.1 To remove the rear upper panel | 4-4 |
| | 4.3.2 To remove the lower access panels | 4-4 |
| | 4.4 How to remove the tabletop | 4-5 |
| | 4.5 Replace pipeline inlet filter | 4-6 |
| | 4.5.1 Replace pipeline inlet check valve | 4-6 |
| | 4.6 Change drive gas | |
| | 4.7 Service the cylinder supply modules | 4-8 |
| | 4.7.1 Tightening procedure for high-pressure tube fittings | 4-8 |
| | 4.7.2 Replace primary regulator module (complete replacement) | |
| | 4.7.3 Replace cylinder inlet filter | |
| | 4.7.4 Replace cylinder check valve | 4-9 |
| | 4.7.5 Replace 3rd-gas cylinder supply module | 4-10 |
| | 4.8 Replace system switch assembly | 4-11 |

| 4.9 Service the flowmeter module | 4-13 |
|--|------|
| 4.9.1 Remove front flowmeter panel shield | 4-13 |
| 4.9.2 Remove flowtubes for cleaning or replacement | 4-13 |
| 4.9.3 Remove complete flowmeter head | 4-15 |
| 4.9.4 Replace flowmeter modules | 4-16 |
| 4.9.5 Replace flowmeter frame | 4-20 |
| 4.9.6 Replace O ₂ supply switch | |
| 4.9.7 Checkout procedure for O ₂ supply switch | 4-21 |
| 4.9.8 Replace secondary regulator manifold or balance regulator manifold $ \ldots $ | 4-22 |
| 4.9.9 Replace 0_2 or $N_2 0$ needle valves (on machines with $N_2 0) \dots$ | 4-23 |
| 4.9.10 Replace an Air needle valve on all machines or an $\rm O_2$ needle valve on mathematical without $\rm N_2O$ | |
| 4.10 Service vaporizer manifold parts | 4-26 |
| 4.10.1 Repair manifold port valve | 4-26 |
| 4.10.2 Checkout procedure for manifold port valve | 4-27 |
| 4.10.3 Replace vaporizer manifold check valve | 4-28 |
| 4.10.4 Replace vaporizer pressure relief valve | 4-30 |
| 4.10.5 Replace vaporizer manifold | 4-31 |
| 4.11 Replace ACGO selector switch | 4-32 |
| 4.12 Clean or replace ACGO port flapper valve | 4-34 |
| 4.13 Reconfigure sample gas return line | 4-35 |
| 4.14 Replace the APL valve | 4-36 |
| 4.15 Replace the bag support arm | 4-37 |
| 4.15.1 Servicing the bag support arm | 4-38 |
| 4.15.2 Replace friction pad in lower bag arm assembly | |
| 4.15.3 Replace bag port housing | 4-40 |
| 4.16 Replace auxiliary O ₂ flowmeter | 4-41 |
| 4.17 Replace the suction control module | 4-42 |
| 4.17.1 Front panel method | 4-42 |
| 4.17.2 Rear panel method | |
| 4.18 Replace ABS breathing system components | |
| 4.18.1 Replace Bag/Vent switch assembly | 4-44 |
| 4.18.2 Replace bellows base latch assembly | |
| 4.19 Replace casters | |
| 4.20 Replace task light and switch | |
| 4.20.1 To replace the task-light switch | |
| 4.20.2 To replace the task-light circuit hoard | 4-47 |

1009-0356-000 07/04 v

| | | 4.21 Replace the display arm or display cables | . 4-48 |
|-----|-----------------|--|--------|
| | | 4.21.1 Cable tie installation | . 4-48 |
| | | 4.21.2 Removing the display arm | . 4-49 |
| | | 4.21.3 Replacing a display cable | . 4-49 |
| | | 4.21.4 Installing the long arm | . 4-50 |
| | | 4.21.5 Installing the short arm | . 4-51 |
| | | 4.22 Replace display and cables in ProTIVA machine | . 4-52 |
| 5 | Maintenance | | |
| | | 5.1 Aespire Planned Maintenance | 5-2 |
| | | 5.1.1 Every twelve (12) months | 5-2 |
| | | 5.1.2 Every twenty-four (24) months | 5-3 |
| | | 5.2 Auxiliary 02 flowmeter tests | 5-4 |
| | | 5.3 Integrated Suction Regulator tests | 5-5 |
| 6 (| Calibration | | |
| | | 6.1 Primary Regulators | 6-2 |
| | | 6.1.1 Test setup | 6-2 |
| | | 6.1.2 Testing Primary Regulators | 6-3 |
| | | 6.1.3 Adjusting Primary Regulators | 6-6 |
| | | 6.2 Secondary Regulators | 6-7 |
| | | 6.2.1 Testing/Adjusting Secondary Regulators or Balance Regulators | 6-7 |
| | | 6.3 Flowmeter Needle Valve Calibration | 6-8 |
| | | 6.3.1 O ₂ Needle Valve Calibration (Minimum Flow) | 6-8 |
| | | 6.3.2 N ₂ O Needle Valve Calibration (Minimum Flow) | . 6-10 |
| | | 6.3.3 Air Needle Valve Calibration (Minimum Flow) | . 6-14 |
| | | 6.3.4 Needle Valve Calibration (Maximum Flow) | . 6-17 |
| | | 6.4 Link system calibration | . 6-18 |
| | | 6.5 O ₂ Flush Regulator | . 6-23 |
| | | 6.6 Airway pressure gauge | . 6-24 |
| | | 6.6.1 Zero the pressure gauge | . 6-24 |
| | | 6.6.2 Checking the pressure gauge accuracy | |
| 7 | Troubleshooting | | |
| | | 7.1 General Troubleshooting | 7-2 |
| | | 7.2 Breathing System Leak Test Guide | 7-4 |
| | | 7.2.1 Breathing system leak test | 7-5 |
| | | 7.2.2 Breathing System Troubleshooting Flowcharts | 7-7 |
| | | 7.2.3 Leak Isolation Tests | . 7-12 |

vi 07/04 1009-0356-000

8 Illustrated Parts

| 8.1 Service tools – Anesthesia machine | 8-3 |
|---|------|
| 8.1.1 Test Devices | 8-3 |
| 8.1.2 Test Tools | 8-4 |
| 8.1.3 Secondary regulator pilot pressure tool | |
| 8.2 External components - front view | 8-6 |
| 8.3 External components - front view references | |
| 8.4 External Components - rear view | 8-8 |
| 8.5 Control module mounting for a ProTIVA machine | 8-9 |
| 8.6 Aespire 100 - exclusive components | 8-10 |
| 8.6.1 AC Inlet (Aespire 100) | 8-12 |
| 8.6.2 Display mount (Aespire 100) | |
| 8.7 Front panel, gauges and system switch | 8-14 |
| 8.8 Rear panel components | 8-15 |
| 8.9 Tabletop components | 8-16 |
| 8.10 Right-side Components | 8-17 |
| 8.11 External components - lower assembly | 8-18 |
| 8.12 Vent Engine Housing | 8-19 |
| 8.13 Display cables, serial board, AGSS flowtube, and sample return | 8-20 |
| 8.14 AC Power cords | 8-21 |
| 8.15 AC Inlet/Outlet Components | 8-22 |
| 8.16 Pipeline inlet fittings | 8-24 |
| 8.17 Cylinder Gas Supplies | 8-25 |
| 8.17.1 Cylinder inlet fittings | 8-26 |
| 8.18 Vaporizer manifold | 8-27 |
| 8.19 Flowmeter components | 8-28 |
| 8.19.1 Flowtube parts | 8-30 |
| 8.19.2 Secondary regulator components | 8-32 |
| 8.20 ABS to machine Interface Components | 8-34 |
| 8.20.1 Flush Regulator, Flush Valve, and ACGO Selector Switch | 8-35 |
| 8 21 Breathing system interface | 8-36 |

1009-0356-000 07/04 vii

| 8.22 Breathing System | 8-37 |
|--|------|
| 8.22.1 APL Valve | 8-37 |
| 8.22.2 Bag/Vent Switch | 8-38 |
| 8.22.3 Absorber canister | 8-39 |
| 8.22.4 Flow Sensor Module | 8-40 |
| 8.22.5 Breathing Circuit Module | 8-41 |
| 8.22.6 Exhalation valve | 8-42 |
| 8.22.7 Bellows | 8-43 |
| 8.22.8 Bellow base | 8-44 |
| 8.22.9 Bag Arms | |
| 8.23 Drawer | 8-46 |
| 8.24 Legris quick-release fittings | 8-47 |
| 8.25 Vent Drive and low-pressure tubing | 8-48 |
| 8.26 Tubing for use with Legris fittings | 8-50 |
| 8.27 Cables and harnesses | 8-52 |
| 8.28 Cables and harnesses (Aespire 100) | 8-54 |
| 8.29 Anesthetic Gas Scavenging System — AGSS | 8-56 |
| 8.29.1 Passive AGSS | 8-56 |
| 8.29.2 Adjustable AGSS | 8-58 |
| 8.29.3 Active AGSS | 8-60 |
| 8.30 Integrated Suction Regulator | 8-62 |
| 8.30.1 Major Components (Continuous and Venturi suction) | 8-62 |
| 8.30.2 Suction Control Module | 8-63 |
| 8.30.3 Venturi assembly | 8-64 |
| 8.31 Auxiliary O ₂ Flowmeter | 8-65 |
| 8.32 Display mounts | 8-66 |
| 8.33 Cable management arm | 8-67 |
| 8.34 Display arm mounting kits for optional equipment | 8-69 |

9 Schematics and Diagrams

viii 07/04 1009-0356-000

In this section

This section provides a general overview of the S/5 Aespire Anesthesia Machine.

| 1.1 | What this manual includes1-2 |
|-----|--|
| 1.2 | 2 Standard service procedures |
| | 1.2.1 User's Reference Manuals |
| | 1.2.2 Technical Reference Manuals1-2 |
| 1.3 | 3 What is an S/5 Aespire |
| | 1.3.1 Aespire 100 |
| | 1.3.2 S/5 ProTIVA1-3 |
| 1.4 | Configuration options1-4 |
| | 1.4.1 Standard configuration |
| | 1.4.2 Options |
| 1.5 | 6 Components |
| 16 | S Symbols used in the manual or on the equipment |

1009-0356-000 07/04 1-1

1.1 What this manual includes

Anesthesia Machine

This manual covers the service information for the S/5 Aespire line of anesthesia machines. It covers the following components:

- · gas delivery components,
- · breathing system components,
- frame component (except those strictly associated with a specific ventilator),
- optional suction regulator and auxiliary O₂ flowmeter.

Ventilator

The ventilator associated with the S/5 Aespire machine has its own Technical Reference Manual:

for the 7100 Ventilator see manual 1006-0836-000.

S/5 ProTIVA

The ProTIVA machine is configured with standard Aespire machine components, with the exception of the vaporizer manifold and the Ventilator/Monitoring Display (7100 Control Module) mounting solution (refer to Section 8.5).

B Braun equipment is not covered in this manual. Refer to B Braun service documentation.

Other equipment

Other equipment may be attached to the system on the display mount, the top shelf, or on the side dovetail rails. Consult separate documentation relative to these items for details.

1.2 Standard service procedures

1.2.1 User's Reference Manuals

Some sections of this manual refer you to the User's Reference Manual for the S/5 Aespire. To expedite repairs, you must have, and be familiar with, the User's Reference Manuals for this product.

Refer to the S/5 Aespire User's Reference Manual if you need further information about the operation of the system.

1.2.2 Technical Reference Manuals

You must first determine where a problem is located before you can determine which Technical Reference Manual to use:

- Use this manual for machine and breathing system related issues.
- Use the 7100 Ventilator manual for ventilator related issues.

1-2 07/04 1009-0356-000

1.3 What is an S/5 Aespire

The S/5 Aespire is a compact, integrated and intuitive anesthesia delivery system. The ventilator portion provides mechanical ventilation to a patient during surgery as well as monitoring and displaying various patient parameters.

The system uses a microprocessor-controlled ventilator with internal monitors, electronic PEEP, Volume Mode, and other optional features. A serial interface permits communication to cardiovascular and respiratory gas monitoring.

Note

Configurations available for this product depend on local market and standards requirements. Illustrations in this manual may not represent all configurations of the product.

The S/5 Aespire is not suitable for use in an MRI environment.

1.3.1 Aespire 100

The Aespire 100 machine is based on the standard S/5 Aespire machine with the following exceptions to available features or options:

- does not include the RS232 Serial Interface
- · does not include the Bi-level LED light strip
- the 2 Vap manifold is standard (does not support 1 Vap manifold)
- not available with AC power outlets (area used by AC Inlet)
- uses 4-inch casters instead of 5-inch casters

1.3.2 S/5 ProTIVA

The S/5 ProTIVA is a special adaptation of the S/5 Aespire machine for use with B Braun intravenous drug delivery components.

1009-0356-000 07/04 1-3

1.4 Configuration options

1.4.1 Standard configuration

The standard configuration includes the following items. Items marked with an asterisk (*) are not included in the Aespire 100 machine.

- 7100 Ventilator
- Advanced Breathing System (ABS)
- Auxiliary Common Gas Outlet (ACGO)
- Serial Interface RS232*
- Bi-level LED light strip*
- Two large drawers

1.4.2 Options

Options include the following items.

Items marked with an asterisk (*) are not available in the Aespire 100 machine).

- · selected software features
- vaporizer manifold (1 Vap* or 2 Vap)
- pipeline configurations (0₂/N₂0, 0₂/Air, or 0₂/N₂0/Air)
- · gas cylinder configurations (two inboard, one outboard)
 - inboard configuration = 0₂/N₂0, 0₂/Air, or 0₂/0₂
 - outboard configuration = N₂O only
- manual bag (on support arm or on tube)
- gas scavenging (active, adjustable, passive, or venturi)
- a suction regulator (pipeline vacuum or venturi vacuum)
- an auxiliary O₂ flowmeter
- localized electrical power outlets* (isolated or non-isolated)
- · various display mounting solutions

1.5 Components

The following figures show the front and rear views of the machine. There are some differences between models.


Figure 1-1 • S/5 Aespire (front view - left side)

Figure 1-2 • S/5 Aespire (front view - right side)

Figure 1-3 • S/5 Aespire (rear view)

Figure 1-4 • S/5 ProTIVA with a typical B Braun fluid manager (fm) system


1-4 07/04 1009-0356-000



- 1. Auxiliary common gas outlet (ACGO) switch
- 2. ACGO
- 3. Inspiratory check valve
- 4. Inspiratory flow sensor or flow port adapter
- 5. Canister (carbon dioxide absorbent)
- 6. Canister release
- 7. Expiratory flow sensor or flow port adapter
- 8. Leak test plug
- 9. Expiratory check valve
- 10. Breathing system release
- 11. Manual bag port
- 12. APL (adjustable pressure-limiting) valve
- 13. Bag/Vent switch
- 14. Bellows assembly
- 15. Pressure gauge (airway)
- 16. Sample gas return port
- 17. Scavenging flow indicator
- 18. RS-232 Serial port (not available in the Aespire 100 machine)

Figure 1-1 • S/5 Aespire (front view - left side)


1009-0356-000 07/04 1-5



- 1. ABS (Advanced Breathing System)
- 2. Flow controls
- 3. Ventilator Display/Control Module
- 4. Dovetail rails
- 5. Vaporizer
- 6. Pipeline pressure gauge(s) (upper row)
- 7. System switch
- 8. Cylinder pressure gauge(s) (lower row)
- 9. 0₂ Flush
- 10. Auxiliary 0₂ flowmeter
- 11. Suction regulator

Figure 1-2 • S/5 Aespire (front view - right side)

1-6 07/04 1009-0356-000




- 1. Circuit Breaker for Electrical Outlets*
- 2. Electrical Outlets*
- 3. Pipeline Connection(s)
- 4. Cylinder Supplies
- 5. System Circuit Breaker (AC Inlet fuses for Aespire 100 machine)
- 6. Mains Inlet
- 7. Equipotential Stud

Figure 1-3 • S/5 Aespire (rear view)

1009-0356-000 07/04 1-7

^{*} Items marked with an asterisk (*) are not available in the Aespire 100 machine.



- 1. Ventilator/Monitoring display
- 2. Syringe pumps (3)
- 3. fm controller
- 4. fm computer
- 5. fm segment

Figure 1-4 • S/5 ProTIVA with a typical B Braun fluid manager (fm) system

1-8 07/04 1009-0356-000

1.6 Symbols in t ini il or on the equi ment

Warnings and Cautions tell you about dangerous conditions that can occur if you do not follow all instructions in this manual.

Warnings tell about a condition that can cause injury to the operator or the patient.

Cautions tell about a condition that can cause damage to the equipment. Read and follow all warnings and cautions.

Other symbols replace words on the equipment or in Datex-Ohmeda manuals. No one device or manual uses all of the symbols. These symbols include:

| Ţ | On (power) | × | Alarm silence button |
|------------|--|---------------------------------|---|
| 0 | Off (power) | $[\mathcal{A}]$ | Alarm silence touch key (Tec 6). |
| (J | Standby | <u>*</u> | Type B equipment |
| Ċ | Standby or preparatory state for part of the equipment | ፟ | Type BF equipment |
| \odot | "ON" only for part of the equipment | • | Type CF equipment |
| Ċ | "OFF" only for part of the equipment | \triangle | Caution, ISO 7000-0434 |
| === | Direct current | $\mathbf{\Lambda}$ | Attention, refer to product instructions, IEC 601-1 |
| \sim | Alternating current | 4 | Dangerous voltage |
| (| Protective earth ground | (| Electrical input |
| Ť | Earth ground | \Longrightarrow | Electrical output |
| 7 | Frame or chassis ground | \leftarrow | Pneumatic input |
| \Diamond | Equipotential | $\qquad \longrightarrow \qquad$ | Pneumatic output |

1009-0356-000 07/04 1-9

| + | Plus, positive polarity | \longrightarrow | Movement in one direction |
|-------------------------|----------------------------------|-----------------------|----------------------------------|
| = | Minus, negative polarity | \longleftrightarrow | Movement in two directions |
| | Variability | T | Read top of float |
| _ 00 | Variability in steps | ○ | Vacuum inlet |
| <u>11</u> | This way up | | Suction bottle outlet |
| - <u>Ö</u> - | Lamp, lighting, illumination | | Cylinder |
| Ī | Lock | 318 | Isolation transformer |
| ī | Unlock | XX | Linkage system |
| Ŧ | Close drain | | Risk of Explosion. |
| ت | Open drain (remove liquid) | | Low pressure leak test |
| 134°C | Autoclavable | } | Mechanical ventilation |
| 13v°C | Not autoclavable | APL | Bag position/ manual ventilation |
| Insp | Inspiratory flow | Exp | Expiratory flow |
| 0 ₂ % | O ₂ sensor connection | 02+ | O ₂ Flush button |
| REF | Stock Number | SN | Serial Number |

1-10 07/04 1009-0356-000



Alarm silence touch key



Volume alarms On/Off touch key



End case touch key



Menu touch key



Circle breathing circuit module



Bain/Mapleson D breathing circuit module



The primary regulator is set to pressure less than 345 kPa (50 psi)



The primary regulator is set to pressure less than 414 kPa (60 psi)



Absorber on



CO₂ Bypass Option



Absorber off (CO₂ Bypass active)



Systems with this mark agree with the European Council Directive (93/42/EEC) for Medical Devices when they are used as specified in their Operation and Maintenance Manuals. The xxxx is the certification number of the Notified Body used by Datex-Ohmeda's Quality Systems.



€ European Union Representative

1009-0356-000 07/04 1-11

Notes

1-12 07/04 1009-0356-000

2 Theory of Operation

| n this section | 2.1 Theory overview | 2-2 |
|----------------|---|------|
| | 2.2 Gas flow through the anesthesia machine | 2-2 |
| | 2.2.1 Overview | 2-2 |
| | 2.2.2 Physical connections | 2-6 |
| | 2.2.3 Venturi Suction | 2-7 |
| | 2.2.4 System switch | 2-8 |
| | 2.2.5 Flow control | 2-9 |
| | 2.3 Flow through the breathing system | 2-11 |
| | 2.3.1 Overview of flow paths | 2-11 |
| | 2.3.2 Manual ventilation | 2-12 |
| | 2.3.3 Mechanical ventilation | 2-15 |
| | 2.3.4 Fresh gas and O _o flush flow | 2-18 |

1009-0356-000 07/04 2-1

2.1 Theory overview

This section describes:

- The flow of gas through the anesthesia machine.
- The flow of gas through the breathing system.
- Electrical signals between the anesthesia machine, including the breathing system, and the ventilator.

2.2 Gas flow through the anesthesia machine

2.2.1 Overview

Refer to Figure 2-1.

Gas supplies

Gas comes into the system through a pipeline (2) or cylinder (4) connection. All connections have indexed fittings, filters, and check valves (one-way valves). Gauges show the pipeline (1) and cylinder (3) pressures.

A primary regulator (5) decreases the cylinder pressures to approximately pipeline levels. A pressure relief valve (6) helps protect the system from high pressures.

To help prevent problems with the gas supplies:

- · Install yoke plugs on all empty cylinder connections.
- When a pipeline supply is adequate, keep the cylinder valve closed.

02 flow

Pipeline or regulated cylinder pressure supplies O_2 directly to the ventilator (**7a** for O_2 drive gas) and the venturi suction (**21a** for O_2 drive gas) supply connection. An additional regulator (**13**) decreases the pressure for the flush valve (**14a**) and the auxiliary flowmeter (**25**).

The flush valve supplies high flows of 0_2 to the fresh gas outlet (**26** or **27**) when you push the flush button. The flush pressure switch (**14b**) monitors activation of the flush valve.

When the system switch (8) is On, O_2 flows to the rest of the system.

A secondary regulator (**10**) supplies a constant O_2 pressure to the O_2 flow control valve (**11**). There is a minimum flow of 25 to 75 mL/min (for dual-tube flowmeters) or 175 to 225 mL/min (for single-tube flowmeters) through the O_2 flowmeter (**12**).

The O_2 pressure switch (9) monitors the O_2 supply pressure. If the pressure is too low, an alarm appears on the ventilator display.

2-2 07/04 1009-0356-000

Air and N₂O flow

Pipeline or regulated cylinder pressure supplies Air directly to the ventilator (**7b** for Air drive gas) and the venturi suction (**21b** for Air drive gas) supply connection.

When the system switch (8) is On, air flows to the rest of the system.

A secondary regulator (18) supplies the Air flow control valve (19). Because there is no balance regulator, air flow continues at the set rate during an O₂ supply failure.

A balance regulator (**15**) controls the N $_2$ O supply pressure to the N $_2$ O flow control valve(**16**). The O $_2$ secondary regulator pressure at a pilot port controls the output of the balance regulator. The N $_2$ O output pressure drops with decreasing O $_2$ supply pressure and shuts off hypoxic gas flow before the O $_2$ supply pressure reaches zero.

A chain link system (Link-25) on the N_2O and O_2 flow controls (**16, 11**) helps keep the O_2 concentration higher than 21% (approximate value) at the common gas outlet.

Mixed gas

The mixed gas goes from the flowmeter outlet, through the vaporizer manifold and vaporizer (23) that is On, to the ACGO selector switch (E). A pressure relief valve (24) sets the maximum outlet pressure.

The ACGO selector switch directs the mixed gas to the selected circuit — to the breathing system (26) or to the ACGO (27).

1009-0356-000 07/04 2-3

Key to Numbered Components

- 1. Pipeline pressure gauge
- 2. Pipeline inlet
- 3. Cylinder pressure gauge
- 4. Cylinder inlet (maximum of 3 cylinders)
- 5. Primary regulator (cylinder pressure)
- 6. High-pressure relief valve (758 kPa / 110 psi)*
- 7. Supply connections for the ventilator
 - a. O2 drive gas
 - b. Air drive gas
- 8. System switch
- 9. Switch for low 02 supply pressure alarm (used with the ventilator)
- 10. O₂ secondary regulator (207 kPa / 30 psi)*
- 11. 02 flow control valve
- 12. 02 flow tube(s)
- 13. O_2 flush and auxiliary flowmeter regulator (241 kPa / 35 psi)*
- 14. 0₂ Flush
 - a. Flush valve
 - b. Pressure switch (used with the ventilator)
- 15. N₂O balance regulator
- 16. N₂O flow control valve
- 17. N₂O flow tube(s)
- 18. Air secondary regulator (207 kPa / 30 psi)*
- 19. Air flow control valve
- 20. Air flow tube(s)
- 21. Supply connection for Venturi suction
 - a. O2 drive gas
 - b. Air drive gas
- 22. Vaporizer port valve
- 23. Vaporizer
- 24. Low-pressure relief valve (38 kPa / 5.5 psi)*
- 25. Auxiliary flowmeter (optional)
- 26. To ABS
- 27. To ACGO
- 28. Test port (primary regulator)
- 29. Test port (secondary/balance regulator)
- * Approximate values

Key to Symbols

K > + <

Pneumatic Connection

 \Diamond

Filter

D

Direction of Flow

0

Check Valve

2-4 07/04 1009-0356-000




Figure 2-1 • Pneumatic circuit

1009-0356-000 07/04 2-5

2.2.2 Physical Figure 2-2 shows the physical path that the gas takes. **connections**




Figure 2-2 • Typical tubing connections - pictorial

2-6 07/04 1009-0356-000

2.2.3 Suction regulators Pipeline vacuum

The suction regulator (shown in Figure 2-2) uses an external vacuum source.

Venturi Drive vacuum

The suction regulator (shown in Figure 2-3) uses an internal, venturi derived vacuum source.

Drive gas (internally plumbed **Air** or **0**₂) enters the Venturi Module (**VM**) at the drive port (**A**). As the drive gas passes through the venturi module, a vacuum is created at port **B**. The drive gas exits the venturi module at port **C** and is exhausted outside the machine through the muffler (**D**).

The control port **(E)** on the venturi module responds to pneumatic signals from the front panel switch on the Suction Control Module **(SCM)** to turn the venturi vacuum drive gas on or off. The check valve **(CV)** helps prevent pressurization of the suction circuitry if the exhaust is occluded or the venturi unit fails.




Figure 2-3 • Venturi suction

1009-0356-000 07/04 2-7

2.2.4 System switch

The system switch has two positions: On and Standby.

In the Standby position


The switch:


- Turns off the ventilator (electrical).
- Stops 0₂ and Air to the flowhead (pneumatic).
- Without O₂ pressure, the N₂O balance regulator stops N₂O.

In the On position

The switch:

- Turns on the ventilator (electrical).
- Supplies O₂ and Air to the flowhead.
- \bullet With adequate ${\rm O}_2$ pressure, the ${\rm N}_2{\rm O}$ balance regulator supplies ${\rm N}_2{\rm O}$.





2-8 07/04 1009-0356-000

2.2.5 Flow control


Needle valves (one for each gas) adjust gas flows. Clockwise rotation decreases flow. Counterclockwise increases flow. Mechanical stops set minimum flows for all gases. The link system sets the maximum ratio of $\rm N_2O$ to $\rm O_2$.

▲ WARNING

The Link 25 Proportioning System sets a minimum $\rm O_2$ concentration in the fresh gas stream when only $\rm O_2$ and $\rm N_2O$ are used. Use of an absorber or another gas can still cause a hypoxic mixture to be delivered to the patient, especially at low $\rm O_2$ flow rates.

Minimum flows


At minimum flow, two tabs prevent clockwise rotation of the valve stem. One tab is on the stop collar; the other is on the valve body.



Link system

The chain link system helps assure an approximate minimum 1 to 3 ratio of flow between $\rm O_2$ and $\rm N_2O$. When engaged (minimum $\rm O_2$ concentration), a tab on the $\rm O_2$ knob is in contact with a tab on the $\rm O_2$ sprocket so that the $\rm O_2$ and $\rm N_2O$ knobs turn together:

- an increase in N₂O flow causes an increase in O₂ flow,
- a decrease in O₂ flow causes a decrease in N₂O flow.



O₂ Knob

1009-0356-000 07/04 2-9

Higher concentrations of O_2 are possible when the link system is not engaged: either by reducing the N_2O flow below the point of engagement or by increasing O_2 flow above the point of engagement.

When the N $_2$ O flow is below the point of engagement, increasing the N $_2$ O flow turns the O $_2$ sprocket without changing the O $_2$ flow. At the point of engagement, the tab on the O $_2$ sprocket makes contact with the tab on the O $_2$ knob. Once the linkage is engaged, turning the N $_2$ O flow control counterclockwise (increase in N $_2$ O flow) also turns the O $_2$ knob counterclockwise (increase in O $_2$ flow) to maintain a nominal 25% minimum O $_2$ concentration.


Decreasing the N_2O flow from the engagement point rotates the tab on the O_2 sprocket away from the tab on the O_2 knob. Increasing the O_2 flow rotates the knob tab away from the sprocket tab. Either action increases the O_2 concentration above 21%. Sufficiently decreasing O_2 flow or increasing the N_2O flow brings the two tabs back into contact and engages the linkage.

The kick-in point is defined as the $\rm N_2O$ flow at which the $\rm N_2O$ valve becomes engaged with the $\rm O_2$ valve flowing at 200 mL/min. This engagement point is an arbitrary benchmark that assists in calibrating the proportioning system. The position of the kick-in is set in the factory. During field calibration, you set the $\rm O_2$ flow to 200 mL/min and the $\rm N_2O$ flow to the kick-in flow (usually in the range of 400 to 700 mL/min) and then install the sprockets with the $\rm O_2$ knob/sprocket engaged.

Maximum flows

All gas flows in Canada require maximum flow stops. A maximum stop collar on the body of the needle valve and a stop collar on the stem of the needle valve set the maximum flow.

At maximum flow, a tab on the stop collar hits the tab on the maximum stop collar and prevents you from turning the knob further counterclockwise. As you decrease the flow, the valve stem moves toward the needle valve assembly and clears the tab.



2-10 07/04 1009-0356-000

2.3 Flow through the breathing system

2.3.1 Overview of flow paths

This section looks at three types of flow paths.

- Ventilation paths: How gas flows from the drive source (bag or bellows) to and from the patient.
- Fresh gas paths: Fresh gas can flow from the machine interface directly to the patient through the inspiratory check valve, or through the absorber into the expiratory flow, or directly to an external circuit through the optional auxiliary common gas outlet.
- Scavenged gas paths: APL or Pop-off.

1009-0356-000 07/04 2-11

2.3.2 Manual ventilation


Manual inspiration

(Figure 2-4)

The Bag/Vent switch closes the ventilator path (B).

Gas flows from the bag (1), through the absorber (2), into the breathing circuit module, and through a unidirectional valve (inspiratory check valve) to the patient (3).

During inspiration, fresh gas (**FG**) flows from the machine into the inspiratory limb, upstream of the inspiratory check valve.



- B Bag/Vent switch to Bag
- FG Fresh Gas
- AP Airway Pressure
- 1 Flow to absorber
- 2 Flow from absorber
- 3 Inspiratory flow

Figure 2-4 • Gas flow during manual inspiration

2-12 07/04 1009-0356-000

Manual expiration (Figure 2-5)

The Bag/Vent switch keeps the ventilator path closed (B).

Gas flows from the patient (4), through a unidirectional valve (expiratory check valve), and into the bag (5).

During exhalation, fresh gas flows backwards through the absorber (**FG**) into the expiratory limb, downstream of the expiratory check valve.

For machines that are plumbed to return sample gas to the breathing system, the returned gas (**SGR**) enters the breathing system after the expiratory check valve (Refer to section 4.13).





Figure 2-5 • Flow during manual expiration

1009-0356-000 07/04 2-13

APL Valve (Figure 2-6)

The APL valve sets a pressure limit for manual ventilation.

As you turn the APL knob, it puts more or less force on the APL disc and seat (**D/S**). If the circuit pressure is too high (**6**), the disc and seat inside the diaphragm opens and vents gas to the scavenging system (**7**).



D/S APL disc and seat

6 APL flow

7 To scavenging

Figure 2-6 • Flow through the APL Valve

2-14 07/04 1009-0356-000

2.3.3 Mechanical ventilation


Mechanical inspiration

(Figure 2-7)

The Bag/Vent switch closes the manual path (\mathbf{V}) . Pilot pressure (\mathbf{P}) closes the exhalation valve.

Drive gas (**D**) pushes down on the bellows. Gas flows from the bellows (**1**), through the absorber (**2**), and through a unidirectional valve (inspiratory check valve) to the patient (**3**).

During inspiration, fresh gas flows into the inspiratory limb, upstream of the inspiratory check valve.



- V Bag/Vent switch to Vent
- P Pilot pressure
- **D** Drive gas
- FG Fresh Gas
- **AP** Airway Pressure
- 1 Flow to absorber
- 2 Flow from absorber
- 3 Inspiratory flow

Figure 2-7 • Mechanical inspiration

Mechanical expiration


(Figure 2-8)

Drive-gas flow stops and the exhalation valve opens. Exhaled gas flows from the patient (4), through a unidirectional valve (expiratory check valve) and into the bellows (5). Residual drive gas (**D**) flows out of the bellows to the scavenging system (6).

If PEEP is selected, static pressure on the pilot port of the exhalation valve sets the PEEP level.

During exhalation, fresh gas flows backwards through the absorber (**FG**) into the expiratory limb, downstream of the expiratory check valve.

For machines that are plumbed to return sample gas to the breathing system, the returned gas (**SGR**) enters the breathing system after the expiratory check valve (Refer to section 4.13).



AP Airway Pressure

D Drive gas

FG Fresh Gas

SGR Sample Gas Return

- 4 Expiratory flow
- 5 Flow to bellows
- 6 To scavenging


Figure 2-8 • Flow through the APL Valve

2-16 07/04 1009-0356-000

Pop-off valve(Figure 2-9)

The pop-off valve limits the pressure inside the bellows to $2.5~{\rm cm}~{\rm H}_2{\rm O}$ above the drive gas pressure. This normally occurs when the bellows reaches the top of the housing at the end of exhalation.

Excess gas (7) vents to the scavenging system (6) through the pop-off valve and the exhalation valve.



- 6 Pop-off flow
- 7 To scavenging

Figure 2-9 • Flow through the pop-off valve

2.3.4 Fresh gas and $\mathbf{0_2}$ flush flow

To ABS breathing system (Figure 2-10)

Fresh gas (1) flows from the vaporizer manifold outlet to the ACGO Selector Switch.

With the ACGO Selector Switch in the ABS position, fresh gas flow is channeled to the breathing system.

The output of the $\rm O_2$ Flush regulator (2) is channeled to the $\rm O_2$ Flush valve. When activated, $\rm O_2$ flush flow joins the fresh gas flow in the ACGO Selector Switch.




Figure 2-10 • Fresh gas and O₂ flush flow (to ABS)

2-18

Auxiliary Common Gas Outlet(Figure 2-11)

Fresh gas (1) flows from the vaporizer manifold outlet to the ACGO Selector Switch.

With the ACGO Selector Switch in the ACGO position, fresh gas flow is channeled to the ACGO outlet.

At the ACGO outlet, a small sample is diverted to the $\rm O_2\,Sensor\,in$ the ABS for $\rm O_2\,$ monitoring.

The output of the $\rm O_2$ Flush regulator (2) is channeled to the $\rm O_2$ Flush valve. When activated, $\rm O_2$ flush flow joins the fresh gas flow in the ACGO Selector Switch.




Figure 2-11 • Fresh gas and O2 flush flow (to ACGO)

Notes

2-20 07/04 1009-0356-000

3 Checkout Procedure

| In this section | 3.1 Inspect the system | |
|-----------------|--|------|
| | 3.2 Pipeline and cylinder tests | |
| | 3.3 Flow control, pressure relief, O ₂ supply alarm, and flush flow tests | 3-4 |
| | 3.3.1 With O ₂ monitoring | |
| | 3.3.2 Without O ₂ monitoring | |
| | 3.3.3 Pressure relief tests | |
| | 3.3.4 O ₂ supply alarm test | |
| | 3.3.5 Flush Flow Test | |
| | 3.4 Vaporizer back pressure test | 3-10 |
| | 3.5 Low-pressure leak test | 3-11 |
| | 3.5.1 Negative low-pressure leak test | 3-11 |
| | 3.5.2 ISO or BSI standard low-pressure leak test | |
| | 3.6 Alarm tests | 3-14 |
| | 3.7 Breathing system tests | 3-16 |
| | 3.8 Auxiliary O ₂ flowmeter tests | 3-18 |
| | 3.9 Integrated Suction Regulator tests | 3-18 |
| | 3.10 Power failure test | 3-19 |
| | 3 11 Flactrical safety tests | 3_10 |

▲ WARNINGS

After any repair or service of the Aespire system, complete all tests in this section.

Before you do the tests in this section:

- Complete all necessary calibrations and subassembly tests. Refer to the individual procedures for a list of necessary calibrations.
- Completely reassemble the system.

If a test failure occurs, make appropriate repairs and test for correct operation.

3.1 Inspect the system

Make sure that:

- The equipment is not damaged.
- · All components are correctly attached.
- Pipeline gas supplies are connected and the pressures are correct.
- A supply of reserve O₂ is provided and connected to the machine.
- · Cylinder valves are closed on models with cylinder supplies.

⚠ CAUTION

The upper shelf weight limit is 34 kg (75 lb).

⚠ WARNING

Do not leave gas cylinder valves open if the pipeline supply is in use. Cylinder supplies could be depleted, leaving an insufficient reserve supply in case of pipeline failure.

- Models with cylinder supplies have a cylinder wrench attached to the system.
- Make sure the casters are not loose and the brake is set and prevents movement.



3-2 07/04 1009-0356-000

3.2 Pipeline and cylinder tests

▲ CAUTION To prevent damage:

- Open the cylinder valves slowly.
- Do not force the flow controls.

If your system does not use cylinder supplies, do not do steps 2 and 3.

- Disconnect the pipeline supplies and close all cylinder valves (if equipped). If the pipeline and the cylinder pressure gauges are not at zero, bleed all gasses from the system.
 - a. Connect an O2 supply.
 - b. Set the system switch to On.
 - c. Set the flow controls to mid range.
 - d. Make sure that all gauges but O2 are at zero.
 - e. Disconnect the O2 supply.
 - f. Make sure that the $\rm O_2$ gauge goes to zero. As pressure decreases, alarms for low $\rm O_2$ supply pressure should occur.
- 2. Make sure that the cylinders are full:
 - a. Open each cylinder valve.
 - Make sure that each cylinder has sufficient pressure. If not, close the applicable cylinder valve and install a full cylinder.
- 3. Test one cylinder at a time for high pressure leaks:
 - a. Set the system switch to Standby, which stops the O2 flow.
 - b. If equipped, turn the auxiliary O₂ flow control fully clockwise (no flow).
 - c. If equipped, turn off venturi derived suction.
 - d. Open the cylinder.
 - e. Record the cylinder pressure.
 - f. Close the cylinder valve.
 - g. Record the cylinder pressure after one minute. If the pressure decreases more than indicated below, there is a leak.

5000 kPa (725 psig) for ventilator drive gas. **690 kPa (100 psig)** for non ventilator drive gas.

Install a new cylinder gasket and do this step again.

h. Repeat step 3 for all cylinders.

⚠ WARNING

Do not leave gas cylinder valves open if the pipeline supply is in use. Cylinder supplies could be depleted, leaving an insufficient reserve supply in case of pipeline failure.

4. Connect the pipeline supplies one at a time and ensure that the corresponding gauge indicates pipeline pressure.

3.3 Flow control, pressure relief, 02 supply alarm, and flush flow tests

If the system includes $\mathbf{0}_2$ monitoring, complete the flow control tests in Section 3.3.1, "With $\mathbf{0}_2$ monitoring".

If the system does not include $\mathbf{0}_2$ monitoring, complete the flow control tests in Section 3.3.2, "Without $\mathbf{0}_2$ monitoring".

3.3.1 With 02 monitoring

⚠ WARNING

Nitrous oxide (N₂O) flows through the system during this test. Use a safe and approved procedure to collect and remove it.

- 1. Set up the gas scavenging system.
 - a. Connect the AGSS to a gas scavenging system.
 - b. Attach a patient circuit and plug the patient port.
 - c. Attach a bag to the bag port (or plug the bag port).
 - d. Set the Bag/Vent switch to Bag.
 - e. Adjust the APL valve to minimum.
- 2. Connect the pipeline supplies or slowly open the cylinder valves.
- 3. Turn all flow controls fully clockwise (minimum flow).
- 4. Set the ACGO selector switch to ABS.
- 5. Turn on the system.
- Confirm that the O₂ sensor measures 21% in room air and 100% in pure
 O₂. If not, calibrate the O₂ sensor.
- 7. Make sure that:
 - For a dual-tube O₂ flowmeter, the O₂ flowtube shows 0.025 to 0.075 L/min.
 - For a single-tube O₂ flowmeter, the O₂ flowtube shows 0.175 to 0.225 L/min.
 - . The other flowtubes show no gas flow.
- 8. Set the flow controls to mid range of each flowtube and make sure that the flowtube floats rotate and move smoothly.

Note

If the system does not include N₂O, skip steps 9 and 10.

- Check the Link proportioning system concentration (increasing N₂O flow).Observe the following precautions:
 - Start with all valves at the minimum setting.
 - b. Adjust only the N₂O flow control.
 - c. Increase the $\rm N_2O$ flow as specified in the following table and make sure the $\rm O_2$ concentration is in range.

3-4 07/04 1009-0356-000

Note: Allow the $\rm O_2$ monitor to stabilize. At the lower flows, the $\rm O_2$ monitor may take up to 90 seconds to stabilize.

d. If you overshoot a setting, turn the $\rm O_2$ flow control clockwise until the $\rm N_2O$ flow decreases to the previous setting before continuing the test.

| Set the N ₂ O flow (L/min) | Measured 0 ₂ |
|---------------------------------------|-------------------------|
| 0.15 (dual flowtubes only) | 21% minimum |
| 0.5 (dual flowtubes only) | 21% minimum |
| 0.8 | 21% to 30% |
| 1.0 | 21% to 30% |
| 2.0 | 21% to 30% |
| 6.0 | 21% to 30% |
| 9.0 | 21% to 30% |

- 10. Check the proportioning system concentration (decreasing $\rm O_2$ flow). Observe the following precautions:
 - a. Start with N2O valve at the maximum setting.
 - b. Adjust only the O2 flow control.
 - c. Decrease the $\rm O_2$ flow as specified in the table and make sure the $\rm O_2$ concentration is in the allowed range.

Note: Allow the $\rm O_2$ monitor to stabilize. At the lower flows, the $\rm O_2$ monitor may take up to 90 seconds to stabilize.

d. If you overshoot a setting, turn the $\rm N_2O$ flow control counterclockwise until the $\rm O_2$ flow increases to the previous setting before continuing the test.

| Set the 0 ₂ flow (L/min) | Measured 0 ₂ |
|-------------------------------------|-------------------------|
| 3.0 | 21% to 30% |
| 2.0 | 21% to 30% |
| 1.0 | 21% to 30% |
| 0.3 | 21% to 30% |

- 11. Check the linearity of the Air flow control.
 - Turn the N₂O flow control fully clockwise to minimum stop.

| Set the 0 ₂ flow (L/min) | Set the Air flow (L/min) | 0 ₂ monitor range |
|-------------------------------------|--------------------------|------------------------------|
| 4.0 | 3.0 | 61% to 71% |
| 3.5 | 6.0 | 45% to 55% |
| 1.5 | 8.0 | 28% to 38% |

3.3.2 Without 0₂ monitoring

⚠ WARNING

The following procedure will test for any significant malfunction of the Link system but it will not confirm proper calibration. Periodic calibration procedures using an accurate and properly calibrated $\rm O_2$ monitor must be performed as recommended in the User's Reference Manual, Part 2, section 3 User Maintenance.

- Nitrous oxide (N₂0) flows through the system during this test. Use a safe and approved procedure to collect and remove it.
 - 1. Set up the gas scavenging system.
 - a. Connect the AGSS to a gas scavenging system.
 - b. Attach a patient circuit and plug the patient port.
 - c. Attach a bag to the bag port (or plug the bag port).
 - d. Set the Bag/Vent switch to Bag.
 - e. Adjust the APL valve to minimum.
 - 2. Connect the pipeline supplies or slowly open the cylinder valves.
 - 3. Turn all flow controls fully clockwise (minimum flow).
 - 4. Set the ACGO selector switch to ABS.
 - 5. Turn on the system.
 - 6. Make sure that:
 - For a dual-tube O₂ flowmeter, the O₂ flowtube shows 0.025 to 0.075 L/min.
 - For a single-tube O₂ flowmeter, the O₂ flowtube shows 0.175 to 0.225 L/min.
 - The other flowtubes show no gas flow.
 - 7. Set the flow controls to mid range of each flowtube and make sure that the flowtube floats rotate and move smoothly.

Note

If the system does not include N2O, skip steps 8 and 9.

- Check the Link proportioning system (increasing N₂O flow).
 Observe the following precautions:
 - a. Start with all valves at the minimum setting.
 - b. Adjust only the N2O flow control.
 - c. Increase the N_2O flow as specified in the following table and make sure the O_2 flow is as specified.

3-6 07/04 1009-0356-000

d. If you overshoot a setting, turn the $\rm O_2$ flow control clockwise until the $\rm N_2O$ flow decreases to the previous setting before continuing the test.

| Set the N ₂ O flow control to (L/min) | The O ₂ flow must be greater than (L/min): |
|--|---|
| 0.8 | 0.2 |
| 2 | 0.5 |
| 4 | 1.0 |
| 10 | 2.5 |

- Check the Link proportioning system (decreasing O₂ flow).
 Observe the following precautions:
 - a. Set the N₂O flow to 9.0 L/min.
 - b. Set the O_2 flow to 3 L/min or higher.
 - c. While reducing the $\rm O_2$ flow, set the $\rm N_2O$ flow to the rates shown in the table. The $\rm O_2$ flow must be greater than the minimum limits.
 - d. If you overshoot a setting, turn the $\rm N_2O$ flow control counterclockwise until the $\rm N_2O$ flow increases to the previous setting before continuing the test.

| Set the N_2 0 flow (using the 0_2 flow control) to (L/min) | The O ₂ flow must be greater than (L/min) |
|--|--|
| 8.0 | 2.0 |
| 4.0 | 1.0 |
| 0.8 | 0.2 |

3.3.3 Pressure relief tests

To check the pressure relief valve (vaporizer manifold outlet).

- 1. Turn all flow controls fully clockwise (minimum flow).
- 2. Set the ACGO selector switch to ACGO.
- Connect a gauge or a digital manometer to the ACGO outlet using the positive pressure leak test adapter.
- 4. Adjust the O₂ flow to 0.5 L/min.
- 5. Verify that the test device reading stabilizes within the following range:
 - 31–60 kPa, 230–450 mm Hg, 4.5–8.5 psi.
- 6. Remove the test device and the adapter.



3.3.4 0₂ supply alarm test

- 1. Set all flow controls to 3 L/min.
- Stop the O₂ supply. (Disconnect the pipeline supply or close the cylinder valve.)
- 3. Make sure that:
 - a. The low ${\rm O}_2$ supply alarm occurs.
 - The N₂O (if equipped) and O₂ flows stop. The O₂ flow stops last.
 - c. Air (if equipped) flow continues.
 - d. Gas supply alarms occur on the ventilator if the ventilator uses \mathbf{O}_2 as the drive gas.
- 4. Turn all of the flow controls fully clockwise (minimum flow).
- 5. Reconnect the pipeline supplies.

3.3.5 Flush Flow Test

With Ventilator

- 1. Set the Bag/Vent switch to Vent.
- 2. Set the system switch to Standby.
- 3. Attach a patient circuit and plug the patient port.
- 4. Set the ACGO selector switch to ABS.
- 5. Ensure that the bellows is completely collapsed.
- Measure the amount of time it takes to fill the bellows when the O₂ Flush button is fully and continuously depressed.
- Repeat the above measurement two more times (deflate bellows by removing the plug from the patient port).
 - The bellows should fill in 1.8 to 2.3 seconds.

3-8 07/04 1009-0356-000

Without Ventilator:


- 1. Set the Bag/Vent switch to Bag.
- 2. Set the system switch to Standby.
- 3. Attach a patient circuit and plug the patient port.
- 4. Attach a 3-liter rebreathing bag on the bag arm or manual bag port (if a 3-liter bag is not available, use the time specified for a 1-liter bag times the volume of the bag used).
- 5. Set the ACGO selector switch to ABS.
- 6. Ensure that the rebreathing bag is completely collapsed.
- 7. Adjust the APL valve to maximum.
- 8. Measure the amount of time it takes to fill the rebreathing bag when the O₂ Flush button is fully and continuously depressed.

Note: When the airway pressure gauge exceeds 10 cm H_20 , the rebreathing bag is full.

- 9. Repeat the above step two more times (deflate the rebreathing bag by removing the plug from the patient port).
 - A 3-liter bag should fill in 3.6 to 5.1 seconds.
 - A 1-liter bag should fill in 1.2 to 1.7 seconds.

Possible Causes of Failure

- · Large leak (if long filling time).
- Flush regulator setting (Section 6.5).
- Flush regulator cross-connection (if long filling time).
- · ACGO selector valve inlet cross-connection (if short filling time).



3.4 Vaporizer back pressure test

⚠ WARNING

Anesthetic agent vapor comes out of the common gas outlet during this test. Use a safe, approved procedure to remove and collect the agent.

- 1. Set up the gas scavenging system.
 - a. Connect the AGSS to a gas scavenging system.
 - b. Attach a patient circuit and plug the patient port.
 - c. Attach a bag to the bag port (or plug the bag port).
 - d. Set the Bag/Vent switch to Bag.
 - e. Adjust the APL valve to minimum.
- 2. Set the ACGO selector switch to ABS.
- 3. Set the system switch to On.
- Set the O₂ flow to 6 L/min.
- 5. Make sure that the ${\rm O}_2$ flow stays constant and the float moves freely.
- 6. Adjust the vaporizer concentration from 0 to 1% one click at a time. The $\rm O_2$ flow must not decrease more than 1 L/min through the full range. If the $\rm O_2$ flow decreases more than 1 L/min:
 - a. Install a different vaporizer and try this step again.
 - b. If the O₂ flow decreases less than 1 L/min with a different vaporizer, the malfunction is in the first vaporizer.
 - c. If the O₂ flow also decreases more than 1 L/min with a different vaporizer, the malfunction is in the Aespire system. Do not use the system until it is serviced (repair vaporizer manifold port valve).
- 7. Complete steps 3 through 5 for each vaporizer and vaporizer position.
- 8. Set the system switch to Standby.

3-10 07/04 1009-0356-000

3.5 Low-pressure leak test

Note

Perform either the "Negative low-pressure leak test" or the "ISO or BSI standard low-pressure leak test. It is not necessary to perform both tests.

⚠ WARNING

Do not use a system with a low-pressure leak. Anesthetic gas will go into the atmosphere, not into the breathing circuit.

3.5.1 Negative lowpressure leak test



- a. Put your hand on the inlet of the leak test device. Push hard for a good seal.
- b. Squeeze the bulb to remove all air from the bulb.
- If the bulb completely inflates in less than 60 seconds, replace the leak test device.
- 2. Set the system switch to Standby.
- 3. Set the ACGO selector switch to ACGO.
- 4. Turn off all vaporizers.




- a. Open the flow controls one and a half turns counterclockwise.
- b. Connect the test device to the ACGO outlet.
- c. Compress and release the bulb until it is empty.
- d. The vacuum causes the floats to move. This is usual. If the bulb completely inflates in 30 seconds or less, there is a leak in the lowpressure circuit.
- 6. Test each vaporizer for low-pressure leaks:
 - a. Set the vaporizer to 1%.
 - b. Repeat step 5.
 - c. Set the vaporizer to OFF.
 - d. Test the remaining vaporizers.
- 7. Disconnect the test device.
- 8. Turn all flow controls fully clockwise (minimum flow). Do not over tighten.

⚠ WARNING

Agent mixtures from the low-pressure leak test stay in the system. Always flush the system with O_2 after the low-pressure leak test (1 L/min for one minute).

Turn off all vaporizers at the end of the low-pressure leak test.

- 9. Flush the system with 0₂:
 - a. Set the system switch to On.
 - b. Set the O₂ flow to 1 L/min.
 - c. Continue the O2 flow for one minute.
 - d. Turn the O₂ flow control fully clockwise (minimum flow).
 - e. Set the system switch to Standby.




3.5.2 ISO or BSI standard low-pressure leak test

⚠ CAUTION

Do the positive pressure leak test at the ACGO outlet only.

- 1. Set the ACGO selector switch to ACGO.
- 2. Turn all flow controls fully clockwise (minimum flow).
- Using the positive pressure leak test adapter, connect the ISO or BSI specific leak test device to the ACGO outlet. Push the adapter into the ACGO outlet throughout the test to get a good seal.



- 4. Keep flowmeter of the test device vertical for accurate results.
- 5. Fully open the needle valve on the test device (counterclockwise).

▲ CAUTION

If the needle valve is not fully open, this test can damage the pressure gauge on the test device.

6. Set the system switch to On.

3-12 07/04 1009-0356-000

- 7. Open the O_2 flow control and set a total flow of 0.4 L/min through the flowmeter on the test device.
- 8. Make sure that the pressure gauge on the test device reads zero and that all other flow controls are fully closed.
- 9. Close the needle valve on the test device until the test gauge reads:

| ISO 5358 | 3 kPa | *** |
|------------|--------|-----|
| BSI 4272.3 | 20 kPa | - |
| | | |

10. If the flow through the test device is less than

0.35 L/min (ISO) or

0.3 L/min (BSI),

there is a low pressure leak in the anesthesia machine.

- Fully open the needle valve on the test device to decrease the back pressure.
- 12. Test each vaporizer for low-pressure leaks:
 - a. Set the vaporizer to 1%.
 - b. Repeat steps 7 through 10.
 - c. Turn the vaporizer OFF.
 - d. Test the remaining vaporizers.
- 13. Remove test tool and adapter.

A WARNING

Agent mixtures from the low-pressure leak test stay in the system. Always flush the system with O_2 after the low-pressure leak test (1 L/min for one minute).

Turn all vaporizers OFF at the end of the low-pressure leak test.

- 14. Flush the system with O₂:
 - a. Set the system switch to On.
 - b. Set the O₂ flow to 1 L/min.
 - c. Continue the O₂ flow for one minute.
 - d. Turn the O₂ flow control fully clockwise (minimum flow).
 - e. Set the system switch to Standby.

3.6 Alarm tests

- Connect a test lung to the patient connection.
- 2. Set the Bag/Vent switch to Vent.
- 3. Set the system switch to On.
- 4. Set the controls:
 - Ventilation Mode: Volume control (select from main menu)
 - · Ventilator:
 - Tidal Vol: 400 ml
 - Rate: 12
 - I:E Ratio:1:2
 - Plimit:40 cm H₂O
 - PEEP: OFF
 - Anesthesia Machine:
 - 02 flow: minimum flow
 - All other gases: OFF
 - ACGO selector switch to ABS
- 5. Push 0₂ Flush to fill the bellows.
- 6. Set the Bag/Vent switch to Bag and back to Vent.
- 7. Make sure that:
 - a. Mechanical ventilation starts.
 - b. A subatmospheric pressure alarm does not occur.

Note: With active gas scavenging, too much scavenging flow can cause subatmospheric alarm.

- c. The ventilator displays the correct data.
- d. The bellows inflate and deflate during mechanical ventilation.
- 8. Set the O₂ flow control to 5 L/min.
- 9. Make sure that:
 - a. The end expiratory pressure is approximately $0\ cm\ H_2 0$. Note: Positive end expiratory pressure when PEEP is off may indicate that the scavenging system is not removing enough gas.
 - b. The ventilator displays the correct data.
 - c. The bellows inflate and deflate during mechanical ventilation.
- 10. Test the low minute volume alarm:
 - a. Go to the alarms menu.
 - b. Set the alarm limit for low minute volume to 6.0 L/min.
 - c. Make sure that a low minute volume alarm occurs.
 - d. Go to the alarms menu.
 - e. Set the low minute volume alarm to OFF.

3-14 07/04 1009-0356-000

11. Test the high airway pressure alarm:

- a. Set P_{limit} to less than the peak airway pressure.
- b. Make sure that the high airway pressure alarm occurs.
- c. Set Plimit to correct level.

12. Test the apnea and low airway pressure alarms:

- a. Turn all flow controls fully clockwise.
- b. Remove the test lung from the patient connection.
- c. Other alarms such as low minute volume can occur.
- d. Make sure that the low airway pressure and apnea alarms occur. The apnea alarm occurs after 30 seconds.

13. Test the sustained airway pressure alarm:

- a. Set the controls:
 - APL valve Closed (70)
 - Bag/Vent switch Bag
- b. Mechanical ventilation stops when the Bag/Vent switch is set to Bag.
- c. Occlude the bag port connector with a test plug.
- d. Close the patient connection using the test plug located on the side of the ABS and push the O_2 Flush button.
- e. Make sure that the sustained pressure alarm occurs after approximately 15 seconds at the sustained pressure limit (6-30 cm H₂O varies with pressure limit).

14. Test the ${\rm O}_2$ monitor and alarms:

- a. Remove the O_2 sensor from the circuit module.
- b. Make sure the sensor measures approximately 21% 0₂ in room air.
- c. Set the low O_2 alarm to 50%. Make sure a low O_2 alarm occurs.
- d. Set the low O_2 alarm back to 21% and make sure that alarm cancels.
- e. Put the O2 sensor back in the circuit.
- f. Remove the test lung from the patient connection.
- g. Set the High O₂ alarm to 50%.
- h. Push the flush button to fill the breathing system.
- i. Set the O₂ flow control to 2 L/min.
- j. Make sure the high O2 alarm comes On.
- k. Set the high 02 alarm back to 100% and make sure that alarm cancels.
- I. After 2 minutes in pure O₂, the O₂ display reads approximately 100%.
- m. Turn the O₂ flow control fully clockwise (minimum flow).
- 15. Set the system switch to Standby.

3.7 Breathing system tests

▲ WARNING

Objects in the breathing system can stop gas flow to the patient. This can cause injury or death:

- Do not use a test plug that is small enough to fall into the breathing system.
- Verify that AGSS is operating. For systems that have a flow indicator on the side, make sure that the flow indicator shows a flow in the green (normal) region.
- 2. Zero the pressure gauge (Section 6.6.1).

Check Valves

- Make sure that the check valves on the breathing circuit module work correctly:
 - The Inspiratory check valve rises during inspiration and falls at the start of expiration.
 - The Expiratory check valve rises during expiration and falls at the start of inspiration.

Ventilator Bellows

- 4. Ventilator bellows test:
 - a. Set the system switch to Standby.
 - b. Set the Bag/Vent switch to Ventilator.
 - c. Set all flow controls to minimum.
 - d. Close the breathing circuit at the patient connection. Use the test plug located on the side of the ABS.
 - e. Push the O₂ flush button to fill the bellows.
 - f. The pressure must not increase to more than 15 cm H₂O on the pressure gauge.
 - g. If the bellows falls more than 100 mL/min (top of indicator), it has a leak.

Service Mode Tests

- Enter the Service Mode: Push and hold the adjustment knob on the ventilator's display and set the system switch to On.
 - a. Select and confirm "Service Modes."
 - Follow the menu structure outline below to reach the adjustment for the inspiratory flow valve. Select and confirm at each step.
 - "Diagnostics Tests/Tools"
 - "Breathing System Leak Test"
 - c. Follow the instructions on the screen.
 - d. The leak rate should be less than 200 mL/min.
 For machines with a single-tube O₂ flowmeter, the pressure reading should reach 30 cm H₂O at minimum flows greater than 200 mL/min.

Note: If test fails, see Section 7.2, "Breathing System Leak Test Guide".

3-16 07/04 1009-0356-000

Bag Circuit 6. Te

- 6. Test the Bag circuit for leaks:
 - a. Set the system switch to On.
 - b. Set the Bag/Ventilator switch to Bag.
 - c. Plug the Bag port (use your hand or the approved test plug).
 - d. Close the APL valve (70 cm H₂0).
 - e. Set the O₂ flow to 0.25 L/min.
 - f. Close the patient connection (using a hand or test plug on the side of the breathing system) and pressurize the bag circuit with the O₂ flush button to approximately 30 cm H₂O.
 - g. Release the flush button. The pressure must not decrease. A pressure decrease large enough to see on the gauge indicates an unacceptable leak.

Note: If test fails, see Section 7.2, "Breathing System Leak Test Guide".

APL Valve

- 7. Test the APL valve:
 - a. Fully close the APL valve (70 cm H₂0).
 - b. Set the total fresh gas flow to approximately 3 L/min and make sure that the value on the inspiratory pressure gauge does not exceed 85 cm H₂O. Some pressure fluctuation is normal.
 - c. Fully open the APL valve (to the MIN position).
 - d. Set O2 flow to 3 L/min. Turn any other gases off.
 - e. Make sure that the value on the inspiratory pressure gauge is less than approximately $5\ cm\ H_2O$.
 - f. Push the O₂ flush button. Make sure that the value on the inspiratory pressure gauge stays less than 10 cm H₂O.
 - g. Set the O₂ flow to minimum and make sure that the value on the inspiratory pressure gauge does not decrease below 0 cm H₂O.
- 8. Remove your hand or the test plug from the patient connection.
- 9. Set the System switch to Standby.

⚠ WARNING

Make sure that there are no test plugs or other objects caught in the breathing system.

3.8 Auxiliary 02 flowmeter tests

- 1. Open the O_2 cylinder valve or connect an O_2 pipeline.
- Rotate the flow control clockwise (decrease) to shut off the flow. The ball should rest at the bottom of the flow tube and not move.
- Rotate the flow control counterclockwise (increase). The ball should rise immediately after rotation is begun. It should rise smoothly and steadily with continued counterclockwise rotation. When a desired flow is set, the ball should maintain in a steady position.
- Occlude the auxiliary O₂ outlet. The ball should rest at the bottom of the flow tube and not move. A ball that does not rest at the bottom of the flow tube indicates a leak and requires service.
- 5. Rotate the flow control clockwise to shut off the flow.

3.9 Integrated Suction Regulator tests


The gauge needle should come to rest within the zero range bracket when no suction is being supplied. Gauges which do not comply may be out of calibration.

- 1. Adjust the regulator setting to minimum.
- 2. Turn the mode selector to I (On).
- 3. Ensure the gauge remains less than 200 mmHg (26 kPa, 0.26 Bar).
- 4. Occlude the inlet.
- 5. Ensure the gauge remains less than 200 mmHg (26 kPa, 0.26 Bar).
- 6. Adjust the regulator in an increasing vacuum level.
- 7. The gauge should rise after rotation has begun. The gauge should rise with continued rotation of the regulator adjustment.
- 8. Adjust the regulator setting to minimum.
- 9. Turn the Mode selector to O (Off).

3-18 07/04 1009-0356-000

3.10 Power failure test

 Connect the power cord to a mains outlet. The mains indicator on the display comes on when AC Power is connected.



If the indicator is not on, the display assembly is not receiving AC power.

- Verify AC power to the machine (reset circuit breaker).
- · Check fuses in display assembly's inlet module.
- · Check fuses in machine's AC inlet assembly.
- 2. Set the system switch to On.
- 3. Unplug the power cord with the system turned on.
- 4. Make sure that the power failure alarm comes on.
- 5. Make sure the following message is on the ventilator display:
 - "On Battery Power OK?"
- 6. Connect the power cable again.
- 7. Make sure the alarm cancels.

3.11 Electrical safety tests

Make sure the system is completely assembled and all accessory devices are connected to electrical outlets.

 Connect an approved test device (e.g. UL, CSA, or AAMI) and verify that the leakage current is less than:

| Voltage | Max. Leakage Current | |
|-------------|----------------------|--|
| 120/100 Vac | $300\mu\text{Amps}$ | |
| 220/240 Vac | 500 μAmps | |

2. Make sure that the resistance to ground is less than 0.2Ω between an exposed metal surface and the ground pin on the power cord.

Notes

3-20 07/04 1009-0356-000

| In this section | This section covers the repair and replacement procedures for components of the Aespire anesthesia machine. |
|-----------------|---|
| | 4.1 Servicing the ventilator |
| | 4.2 How to bleed gas pressure from the machine |
| | 4.3 How to remove the rear panels |
| | 4.3.1 To remove the rear upper panel |
| | 4.3.2 To remove the lower access panels |
| | 4.4 How to remove the tabletop |
| | 4.5 Replace pipeline inlet filter |
| | 4.5.1 Replace pipeline inlet check valve4-6 |
| | 4.6 Change drive gas |
| | 4.7 Service the cylinder supply modules |
| | 4.7.1 Tightening procedure for high-pressure tube fittings4-8 |
| | 4.7.2 Replace primary regulator module (complete replacement)4-8 |
| | 4.7.3 Replace cylinder inlet filter4-9 |
| | 4.7.4 Replace cylinder check valve |
| | 4.7.5 Replace 3rd-gas cylinder supply module |
| | 4.8 Replace system switch assembly |
| | 4.9 Service the flowmeter module |
| | 4.9.1 Remove front flowmeter panel shield |
| | 4.9.2 Remove flowtubes for cleaning or replacement |
| | 4.9.3 Remove complete flowmeter head |
| | 4.9.4 Replace flowmeter modules |
| | 4.9.5 Replace flowmeter frame |
| | 4.9.6 Replace O ₂ supply switch 4-21 |
| | 4.9.7 Checkout procedure for O ₂ supply switch |
| | 4.9.8 Replace secondary regulator manifold or balance regulator manifold 4-22 |
| | 4.9.9 Replace O ₂ or N ₂ O needle valves (on machines with N ₂ O) |
| | 4.9.10 Replace an Air needle valve on all machines or an O_2 needle valve on machines without N_2O |

| 4.10 Service vaporizer manifold parts | 4-26 |
|---|------|
| 4.10.1 Repair manifold port valve | 4-26 |
| 4.10.2 Checkout procedure for manifold port valve | 4-27 |
| 4.10.3 Replace vaporizer manifold check valve | 4-28 |
| 4.10.4 Replace vaporizer pressure relief valve | 4-30 |
| 4.10.5 Replace vaporizer manifold | |
| 4.11 Replace ACGO selector switch | 4-32 |
| 4.12 Clean or replace ACGO port flapper valve | 4-34 |
| 4.13 Reconfigure sample gas return line | 4-35 |
| 4.14 Replace the APL valve | 4-36 |
| 4.15 Replace the bag support arm | 4-37 |
| 4.15.1 Servicing the bag support arm | 4-38 |
| 4.15.2 Replace friction pad in lower bag arm assembly | 4-39 |
| 4.15.3 Replace bag port housing | 4-40 |
| 4.16 Replace auxiliary O ₂ flowmeter | 4-41 |
| 4.17 Replace the suction control module | 4-42 |
| 4.17.1 Front panel method | 4-42 |
| 4.17.2 Rear panel method | 4-43 |
| 4.18 Replace ABS breathing system components | 4-44 |
| 4.18.1 Replace Bag/Vent switch assembly | 4-44 |
| 4.18.2 Replace bellows base latch assembly | 4-45 |
| 4.19 Replace casters | 4-46 |
| 4.20 Replace task light and switch | 4-47 |
| 4.20.1 To replace the task-light switch | 4-47 |
| 4.20.2 To replace the task-light circuit board | 4-47 |
| 4.21 Replace the display arm or display cables | 4-48 |
| 4.21.1 Cable tie installation | 4-48 |
| 4.21.2 Removing the display arm | 4-49 |
| 4.21.3 Replacing a display cable | 4-49 |
| 4.21.4 Installing the long arm | 4-50 |
| 4.21.5 Installing the short arm | 4-51 |
| 4.22 Replace display and cables in ProTIVA machine | 4-52 |

4-2 07/04 1009-0356-000

⚠ WARNING

To prevent fires:

- Use lubricants approved for anesthesia or O₂ equipment, such as Krytox.
- Do not use lubricants that contain oil or grease; they burn or explode in high O₂ concentrations.
- All covers used on the system must be made from antistatic (conductive) materials. Static electricity can cause fires.
- Obey infection control and safety procedures. Used equipment may contain blood and body fluids.
- A movable part or a removable component may present a pinch or a crush hazard. Use care when moving or replacing system parts and components.
- Some internal parts have sharp edges and can cause cuts or abrasions. Use care when servicing internal components.
- After repairs are completed, always perform the checkout procedure. Refer to Section 3 of this manual.

4.1 Servicing the ventilator

The Aespire anesthesia machine is configured with the 7100 Ventilator.

Service information for this ventilator is provided in a separate manual as detailed in Section 1.2.2.

The 7100 Ventilator comprises the following components:

- The Display/Control module.
- The Serial Isolation and Connector Board, located in the upper structure.
- The Vent Engine, located in the rear portion of the breathing system.
- The Ventilator Monitoring Board, located under the tabletop.

Common machine components such as switches and sensors which provide input to the ventilator are covered in this manual.

4.2 How to bleed gas pressure from the machine

Before disconnecting pneumatic fittings, bleed all gas pressure from the machine.

- 1. Set the system switch to On.
- Close all cylinder valves and disconnect all pipeline supplies from the source.
 Note: If the machine includes N₂O, do not disconnect the O₂ pipeline.
 If pipeline O₂ is not available, open the O₂ cylinder valve.
- Turn the flow controls for all gases (except O₂) at least one turn counterclockwise.
- 4. Ensure that all cylinder and pipeline gauges read zero before proceeding.
 - For machines with N₂O, disconnect the O₂ pipeline supply from the source (or close the O₂ cylinder valve).
 - Press the O₂ flush button to bleed O₂ from the system.
- 5. Set the system switch to Standby.

4.3 How to remove the rear panels

You must remove the rear upper panel to repair or replace many of the machine's components. To access the 3rd cylinder supply (if equipped) you must remove the lower access panels.

4.3.1 To remove the rear upper panel

- 1. Bleed all gas pressure from the machine (Section 4.2).
- 2. Ensure that all cylinder and pipeline gauges read zero before proceeding.
- 3. Disconnect all electrical cables.
- 4. To remove the rear panel, fully loosen the three captive screws that hold the panel in place. Remove the panel.
 - If the machine includes integrated suction, disconnect the two tube fittings from the overflow safety trap manifold.
 - If the machine includes electrical outlets, lower the panel and place it so that it does not stress the power cable.

4.3.2 To remove the lower access panels

- 1. Bleed all gas pressure from the machine (Section 4.2).
- 2. Ensure that all cylinder and pipeline gauges read zero before proceeding.
- 3. If present, remove the inboard cylinders.
- Remove (two screws) the small upper access panel to access the 3rd cylinder regulator (N₂O) test port.
- Remove (six screws) the large lower access panel to access the 3rd cylinder regulator.

4-4 07/04 1009-0356-000

4.4 How to remove the tabletop

The tabletop is held in place with five captive screws along the periphery of the pan assembly (accessed from below the rim of the tabletop).

- One screw (A) is in a deep recess at the right-rear corner of the tabletop.
- Two screws (B) are at the front of the tabletop: one screw is at the right corner of the tabletop, one is near the O₂ Flush button.
- To access the remaining two screws (**C**), you must remove the ABS: one screw is at the left corner of the tabletop, one is near the APL Valve.



4.5 Replace pipeline inlet filter


- 1. Remove the pipeline inlet fitting.
- 2. Pull the pipeline inlet filter out of the fitting. The o-ring should come out with the filter.




3. Install the new pipeline inlet filter in the pipeline inlet fitting. The new filter comes with an o-ring.

4.5.1 Replace pipeline inlet check valve

- 1. Remove the rear panel (Section 4.3).
- 2. Remove the pipeline inlet fitting.
- 3. The Air and O₂ pipeline manifolds include a drive gas connection at the back of the manifold. Remove the drive gas tube or plug to access the check valve.
- 4. From the back of the pipeline manifold, use a thin tool to push out the check valve. (For an N₂O manifold, you will have to carefully apply pressure at the outlet of the manifold – with a syringe for example – to gently force the check valve out of the manifold).



Push the new check valve into the opening, using the same thin tool. The new check valve includes an o-ring — orient it toward the pipeline inlet. Note: Make sure to push the new check valve all the way back into the opening until it bottoms out on the shoulder.



- Install the pipeline inlet fitting.
- 7. Perform the checkout procedure (Section 3).

4-6 07/04 1009-0356-000

4.6 Change drive gas

⚠ CAUTION

If you change the drive gas, you must also change the drive gas selection on the ventilator service setup screen. Refer to Section 4 of the ventilator Technical Reference manual.

 If the drive gas selection and the actual drive gas do not agree, volumes will not be correct.

The ventilator will alarm with the message "Low Drive Gas Press" if the selected drive gas pressure, either O_2 or Air, is lost.

1. Remove the rear panel (Section 4.3).

Note:

The $\rm O_2$ and Air pipeline manifolds have a drive gas connection at the back. The connection not in use is plugged.

- 2. Remove the plug from the new connection.
- 3. Disconnect the drive gas hose from the present connection.
- 4. Install the plug in this connection (pull on the plug to ensure that it is locked into the fitting).
- 5. Reroute the drive gas hose so that it does not cause kinks in other tubing.
- 6. Connect the drive gas hose to the new connection (pull on the hose connector to ensure that it is locked into the fitting).
- 7. Do a high-pressure leak test (Section 3.2).
- 8. Enter the service mode and select the correct drive gas.
- 9. Test the primary regulator. Verify that it functions within specifications now that it will be supplying drive gas to the ventilator (Section 6.1).
- 10. Perform the checkout procedure (Section 3).

4.7 Service the cylinder supply modules

△ WARNING

Be careful not to expose internal components to grease or oil (except Krytox or equivalent).

4.7.1 Tightening procedure for high-pressure tube fittings

The cylinder pressure gauge is connected to the cylinder supply through a copper tube with fittings at both ends. Use the following tightening procedure whenever you are replacing a cylinder supply or a cylinder pressure gauge.


- 1. Insert the tubing into the fitting until the ferrule seats in the fitting.
- 2. Tighten the nut by hand.
- 3. Continue tightening the nut with a wrench until it reaches the original position (about 1/4 turn). You will feel an increase in resistance at the original position.
- 4. After reaching the original position, tighten the nut just slightly.

Note

If you are installing a new tube that has not been tightened before, tighten the nut with a wrench an additional 3/4 of a turn after the nut is finger tight.

4.7.2 Replace primary regulator module (complete replacement)

- 1. Bleed all gas pressure from the machine (Section 4.2).
- 2. Ensure that all cylinder and pipeline gauges read zero before proceeding.
- 3. Remove the rear panel (Section 4.3).
- 4. Disconnect the high-pressure cylinder gauge fitting.
- 5. Disconnect the output tube fitting.
- Remove the three mounting screws and lockwashers.
- To reassemble, perform the previous steps in reverse order.
 - Tighten the high-pressure fitting as detailed in Section 4.7.1
 - Pull on the cylinder output fitting to ensure it is locked in place.
- Check the output of the regulator BEFORE you install the rear panel. Adjust if necessary (Section 6.1).
- Perform the checkout procedure (Section 3).



4-8 07/04 1009-0356-000

4.7.3 Replace cylinder inlet filter


- 1. Open the cylinder yokes.
- 2. Remove the inlet adapter from the cylinder yoke, using a 4 mm hex wrench.

Note: A brass retaining ring keeps the filter inside the inlet adapter.

3. Thread a 6-mm screw (two turns only) into the brass retaining ring and pull it out.

△ CAUTION

Be careful not to crush the filter. Do not thread in the screw more than two full turns.




- 4. Remove the filter.
- 5. Install the new filter and brass retaining ring.
- 6. Install the inlet adapter in the cylinder yoke.
- 7. Perform the checkout procedure (Section 3).

4.7.4 Replace cylinder check valve

The cylinder check valve is not a replaceable item. If the check valve is defective, you must replace the complete cylinder supply module.

4.7.5 Replace 3rd-gas cylinder supply module

- 1. Remove the lower rear cover (Section 4.3.2).
- 2. Disconnect the high-pressure cylinder gauge fitting.
- 3. Disconnect the output tube fitting.
- 4. Remove the three mounting screws and lockwashers.




- $5. \quad \text{To reassemble, perform the previous steps in reverse order.} \\$
 - Tighten the high-pressure fitting as detailed in Section 4.7.1
 - Pull on the cylinder output fitting to ensure it is locked in place.
- Check the output of the regulator BEFORE you install the rear cover. Adjust if necessary (Section 6.1).
- 7. Perform the checkout procedure (Section 3).

4-10 07/04 1009-0356-000


4.8 Replace system switch assembly

- 1. Bleed all gas pressure from the machine (Section 4.2).
- 2. Ensure that all cylinder and pipeline gauges read zero before proceeding.
- 3. Remove the tabletop (Section 4.4).
- 4. Remove the gauge panel mounting screws and move the panel forward to access the system switch.



- 5. Disconnect the wires from the electrical switch.
- 6. Back out the system switch mounting screws just enough to allow the knob collar to be released.
- 7. While holding the switch assembly, push in the knob and turn it counterclockwise.
- 8. Pull the knob and collar out from the front and remove the switch assembly.

- 9. Install the replacement switch assembly:
 - a. Loosen the two outside screws on the electrical module.
 - b. Insert the wires in the electrical module and tighten the screws.
 - c. Pull the wires on the electrical module to ensure that there is a good connection.
 - d. Turn back the system switch mounting screws until their tips recede.
 - e. Orient the switch assembly with the Air fittings toward the right and the $\rm O_2$ fittings toward the left.
 - f. Install the switch assembly through the gauge panel.
 - g. Push the knob collar in with the indicator up and turn it clockwise until it locks.
 - h. Tighten the mounting screws. Make sure that the top edge of the switch assembly is parallel to the top edge of the gauge panel.
 - Transfer the tubing from the old system switch to the new system switch on the pneumatic module (pull on the tubing to ensure that it is locked into the module).




- 10. Test the replacement switch assembly:
 - a. Connect Air and O2 supplies.
 - b. Connect the power cable to an electrical outlet.
 - c. Set the system switch to On.
 - d. Increase the O₂ and Air flow. Make sure that gas flows.
 - e. Make sure that you do not feel or hear any leaks.
 - f. Make sure that the display comes On.
 - g. Set the system switch to Standby.
 - h. Make sure all gas flow stops and the display turns Off.
- 11. Reinstall the gauge panel and the tabletop.
- 12. Perform the checkout procedure (Section 3).

4-12 07/04 1009-0356-000

4.9 Service the flowmeter module

4.9.1 Remove front flowmeter panel shield


- 1. Bleed all gas pressure from the machine (Section 4.2).
- 2. Ensure that all cylinder and pipeline gauges read zero before proceeding.
- The flowmeter panel is held in place with two latching tabs at the right side. To remove the panel, release each latch by pushing it toward the center of the panel with a thin rod (3-mm hex wrench) through the access hole in the shroud.



- 4. Remove the panel.
- 5. To reinstall the panel, engage the retaining tabs on the left side and press the right side against the shroud to latch it in place.

4.9.2 Remove flowtubes for cleaning or replacement


- 1. Remove the front flowmeter panel shield (Section 4.9.1).
- To remove a flowtube assembly, push up on the tube just enough to clear the bottom seal, pull out from the bottom until the tube clears the flowtube module, then pull down slowly to release it from the module.



⚠ WARNING

Floats are calibrated to a specific tube. Keep each float with its tube. Replace tube and floats together. Interchanging floats can cause incorrect readings. Disassemble the flowtube assemblies only when service is required. Excessive cleaning can remove the antistatic coating from inside the tube. Damage to the float requires replacement of the entire flowtube.

3. Disassemble the flowtube assembly.



- Completely clean, rinse, and dry the flowtube. Use hospital grade alcohol and a flowtube brush.
- 5. Replace stops, filter and o-rings, as necessary. Lightly coat all o-rings with Krytox. Be careful to not get any Krytox on the bottom of small flowtube float stops.
- 6. Reassemble the flowtube assembly.
- 7. Insert the flowtube, spring side first, into the top of the module with the scale oriented forward.
- 8. Push up and slide the bottom of the flowtube into place on the bottom o-ring. It may be necessary to rotate the tube to engage the index tabs.


Note: Be sure o-rings are inserted completely into the collar.

- 9. Push down on the tube to seat the bottom o-ring.
- Reinstall the front flowmeter panel shield.
- 11. Perform the checkout procedure (Section 3).

4-14 07/04 1009-0356-000

4.9.3 Remove complete flowmeter head


- 1. Remove the rear panel (Section 4.3).
- Disconnect the tubing at the rear of each gas module. The following example is a back view of the flowmeter head.



3. Disconnect the O₂ supply switch harness. Note position of switch connections so that you can reassemble correctly later.




- 4. Remove the front flowmeter panel shield (Section 4.9.1).
- 5. Remove the mounting screw from each regulator manifold.



- 6. Pull the flowhead forward.
- 7. To reinstall, perform the previous steps in reverse order (pull on the tubing to ensure it is locked into the fittings).
- Check for proper alignment of the front flowmeter panel. If any of the needle valve knobs rub against the flowmeter panel, reposition the flowhead to allow for proper clearance.
- 9. Perform the checkout procedure (Section 3).

4.9.4 Replace flowmeter modules

- 1. Remove the complete flowmeter head (Section 4.9.3).
- Refer to the following illustrations. Note that these illustrations show ANSI flowmeter module positions. The order is reversed in ISO machines.



0₂/Air modules


 O_2 /Air flowmeter modules are connected at the bottom with a long screw (**A**) and nut (**B**) that is recessed. To remove, retain nut while loosening screw.



Note

The flowmeter modules are held together at the top with a u-clip. To separate the modules, pivot the modules (front to back) 45 degree. The u-clip will disengage and allow the modules to separate.


The flowmeter modules are interconnected at the top by a spacer tube. The o-ring on the spacer tube makes a leak-tight seal.



4-16 07/04 1009-0356-000

Single-tube flowhead

The outlet fitting for a single-tube flowhead is not an integral component of the $\rm O_2$ flowmeter. The outlet fitting is a separate component that includes an o-ring seal and is held in place with two u-clips.



3. To remove the Air flowmeter module:


- a. Remove the screw located on the side of the Air flowmeter module.
- b. Hold the flowmeter module with the flowtubes facing you.
- c. Grasp the outer modules at the bottom of the regulator manifold and push the left module away from you until the u-clip pops off and the module separates from the other assemblies.
- Pull the modules sideways to separate them at the top.
 Save the u-clip, spacer tube, and the o-ring for reassembly.

4. To remove the O2 or N2O flowmeter module:


- a. Set the O_2 and N_2O needle valves to their maximum position (counterclockwise).
- b. Loosen the set screws on the N₂O knob, then remove the knob.
- c. Loosen the set screws on the N_2O sprocket and the O_2 knob.
- d. To remove, grasp the $\rm O_2$ knob/sprocket, $\rm N_2O$ sprocket, and chain as an assembly. Remove as an assembly.



- e. Remove the spacer from the N₂O needle valve spindle.
- f. Remove the pilot tube going to the balance regulator.



- g. The $\rm O_2$ and $\rm N_2O$ flowmeter modules are held together by a single screw. Remove the screw located on the side of the $\rm O_2$ flowmeter module.
- h. Hold the flowmeter modules with the flowtubes facing you.
- Grasp the modules at the bottom of the regulator manifolds and push the left module away from you until the N₂O module separates from the O₂ module.
- Pull the modules sideways to separate them at the top.
 Save the u-clip, spacer tube, and the o-ring for reassembly.
- To reassemble the flowmeter modules, perform the previous steps in reverse order.Note: The u-clips must be reinstalled with the barbed leg to the left as viewed from the front.




- Install the screw through the O₂ flowmeter module to locks the O₂ and N₂O flowmeter modules together.
- 7. Reattach the pilot tube that goes to the balance regulator.
- Install the flowhead into the machine. Reconnect the tubing and the O₂ supply switch harness.
- 9. Confirm needle valve calibration (Section 6.3).
- 10. Install the spacer on the N2O needle valve stem.
- 11. Install the chain on the O2 knob/sprocket assembly and the N2O sprocket.
- 12. Install the chain and sprockets on the needle valve stems as an assembly. Do not tighten the set screws.

4-18 07/04 1009-0356-000








- 13. Install the ${\rm N_2O}$ knob. Snug one set screw to hold the knob in place.
- 14. Perform the link system calibration (Section 6.4).
- 15. Install the flowmeter panel shield.
- 16. Perform the checkout procedure (Section 3).

4.9.5 Replace flowmeter frame

- 1. Remove the front flowmeter panel shield (Section 4.9.1).
- 2. Remove the complete flowmeter head (Section 4.9.3).
- 3. Separate the flowmeter modules as required (Section 4.9.4).
- 4. Remove the flowtubes (Section 4.9.2). Keep all the parts for reassembly.
- 5. Remove the gas identification panel by removing the two screws at the back of the frame. Keep all the parts for reassembly.



Remove the flowmeter frame by loosening the four mounting screws at the back of the regulator manifold.

Note: There is a retainer in each screw location that keeps each screw within the manifold.

- 7. To reassemble, perform the previous steps in reverse order.
- 8. If replacing O_2 or N_2O frames, perform the link system calibration (Section 6.4).
- 9. Perform the checkout procedure (Section 3).


4-20 07/04 1009-0356-000

4.9.6 Replace 0₂ supply switch

The O_2 supply switch is located on the O_2 flowmeter's regulator module.

- 1. Remove the upper rear panel (Section 4.3).
- 2. Remove the two mounting screws from the O_2 supply switch.





- 3. Pull the O₂ supply switch out of the regulator manifold.
- 4. Install the replacement O2 supply switch.
- 5. Tighten the screws.
- 6. Disconnect the leads from the old switch and reconnect them to the new switch.
- Adjust the alarm threshold for the new O₂ supply switch, as explained in the checkout procedure below (Section 4.9.7).
- 8. Replace the rear panel.
- 9. Perform the checkout procedure (Section 3).

4.9.7 Checkout procedure for O₂ supply switch

- 1. Remove the upper rear panel (Section 4.3).
- 2. Attach a gauge to the $\rm O_2$ primary regulator test port. (On pipeline only machines, attach the gauge to a 6-mm $\rm O_2$ port).
- 3. Adjust the O₂ flow control to minimum flow (clockwise).
- Install an O₂ cylinder and open the cylinder valve (for pipeline only, connect O₂ pipeline source).
- 5. Turn the system on.
- Close the cylinder valve (disconnect pipeline from source) and watch the test gauge as the O₂ pressure bleeds down slowly.

Note: The "No 0_2 pressure" alarm should occur between descending pressure of 221-193 kPa (32-28 psi).

- 7. If adjustment is required, set the adjustment screw so that the "No O_2 pressure" alarm occurs at 207 \pm 7 kPa (30 \pm 1 psi).
- 8. Disconnect the gauge and plug the test port
- 9. To reassemble, perform the previous steps in reverse order.
- 10. Perform the checkout procedure (Section 3).

4.9.8 Replace secondary regulator manifold or balance regulator manifold

- 1. Remove the front flowmeter panel shield (Section 4.9.1).
- 2. Remove the complete flowmeter head (Section 4.9.3).
- 3. Separate the flowmeter modules (Section 4.9.4).
- Remove the flowmeter frame from the regulator manifold by removing the four screws at the rear of the regulator manifold (no need to remove flowtubes).
- 5. Remove the needle valve from the regulator manifold:
 - a. Unscrew the complete assembly together (stop collar, needle valve).
 - b. Replace the o-ring if necessary.
- 6. Screw the needle valve into the new regulator manifold.
- 7. Remove the plugs and balance regulator elbow fitting (and $\rm O_2$ supply switch if an $\rm O_2$ module) from the old regulator manifold.
- Install the plugs and balance regulator elbow fitting (and O₂ supply switch if an O₂ module) into the new regulator manifold (pull on the plugs and fittings to ensure that they are locked into the manifold).
- 9. Reinstall the flowmeter frame to the regulator manifold.
- 10. Reinstall all the flowmeter modules to the flowmeter head.
- 11. Reinstall the flowmeter head (Section 4.9.4).
- 12. Do the necessary calibrations (Section 6).

| Necessary calibrations | Section |
|-----------------------------|---------|
| Secondary Regulator | 6.2 |
| O ₂ minimum flow | 6.3.1 |
| Maximum flow | 6.3.4 |
| Link system | 6.4 |

- 13. Reinstall the front flowmeter panel shield.
- 14. Perform the checkout procedure (Section 3).

4-22 07/04 1009-0356-000

4.9.9 Replace 0₂ or N₂O needle valves (on machines with N₂O)


For machines without N_2O , refer to Section 4.9.10 for replacing the O_2 needle valve.


- 1. Bleed all gas pressure from the machine (Section 4.2).
- 2. Ensure that all cylinder and pipeline gauges read zero before proceeding.
- 3. Remove the front flowmeter panel shield (Section 4.9.1).
- 4. Set the O₂ and N₂O needle valves to their minimum position.
- Loosen the set screws on the N₂O knob, then remove the knob.
- Loosen the set screws on the N₂O sprocket and the O₂ knob.
- To remove, grasp the O₂ knob/sprocket, N₂O sprocket, and chain as an assembly.
 Remove as an assembly.
- 8. Remove the spacer from the N₂O needle valve spindle.
- Loosen the set screws on the needle valve stop collar for the needle valve that is being replaced.
- 10. Remove the stop collar.
- To remove the needle valve from the flowmeter block, turn the needle valve counterclockwise with a 16-mm wrench.
- 12. To install the new needle valve, turn it clockwise and tighten it with the wrench.

Note: Be sure the o-ring is properly located on the tip of the needle valve.

⚠ WARNING

The $\rm O_2$ and $\rm N_2O$ needle valves are not the same. Patient injury can result if the wrong needle valve is installed in the flowmeter block. You can identify the $\rm N_2O$ needle valve by a groove located just below the top brass hex.






- 13. Install the stop collar on the new needle valve. Do not tighten the set screws.
- 14. Perform the needle valve calibration (Section 6.3).

15. After calibrating the needle valve, put the spacer the $N_2\mathrm{O}$ needle valve spindle.



- 16. Put the chain on the O_2 knob/sprocket assembly and the $\mathrm{N}_2\mathrm{O}$ sprocket.
- 17. Install the chain and sprockets on the needle valve spindles as an assembly. Do not tighten the set screws.





- 18. Install the N_2O knob. Do not tighten the set screws.
- 19. Perform the link system calibration (Section 6.4).
- 20. Install the flowmeter panel shield.
- 21. Perform the checkout procedure (Section 3).

4-24 07/04 1009-0356-000

4.9.10 Replace an Air needle valve on all machines or an O_2 needle valve on machines without N_2O

For machines with N_2O , refer to Section 4.9.9 for replacing the O_2 needle valve.

- 1. Bleed all gas pressure from the machine (Section 4.2).
- 2. Ensure that all cylinder and pipeline gauges read zero before proceeding.
- 3. Remove the front flowmeter panel shield (Section 4.9.1).
- 4. Loosen the set screws on the flow control knob and remove the knob.
- 5. Loosen the set screws on the stop collar and remove the collar.
- 6. If equipped, remove the maximum flow stop collar.
- 7. Using a 16-mm wrench, remove the needle valve by turning it counterclockwise.
- 8. Install the new needle valve and tighten.


Note: Be sure the o-ring is properly located on the tip of the needle valve.

- 9. If equipped, install the maximum flow stop collar (do not tighten).
- 10. Install the stop collar (do not tighten the screws).
- 11. Install the flow control knob on the shaft. Tighten one set screw to snug.
- 12. Reconnect the gas supplies.
- 13. Perform the flow control stop procedures explained in:
 - Section 6.3.1 for 0₂.
 - · Section 6.3.3 for Air.
 - · Section 6.3.4 for maximum flow.
- 14. Install the flowmeter panel shield.
- 15. Perform the checkout procedure (Section 3).

4.10 Service vaporizer manifold parts


4.10.1 Repair manifold port valve

- 1. Set the system switch to Standby.
- 2. Remove the vaporizers from the vaporizer manifold.
- 3. Using a 14-mm wrench, carefully remove the valve nipple (threaded).



 Disassemble as necessary to replace parts. The following illustration shows the parts.

Note: The port valve replacement kit includes the valve cartridge assembly and the seal. The kit does not include o-rings.



- 5. When installing a new valve cartridge assembly into the vaporizer manifold, put a light coat of Krytox on the bottom portion of the cartridge. The bottom portion of the cartridge is defined as the brass surface that is inserted in the lower spring. Note: Do not apply Krytox to the valve seal.
- 6. Verify that the parts are free of dust and dirt.
- 7. To reassemble, perform the previous steps in reverse order.
- 8. Complete the port valve checkout procedure described below (Section 4.10.2).

4-26 07/04 1009-0356-000

4.10.2 Checkout procedure for manifold port valve

Use the Vaporizer Manifold Valve Test Tool to perform the checkout procedure for the manifold port valve. This tool and test procedure are intended for use only when the valve cartridge assembly is replaced.

Note

This replacement and test procedure is a service action and is not part of the maintenance program.


- 1. Set the system switch to Standby.
- 2. After replacing the valve cartridge assembly, remove the vaporizer port o-ring.
- 3. Attach the valve tester to the top of the valve by sliding the bottom of the tester onto the o-ring groove.
- 4. Tighten the tester screw down onto the valve until the screw bottoms out on the top of the valve. The tester o-ring should create a seal with the top of the valve.
- 5. Remove the upper rear panel (Section 4.3).
- 6. Remove the inlet tube from the vaporizer manifold.
- 7. Set the SCGO Selector switch to ACGO.
- 8. Test the negative low-pressure leak-test device:
 - a. Put you hand on the inlet of the leak-test device. Push hard for a good seal.
 - b. Remove all air from the bulb.
 - c. The bulb should not inflate in less than 60 seconds.
- 9. Attach the negative low-pressure leak-test device to the ACGO outlet.
- 10. Remove all air from the bulb. The bulb should not inflate in less than 45 seconds.
- 11. Remove the valve tester.
- 12. Reassemble the inlet tube, vaporizer port o-ring, and the upper rear panel.
- 13. Conduct a negative low-pressure leak test on the system (Section 3.5.1).

⚠ WARNING


If the valve test tool is not removed before flowing gas through the system, pneumatic head damage could result.

4.10.3 Replace vaporizer manifold check valve

- 1. Set the system switch to Standby.
- 2. Remove the vaporizers from the vaporizer manifold.
- 3. Remove the upper rear panel.
- 4. Disconnect the tubing from the valve block.




- 5. Remove the valve block.
 - To access the left-hand mounting screw (A), remove the right (viewed from front) side panel (B).
 - The right-side extrusion include an access hole (C) for removing the left-hand mounting screw.
 - Note: For early production machines that do not have an access hole in the extrusion, you must remove the vaporizer manifold to remove the valve block.



4-28

Note The valve body, o-ring, and flapper do not come out with the block. They stay intact at the bottom of the vaporizer manifold.

6. Pull the flapper out of the valve body.



- 7. Using a hex wrench, put the wrench through one of the holes of the valve body and pull down to remove the valve body and o-ring.
- 8. Verify that parts are free of dust and dirt.
- 9. Replace the flapper by inserting the flapper stem and gently pulling the stem until the flapper secures to the valve body.
- 10. Lightly lubricate the o-ring with Krytox.
- 11. Place the lubricated o-ring on the valve body port at the bottom of the manifold.
- 12. Gently install the valve body in the manifold:
 - Check that the o-ring makes a good seal between the manifold and the valve body.
 - Check that the flapper valve makes solid contact with the valve body.
- 13. Install the valve block.
- 14. Reconnect the tubing to the valve block. Pull on the tube to ensure that it is locked in the fitting.
- 15. Install the vaporizer front panel.
- 16. Perform the checkout procedure (Section 3).

4.10.4 Replace vaporizer pressure relief valve

- 1. Set the system switch to Standby.
- 2. Remove the vaporizers from the vaporizer manifold.
- 3. Remove the upper rear panel (Section 4.3).
- 4. Using a 13mm open ended wrench, remove the vaporizer pressure relief valve by turning counterclockwise.




- 5. Verify that the parts are free of dust and dirt.
- 6. Install a new vaporizer pressure relief valve.
- 7. To reassemble, perform the previous steps in reverse order.
- 8. Perform the checkout procedure (Section 3).

4-30 07/04 1009-0356-000

4.10.5 Replace vaporizer manifold

- 1. Remove the upper rear panel (Section 4.3).
- 2. Remove the front flowmeter shield (Section 4.9.1).
- 3. Remove the right side panel (A).
- 4. From the front of the machine, remove the screw (**B**) at the right upright of the flowhead bezel.
- From the back of the machine, remove the two screws (C) that hold the vaporizer manifold vertical support to the flowhead bracket.
- From the back of the machine, remove the two screws (D) that hold the vertical support to the vaporizer manifold.
- 7. Remove the vertical support from the machine.
- 8. While holding the vaporizer manifold, remove the two screws (**E**) at the right-hand extrusion to release the manifold.




- Install the new vaporizer manifold in reverse order. Do not fully tighten the screws until they are all in place.
 - Attach the new manifold to the right-hand extrusion (E).
 - Attach the vertical support to the vaporizer manifold (D).
 - Attach the bottom of the vertical support to the flowhead bracket (C).
 - Attach the vertical support to the front bezel (B).
- 10. Tighten the mounting screws in the following order: E, D, C, B.
- 11. Reassemble the machine.
- 12. Perform the checkout procedure (Section 3).

4.11 Replace ACGO selector switch

Removal

- 1. Remove the tabletop (Section 4.4).
- 2. Clip the tie wraps (1) from the outlet barb fittings at the side of the switch.



- 3. Disconnect the fresh gas (2) and flush (3) tubes at the back of the switch.
- Disconnect the wires from the ACGO mode microswitch (4) at the back of the selector switch.
- 5. Disconnect the wires from the flush pressure switch (5) on top of the selector switch.
- 6. Set the ACGO selector switch to ABS.
- 7. Back out the selector switch mounting screws (6) until the tips are flush with the face of the mounting casting.
- While pushing the selector knob toward the machine and holding it steady, push the valve body toward the knob and rotate it counterclockwise to separate the valve body from the knob assembly.
- 9. Remove the knob assembly and protective shroud from the machine.
- 10. Remove the valve from the silicone output tubes.

Replacement

- Remove the knob assembly from the valve body.
- Back out the selector switch mounting screws until the tips are flush with the face of the mounting casting.
- 3. Guide the outlet fittings of the valve body into their respective silicone tubes.
- Hold the selector knob with the indicator mark facing down. Turn the chrome collar to its maximum counterclockwise position (as viewed from the front).





4-32 07/04 1009-0356-000

- Place the shroud over the knob and guide the assembly into the pan opening.
- Ensure that the indicators on the shroud align with label on the pan and the alignment tab mates with the alignment hole in the pan.
- While holding the knob assembly steady against the pan, place the valve assembly over the knob actuator. Using moderate force press the two assemblies together. The knob should rotate to the ACGO position.
- While continuing to force the assemblies together, rotate the knob assembly to the ABS position. The assemblies should snap into place.
- 9. Verify proper alignment of the knob with the setting indicators. Tighten the mounting screws evenly to secure the switch assembly to the pan.
- Secure the outlet tubing with tie wraps.
- Connect the fresh gas and flush gas tubing. Pull on the tubing to ensure that it is locked in the fitting.
- 12. Reconnect the wires to the ACGO mode microswitch at the back of the valve (top two terminals).
- Reconnect the wires to the flush pressure switch at the top of the valve (upper and lower terminals).
- 14. Replace the tabletop.

Test procedure


- 1. Confirm that flush flow and 5 L/min fresh gas flow are diverted to the ACGO port and the ABS in the respective knob positions.
- 2. Confirm that the ventilator display indicates ACGO mode when the valve is set in the ACGO position.
- Test the function of the flush pressure switch (Service Mode - "Display Discrete I/O Signals").
- 4. Perform the low-pressure leak test (Section 3.5).
- 5. Perform the checkout procedure (Section 3)







4.12 Clean or replace ACGO port flapper valve

- 1. Remove the tabletop (Section 4.4).
- Remove the ACGO cap mounting screws.
- 3. Remove the cap.
- Examine the flapper and disk for obstructions or debris. Clean with isopropyl alcohol if necessary; retest.
- 5. If leak persists, replace the flapper.
 - Remove the flapper from the check valve disk.
 - Clean the new flapper with isopropyl alcohol.
 - Apply a drop of isopropyl alcohol to the center hub of the new flapper.
 - Before the alcohol evaporates, align the center hub of the new flapper with the center hole of the check valve disc.
 - While pressing the flapper against the disc, use you fingernail to help pull the hub through the disc from the other side.







- 6. Lubricate the o-ring sparingly with Krytox (do not get Krytox on the flapper).
- 7. Insert the flapper assembly into the ACGO outlet with the flapper up.
- 8. Replace the cap.

4-34 07/04 1009-0356-000


4.13 Reconfigure sample gas return line

Sample gas return is directed to the scavenging system as a factory default. Perform the following to reroute the sample gas back to the breathing system. Refer to "Tubing" on page 9-8.

- 1. Remove the tabletop (Section 4.4).
- Port 4 (A) of the ABS breathing system is connected to the expiratory circuit, downstream of the expiratory check valve. As a factory default, Port 4 is plumbed with a length of tubing that is plugged (B) at the far end.
- 3. Remove the plug from the tube.
- 4. Find the sample return line at the left-rear corner of the pan assembly. The sample return line includes an inline connector (C) at the point where the sample line goes down into the vent engine housing.
- A
- Separate the scavenging tube, removing the inline connector from the portion of the tube that extends into the vent engine housing. Plug the open end of the scavenging tube with the plug removed above.
- 6. Insert the inline connector from the sample return port into the open tube to Port 4. Pull on the connector to ensure that it is securely connected.
- 7. Replace the tabletop.
- 8. Perform the checkout procedure (Section 3).

4.14 Replace the APL valve

- 1. Remove the ABS breathing system.
- The APL valve is held in place with a spring and a retainer (A) that snaps into a recess in the lower body of the APL valve. To release the retainer, place an appropriately sized straight blade screwdriver into the housing cutout (B). Twist the screwdriver to release the retainer.
- Place the new APL valve into position with the setting indicator facing forward.
- 4. Place the spring into the retainer.
- While holding the APL valve tight to the housing, Snap the spring and retainer onto the valve body from below.
- 6. Reinstall the ABS breathing system.
- 7. Perform the checkout procedure (Section 3).



4-36 07/04 1009-0356-000

4.15 Replace the bag support arm

- 1. Remove the ABS breathing system from the machine.
- 2. From the underside of the casting, remove the two screws/ lockwashers (A) that hold the arm in place.
 - If either of the pins (see below) remain in the casting, remove them from the casting.



- 3. Install the new bag support arm assembly.
 - Position the bag arm over mounting pattern of 4 small holes in the support casting. The arm should extend towards the front of the machine. Align the two pins (B) extending from the base of the bag arm assembly, with two of the small holes in the casting that are in line with the APL valve.



- Lower the bag arm, pushing the two pins into the holes.
- From the underside of the casting, secure the bag arm with two M3x16 screws and lockwashers.
- 4. Test the force required to swing the bag arm from side to side and adjust if necessary.

The force is adjusted by turning the lock nut (8-mm socket) which is accessible from underneath the support casting. Turn clockwise to increase the force and counterclockwise to reduce the

- force. Swing the bag arm sideways
- through the 90 degree arc permitted by its internal stop.
- Adjust to just enough friction to prevent the bag arm from swinging sideways as the bag height is being changed. The bag arm height is changed by squeezing the lock release lever (C) at the free end of the bag arm and rotating it to the desired position.
- 5. Replace the ABS breathing system.


Note: The adjustment nut is initially set so that 5-mm of exposed thread extends from the adjusting nut. With use, the force required to move the arm increases and may require readjustment.


4.15.1 Servicing the bag support arm

Service parts for the bag support arm include the upper and lower assemblies.


To replace either assembly:

- 1. Remove the bag support arm from the machine (Section 4.15).
- To separate the upper assembly from the lower assembly, use a small (2.5-mm) pin punch from the bottom to drive the dowel pin up and out.
- To assemble the bag arm, apply a light coat of Krytox to the area of the upper arm (A) that extends into the lower arm (including the dowel pin groove).





- 4. Insert the upper assembly into the lower assembly. Align the surface (**B**) of the upper assembly with the surface (**C**) of the lower assembly.
- Insert the dowel pin into the hole (from the top side as shown). Drive the dowel pin into the bag arm until it is flush with the top surface.




4-38 07/04 1009-0356-000


4.15.2 Replace friction pad in lower bag arm assembly


- 1. Remove the ABS breathing system from the machine.
- Using an 8-mm socket, remove the nut (A), shoulder washer, and spring from the lower assembly.
- Lift the bag support arm off of the swivel post. Remove the old friction pad.
- 4. Wipe any residue and friction particles from the post.
- Insert a new friction pad into the base. Keep approximately
 1 mm of space between the end of the pad and the bottom of the base.

Note: Align the friction pad gap with the seam (**B**) in the base. Position the retaining screw so the pin (**C**) at the base is perpendicular to the seam.

- With the bag support arm facing forward, place the base of the arm over the swivel post. Ensure that the slot in the base (D) engages the tab (E) on the swivel post.
- Replace the spring, shoulder washer and nut. Tighten the nut until 5 mm of thread extends beyond the nut.
- 8. Follow the procedure in Section 4.15 to adjust the force required to swing the bag arm from side to side.









4.15.3 Replace bag port housing

- Remove the bag support arm cover (A) – screw and lockwasher from below.
- Remove nut (B) to remove the release lever (C).
- 3. Remove the retaining ring (D).
- Slide the bag port housing (E) off the end of the bag support arm.
- Before installing the new bag port housing, clean and lubricate sparingly with Krytox the exposed metal end (F) and the guide slot (G) of the bag support arm.






- Slide the new bag port housing onto the bag arm. Secure it with the retaining ring.
- Lubricate sparingly with Krytox the pivot boss (H) before replacing the release lever.
- After replacing the release lever, adjust the mounting nut so that a 2-mm gap remains between the lever and housing when the release lever is fully depressed.
- 9. Replace the bag arm cover.



4-40 07/04 1009-0356-000

4.16 Replace auxiliary 02 flowmeter

- 1. Bleed all gas pressure from the machine (Section 4.2).
- 2. Ensure that all cylinder and pipeline gauges read zero before proceeding.
- 3. Remove the tabletop (Section 4.4).
- 4. Remove the adjustment knob from the flowmeter; pull forward.
- 5. Remove the gauge panel mounting screws and move the panel forward to access the flowmeter.




- 6. Disconnect the inlet tube fitting (A).
- 7. Disconnect the tube (B) from the outlet fitting.
- 8. Remove the four screws (C) that hold the flowmeter mounting bracket to the front panel.
- 9. Transfer the mounting bracket to the new flowmeter.
- 10. Reassemble in reverse order.
- 11. Perform the checkout procedure (Section 3).


4.17 Replace the suction control module

The suction control module can be replaced by removing the front panel, along with the ABS and the tabletop, to gain access. Alternatively, if the situation warrants, the suction control module can be accessed by removing the rear panel.

4.17.1 Front panel method

- 1. Bleed all gas pressure from the machine (Section 4.2).
- 2. Ensure that all cylinder and pipeline gauges read zero before proceeding.
- 3. Remove the tabletop (Section 4.4).
- 4. Remove the gauge panel mounting screws and move the panel forward to access the suction control module.







- 5. Disconnect the tubing from the suction control module.
 - Tube (A) from overflow safety trap (suction).
 - Tube (B) from vacuum source.
 - If venturi drive, tube (C) from pilot valve adapter.
- Remove the two mounting screws that hold the suction control module to the mounting bracket.
- 7. Transfer the mounting bracket to the new suction control module.
- 8. Reassemble in reverse order.
- 9. Perform the checkout procedure (Section 3).

4-42 07/04 1009-0356-000

4.17.2 Rear panel method

- 1. Lower the upper rear panel (Section 4.3).
- 2. Disconnect the white (A) vacuum and black (B) suction fittings from the rear panel. Do not remove the tubing from the regulator.
- 3. If you are replacing a venturi drive suction control module, disconnect the tube (**C**) from the pilot valve adapter.






- Remove the two screws (D) that hold the suction control module to the mounting bracket.
- 5. Remove the regulator assembly from the front panel.
- 6. Transfer the tubing to the new regulator:
 - Attach the vacuum source tube (white fitting) to the lower connector (E).
 - Attach the suction tube (black fitting) to the upper connector (F).
- 7. Guide the tubes into the front panel opening.
- 8. While holding the regulator assembly against the front panel, attach the retaining bracket to the regulator. Tighten the screws to secure the regulator assembly.
- 9. If applicable, attach the control port tube to the pilot valve adapter (C).
- 10. Attach the vacuum and suction fitting to the rear panel manifold.
- 11. Replace the rear panel.
- 12. Perform the checkout procedure (Section 3).


4.18 Replace ABS breathing system components

4.18.1 Replace Bag/Vent switch assembly

- 1. Remove the ABS breathing system.
- 2. From the underside, remove the bellows base manifold (A) and fully loosen the two captive screws (B) at the bag port side of the APL/BTV manifold.



3. From the topside, rotate the Bag/Vent switch cartridge counterclockwise until the Bag/vent switch outlet port (**C**) clears the bellows housing.




- 4. Lift out the Bag/Vent switch cartridge from the housing.
- 5. Replace the Bag/Vent switch cartridge in reverse order.
- 6. Reinstall the ABS breathing system.
- 7. Perform the checkout procedure (Section 3).

4-44 07/04 1009-0356-000


4.18.2 Replace bellows base latch assembly


To replace the latch assembly, you must disassemble the bellows base assembly to the point where you can remove the guide (A) and latch assembly (B) as a unit.

- 1. Remove the Bag/Vent switch cartridge (Section 4.18.1).
- 2. Remove the two remaining screws (**C**) that hold the APL/BTV manifold to the bellows base assembly. Remove the APL/BTV manifol.



- To remove the guide/latch assembly, remove two mounting screws (D) from the underside. Remove two additional mounting screws from the topside. Remove the guide/latch assembly from the bellows base assembly.
- 4. Separate the latch assembly from the guide assembly.
- 5. To install the new latch assembly, put the spring (**E**) into place in the guide assembly (long leg down).
- 6. Place the latch assembly on the guide assembly so that the latch engages the short leg of the spring. Secure the latch assembly (**F**) to the guide assembly.





- 7. Mount the guide/latch assembly into the bellows base assembly.
- 8. Reassemble the breathing system in reverse order.
- 9. Perform the checkout procedure (Section 3).

4.19 Replace casters

⚠ WARNING

Replacing a caster requires at least two people to maneuver and tip the machine. Personal injury and/or machine damage is possible if one person attempts this procedure alone.

 Disconnect all pipeline hoses from the wall and the machine, close all gas cylinders, unplug the power cord, and set the system switch to standby.

△ CAUTION

Remove the vaporizers before tipping the machine. If a vaporizer is inverted, it must be set to 5% and purged for 30 minutes with a 5 L/min flow. The interlock system prevents purging more than one vaporizer at a time.

2. Remove the absorber, the vaporizers, gas cylinders, drawers and all auxiliary equipment.

⚠ CAUTION

To prevent damage, do not tip the Aespire machine more than 10 degrees from vertical.

3. Block the opposite wheels; then, block up the machine until there is enough room to remove the defective caster.

To block up the machine, tip and slide blocks under the caster base. Raise both sides evenly until the unit is high enough to remove the caster.


- The casters are threaded into the base and held with a Loctite compound.
 Remove the caster with an appropriately sized open-end wrench.
- If required, clean the treads of the new caster with denatured alcohol.
- Apply Loctite 242 to the threads of the new caster. Install the caster securely into place.
- Make sure the caster turns freely.
- 8. Carefully lower the machine to the floor.
- 9. Perform the checkout procedure (Section 3).



4-46 07/04 1009-0356-000

4.20 Replace task light and switch

Remove the four screws (A) that hold the task-light lens to the upper shelf.



4.20.1 To replace the task-light switch

- 1. Using a small needle-nose pliers, disconnect the switch harness from the task-light circuit board connector (B).
- 2. Remove the two screws (C) that hold the switch retainer plate to the upper shelf.
- 3. Transfer the switch retainer plate to the new switch, counter-sunk side to the outside.
- 4. Mount the switch to the upper shelf.
- Remount the task-light assembly. Ensure that the switch harness and the task-light harness wires are positioned in their respective recesses and are not pinched under the task-light lens.

4.20.2 To replace the task-light circuit board

- 1. Using a small needle-nose pliers,
 - disconnect the switch harness from the task light circuit board connector (B).
 - disconnect the task-light harness from the task light circuit board connector (D).
- 2. Slide the task-light circuit board out of the lens.
- Slide the new task light into the lens, ensuring that the connectors are aligned with the lens cutouts.
- Plug the task-light harness and the switch harness into their respective connectors on the task-light circuit board. Use a small screwdriver to push the connectors securely into place.
- Remount the task-light assembly. Ensure that the switch harness and the task-light harness wires are positioned in their respective recesses and are not pinched under the task-light lens.

4.21 Replace the display arm or display cables

Cable replacement requires that you first remove the display arm from the dovetail extrusion.

Before replacing the display arm, note the routing of the cables.

After replacing the display arm, ensure that the cables are dressed properly and do not interfere with the motion of the display arm.

Follow the procedure in Section 4.21.1 for the recommended use of cable ties.

4.21.1 Cable tie installation

 Wrap the cable tie around the cables in the indicated position. Start the tail through the clamp as shown.



2. Feed the tail between the cables.



Pull the tie tight and cut off the tail below the outer surface of the cables. Do not leave a sharp edge exposed.




4-48 07/04 1009-0356-000

4.21.2 Removing the display arm

If equipped, remove additional equipment from the arm before removing the arm.


- 1. Disconnect the cables from the display.
- Remove the display from the display arm.
- 3. Remove the cables from the cable clamps.
- 4. Loosen the screws that secure the display arm in the dovetail.
- If required, use a rubber mallet to tap the display arm out of the dovetail.



4.21.3 Replacing a display cable

- Remove the three screws that hold the dovetail extrusion to the upright. Remove the extrusion to allow cable replacement.
- 2. As require to access the particular cable routing for replacement, remove either (or all):
 - the rear cover (Section 4.3),
 - the tabletop (Section 4.4),
 - or the AC Inlet module.
- After replacing the cable, place the cables in the notches in the order shown.
 - A Serial Isolation
 - · B Monitoring
 - · C Vent Engine
 - D Power cable
- 4. Install the dovetail extrusion loosely to hold the cables in place.
- Adjust the cable length outside the machine to approximately 66 cm.
- 6. Securely tighten the extrusion mounting screws.






4-49

1009-0356-000 07/04


4.21.4 Installing the long arm

- 1. Place the arm into the extrusion.
- Use a rubber mallet to tap the arm into place. Leave a 12-mm gap between the lower edge of the arm mounting plate and the end of the dovetail.
- 3. Tighten the mounting screws to secure the display arm.
- 4. Remount the display.
- 5. Route the display cables neatly through the cable clamps.
- 6. Attach the cables to the display.
- 7. If required, install cable tie in the locations shown.
 - Place one cable tie close to the pivot of the arm.
 - Place the second cable tie near the display as shown.
- Ensure that the cables are secured so that they do not interfere with the display arm through the entire range of motion.












4-50 07/04 1009-0356-000

4.21.5 Installing the short arm

- 1. Place the arm into the extrusion.
- Use a rubber mallet to tap the arm into place. Leave a 12-mm gap between the lower edge of the arm mounting plate and the end of the dovetail.
- 3. Tighten the mounting screws to secure the display arm.
- 4. Remount the display.
- 5. Route the display cables neatly through the cable clamps.
- 6. Attach the cables to the display.
- 7. If required, install cable tie in the locations shown.
 - Place one cable tie close to the pivot of the arm.
 - Place the second cable between the two straps.
 - Place the third cable tie near the display as shown.
- Ensure that the cables are secured so that they do not interfere with display through the entire range of motion.









4.22 Replace display and cables in ProTIVA machine


- Remove two thumbscrews that secure the display mounting bracket to the no-vap manifold.
- 2. Lift the display slightly to disengage the mounting pins.
- Lower the display face down on the worksurface.
- Replace the display and reassemble in reverse order.
- To replace a cable, remove the retaining block and feed the cable to the back of the machine.
- As require to access the particular cable routing for replacement, remove either (or all):
 - the rear cover (Section 4.3),
 - the tabletop (Section 4.4),
 - · or the AC Inlet module.
- 7. Attach the cable to the display.
- 8. When replacing the cable retainer, ensure:
 - that the cables fit properly in their feed-through slots (as detailed below),
 - · and that they extend out of the machine with minimal slack.
- 9. From the back of the machine, ensure that the cables are routed such that they do not interfere with the replacement of the rear cover.





- A Serial Isolation
- B Monitoring board
- C Vent Engine
- D Power cable





4-52 07/04 1009-0356-000

5 Maintenance

In this section

This section covers the regular maintenance procedures (minimum requirements) needed to make sure that the Aespire Anesthesia Machine — including the ventilator — operates to specifications.

| 5.1 Aespire Planned Maintenance | 5-2 |
|--|-----|
| 5.1.1 Every twelve (12) months | 5-2 |
| 5.1.2 Every twenty-four (24) months | 5-3 |
| 5.2 Auxiliary 0 ₂ flowmeter tests | 5-4 |
| 5.3 Integrated Suction Regulator tests | 5-5 |

⚠ WARNINGS

Do not perform testing or maintenance on the Aespire Anesthesia Machine while it is being used on a patient. Possible injury can result.

Items can be contaminated due to infectious patients. Wear sterile rubber gloves. Contamination can spread to you and others.

Obey infection control and safety procedures. Used equipment may contain blood and body fluids.

5.1 Aespire Planned Maintenance

| Serial Number: | Date: (YY/MM/DD) / / | |
|-----------------------------------|---|--|
| Hospital: | Performed by: | |
| □ 12 months □ 24 month | Δ | |
| 5.1.1 Every twelve (12) months | Perform the following steps every 12 months. | |
| Machine Parts Replacement | Refer to the listed section in this manual. Perform the following step: | |
| | Replace the vaporizer port o-rings (Section 4.10.1) (Kit Stock Number 1102-3016-000) | |
| Machine Checks and Tests | Refer to the Aespire User's Reference Manual, Part 2. Perform the following steps: | |
| | User maintenance listed below. Including disassembly, inspection, cleaning and parts replacement as required (Section 3 and Section 2) | |
| | AGSS Maintenance: Empty any condensate from the reservoir (disposable item). Inspect air brake for occlusion. Inspect, clean or replace filter on active AGSS. Breathing Circuit Maintenance Bellows Assembly Maintenance Bellows Assembly Tests O₂ Sensor Calibration Flow Sensor Calibration | |

Refer to listed sections in this manual.

Perform the following steps:

| 2. | Inspect the system (Section 3.1) |
|----|--|
| 3. | Pipeline and cylinder tests (Section 3.2) |
| 4. | Flow control and pressure relief tests (Section 3.3) |
| 5. | Vaporizer back pressure test (Section 3.4) |
| 6. | Low-pressure leak test (Section 3.5) |
| 7. | Airway pressure gauge accuracy check (Section 6.6.2) |
| 8. | Alarm tests (Section 3.6) |
| 9. | Breathing systems tests (Section 3.7) |
| 10 | . Auxiliary ${\rm O}_2$ flowmeter tests, if equipped with option (Section 5.2) |
| | |

5-2 07/04 1009-0356-000

| | 11. Integrated suction regulator tests, if equipped with option (Section 5.3) |
|--|---|
| | 12. Power failure test (Section 3.10) |
| | 13. Electrical safety tests (Section 3.11) |
| 7100 Ventilator Checks, Tests and Calibrations | Refer to the listed sections in the Aestiva 7100 Ventilator Service Manual. Perform the following steps: |
| | 1. MOPV pressure relief valve test (Section 6.3) |
| | 2. From the Ventilator Service Mode menu, perform the following: |
| | Display Discrete I/O Signals. Verify proper operation of all switches. (Section 4.10.2) |
| | Display Error Log. If any error codes have been logged follow the appropriate troubleshooting procedures. Clear the error log. (Section 4.5) |
| | Adjust Drive Gas Regulator (Section 4.9.3) |
| | Airway Sensor Span (Section 4.9.4) |
| | PEEP Valve Calibration (Section 4.9.5) |
| | Inspiratory Valve Calibration (Section 4.9.6) |
| | Pressure Sensitivity (Section 4.9.7) |
| 5.1.2 Every twenty-four (24) months | In addition to the 12-month requirements, replace the following parts every 24 months. All machine and ventilator parts should be replaced before performing the checks, tests, and calibrations. |
| 7100 Ventilator Parts Replacement | Refer to the listed sections in the Aestiva 7100 Ventilator Service Manual. Perform the following steps: |
| | Replace the internal backup battery (Section 7.2.3) (Stock Number 1504-3505-000). |
| | 2. Replace the free breathing flapper valve (Section 6.2) (Stock Number 0211-1454-100). |
| | 3. Replace the free breathing valve o-ring (Section 6.2) (Stock Number 1503-3208-000). |

5.2 Auxiliary 02 flowmeter tests

- 1. Open the O₂ cylinder valve or connect an O₂ pipeline.
- Rotate the flow control clockwise (decrease) to shut off the flow. The ball should rest at the bottom of the flow tube and not move.
- Rotate the flow control counterclockwise (increase). The ball should rise immediately after rotation is begun. It should rise smoothly and steadily with continued counterclockwise rotation. When a desired flow is set, the ball should maintain in a steady position.
- 4. Rotate the flow control clockwise to shut off the flow.

Flow Accuracy Test

Note: To check flow accuracy, be sure that the flow test device is capable of measuring 0 to 15 L/min with an accuracy of $\pm 2\%$ of reading.

- 1. Connect the flowmeter outlet to the flow test device.
- Adjust the flowmeter so the center of the ball aligns with the selected test point (observe that the ball maintains a steady position for 10 seconds).
- The test device reading should be between the limits shown for each of the selected settings in the table below.

Flow Tester Reading

| Flowmeter Setting L/min | Lower Limit L/min | Upper Limit L/min |
|-------------------------------|----------------------|----------------------|
| 1 | 0.52 | 1.48 |
| 3 | 2.56 | 3.44 |
| 5 | 4.60 | 5.40 |
| 10 | 9.70 | 10.30 |
| maximum (valve fully open) | 12.00 | |

- 4. Rotate the flow control clockwise to shut off the flow.
- 5. Close the O₂ cylinder valve or disconnect the O₂ pipeline.

5-4 07/04 1009-0356-000

5.3 Integrated Suction Regulator tests

Note

There are two types of integrated suction systems for the Avance anesthesia machine:

- Continuous Vacuum Regulator, Three-Mode, Pipeline Vacuum
- Continuous Vacuum Regulator, Three-Mode, Venturi Derived Vacuum

For Pipeline Vacuum systems,

a vacuum source of at least 500 mm Hg (67 kPa or 20 in Hg) is required for testing. The supply open flow must be a minimum of 50 L/min.

For Venturi Derived Vacuum systems,

an O2 or Air source of at least 282 kPa (41 psi) is required for testing.

Gauge Accuracy

The gauge needle should come to rest within the zero range bracket when no suction is being supplied. Gauges which do not comply may be out of tolerance.

Note

To check gauge accuracy, be sure that the test gauge is capable of measuring 0 to 550 mm Hg with an accuracy of $\pm 1\%$ of reading.

- 1. Connect the suction patient port to the test gauge.
- 2. Turn the mode selector switch to I (ON).
- Ensure that the vacuum test gauge is in agreement with the suction vacuum gauge ± 38 mm Hg/5 kPa at the following test points.

| Test points | |
|----------------------|--------------------------------|
| Suction vacuum gauge | Test gauge tolerance |
| 100 mm Hg (13.3 kPa) | 62-138 mm Hg (8.3-18.4 kPa) |
| 300 mm Hg (40 kPa) | 262-338 mm Hg (35-45 kPa) |
| 500 mm Hg (66.7 kPa) | 462-538 mm Hg (61.6-71.7) kPa) |

Flow Test

Note: To check flow accuracy, be sure that the flow test device is capable of measuring 0-30 L/min.

- 1. Connect the patient port of the suction regulator to the flow test device.
- 2. Rotate the suction control knob fully clockwise (increase).
- 3. Turn the mode selector switch to I (ON) and verify that the flow rate is:
 - at least 20 L/min.
- 4. Disconnect the test flowmeter.

(Tests continue on next page.)

Regulation Test

- Turn the mode selector switch to I (ON).
- 2. Occlude the patient port of the suction regulator.
- 3. Set the vacuum regulator gauge to 100 mm Hg/13 kPa.
- 4. Open and close the patient port several times.
- 5. With the patient port occluded, the gauge should return to 100 mm Hg/13 kPa within a tolerance of $\pm 10 \text{ mm Hg}/1.3 \text{ kPa}$.

Vacuum Bleed Test

- 1. Occlude the patient port of the suction regulator.
- 2. Set the vacuum regulator gauge to 100 mm Hg/13 kPa.
- 3. Turn the mode selector switch to 0 (OFF) and observe the gauge needle. It must return to the zero range bracket or stop pin within 10 seconds.

Vacuum Leak Test

- 1. Turn the mode selector switch to 0 (OFF).
- Rotate the suction control knob a minimum of two full turns in the clockwise direction (increase suction) to ensure its setting is not at the off position.
- 3. Occlude the patient port of the suction regulator.
- 4. Observe the suction gauge, the needle should not move.
- 5. Rotate the suction control knob fully counterclockwise to ensure its setting is at the fully off position.
- 6. Turn the mode selector switch to I (ON).
- 7. Observe the suction gauge, the needle should not move.

5-6 07/04 1009-0356-000

6 Calibration

| After adjustments and calibration are completed, always perform the checked procedure. Refer to Section 3 of this manual. | out |
|---|--|
| This section covers calibration procedures for components of the Aespire anesthesi machine. | a |
| 6.1 Primary Regulators | .6-2 |
| 6.1.1 Test setup | .6-2 |
| 6.1.2 Testing Primary Regulators | .6-3 |
| 6.1.3 Adjusting Primary Regulators | .6-6 |
| 6.2 Secondary Regulators | .6-7 |
| 6.2.1 Testing/Adjusting Secondary Regulators or Balance Regulators | .6-7 |
| 6.3 Flowmeter Needle Valve Calibration | .6-8 |
| 6.3.1 O ₂ Needle Valve Calibration (Minimum Flow) | 6-8 |
| 6.3.2 N ₂ O Needle Valve Calibration (Minimum Flow) | 3-10 |
| 6.3.3 Air Needle Valve Calibration (Minimum Flow) | 3-14 |
| 6.3.4 Needle Valve Calibration (Maximum Flow) | 5-17 |
| 6.4 Link system calibration | 3-18 |
| 6.5 O ₂ Flush Regulator | 3-23 |
| 6.6 Airway pressure gauge | 3-24 |
| 6.6.1 Zero the pressure gauge | 3-24 |
| 6.6.2 Checking the pressure gauge accuracy 6 | 3-25 |
| | This section covers calibration procedures for components of the Aespire anesthesi |

6.1 Primary Regulators

Follow the procedure in Section 6.1.1 to gain access to the regulators. Then, in Section 6.1.2, select the test that is appropriate for the regulator you are testing.

⚠ WARNING

When testing/adjusting N₂O regulators, nitrous oxide flows through the system. Use a safe and approved procedure to collect and remove it.

6.1.1 Test setup


⚠ WARNING

Wear safety glasses while test device is connected to the test port.

⚠ CAUTION

Be careful not to plug the output of the primary regulator without having a pressure relief valve in the output circuit.

- 1. Set the system switch to Standby.
- 2. Disconnect all pipeline supplies.
- Remove the upper rear panel (Section 4.3).For 3rd gas cylinder regulators, also remove the lower rear panel.
- 4. If equipped, turn the auxiliary O₂ flowmeter control fully clockwise (no flow).
- 5. Install a full cylinder in the cylinder supply to be tested. It is essential that the cylinder be within 10% of its full pressure.
- 6. Remove the plug from the test port and connect a test device capable of measuring 689 kPa (100 psi).



6-2 07/04 1009-0356-000

6.1.2 Testing Primary Regulators

There are two variations of the test procedure for the primary regulators:

- Test A For primary regulators that supply drive gas to the ventilator.
- Test B − For all gases not used to supply drive gas to the ventilator.

Test A For primary regulators that supply drive gas to the ventilator (02 or Air):

Under low flow conditions, the output pressure of a properly adjusted/functioning regulator should fall within specifications listed in step 4. Under high flow conditions, the output pressure should not drop below the specifications in step 12.

- 1. Remove the bellows assembly.
- 2. Slowly open the cylinder valve.
- 3. Set the system switch to On.
- 4. Low Flow Test: Set the fresh gas flow to 0.05 L/min (or minimum flow for O₂). When checking an Air regulator on systems that have a single flowtube, open the needle valve 1/8 turn from the minimum stop to achieve a flow close to 0.05 L/min.
 - Close the cylinder valve and allow the pressure to decay to 2068 kPa (300 psi) as indicated on the cylinder gauge (upper limit of the red band).
 The flow may be temporarily increased to facilitate the decay.
 - At the time that the cylinder pressure reaches 2068 kPa (300 psi), set the system switch to Standby.
 - Within one minute, the test device must stabilize between:
 - (60) DIN 372-400 kPa (54-58 psi)
 - (50) Pin Indexed 310-341 kPa (45.0-49.5 psi).
 - If the test device pressure does not stabilize within one minute, replace the cylinder supply.
 - If the test device stabilizes within one minute, but the readings are not within specifications, readjust the regulator (Section 6.1.3).
- 5. Slowly open the cylinder valve.
- Enter the Service Mode: (Push and hold the adjustment knob on the ventilator's display and set the system switch to On.)
- Select and confirm "Service Modes."
- Follow the menu structure outline below to reach the adjustment for the inspiratory flow valve. Select and confirm at each step.
 - "Diagnostics Tests/Tools"
 - · "Valves Test Tool"
 - · "Set Inspiratory Valve"

- High Flow Test: Rotate adjustment knob counterclockwise to obtain 65 (L/min):
 - · While watching the test device press confirm.
 - After 2 seconds, select "Go to Diagnostic Tests/Tools Menu" and press confirm to stop the gas flow.
 - The minimum test device reading observed must be greater than:

(60) DIN 221 kPa (32 psi)

(50) Pin Indexed 207 kPa (30 psi)

Repeat this step three times.

If the test device reading under "high flow" conditions is less than specified, readjust the regulator per the procedure in Section 6.1.3; however, set the regulated pressure higher by the difference you noted in this step plus 7 kPa (1 psi). This adjusts the "low flow" regulated output to the high side of the specification so that the "high flow" regulated pressure can fall within the specification.

If the regulator subsequently fails the "low flow" specification (step 4) because the reading is too high, replace the cylinder supply.

- Set the system switch to Standby.
- 11. Close the cylinder valve.
- 12. Bleed the system of all pressure.
- 13. Disconnect the test device and plug the test port (pull on the plug to ensure it is locked in the fitting).
- 14. Replace the bellows assembly.
- 15. Replace the rear panel(s).
- 16. Perform the checkout procedure (Section 3).

Test B For all gases not used to supply drive gas to the ventilator:

Under low flow conditions, the output pressure of a properly adjusted/functioning regulator should fall within specifications listed in step 4. Under high flow conditions, the output pressure should not drop below the specifications in step 7.

- If the cylinder supply being tested is N₂O, connect a source of O₂ and set the O₂ flow control to the minimum stop (pilot pressure for secondary regulator).
- 2. Slowly open the cylinder valve for the regulator being tested.
- 3. Set the system switch to On.

6-4 07/04 1009-0356-000

- 4. Low Flow Test: Set the flow of the gas being tested to 0.05 L/min (or minimum flow for O₂). When checking a regulator on systems that have a single flowtube, open the needle valve 1/8 turn from the minimum stop to achieve a flow close to 0.05 L/min.
 - Close the cylinder valve and allow the pressure to decay to 2068 kPa (300 psi) as indicated on the cylinder gauge (upper limit of the red band).
 The flow may be temporarily increased to facilitate the decay.
 - At the time that the cylinder pressure reaches 2068 kPa (300 psi), set the system switch to Standby.
 - Within one minute, the test device must stabilize between:
 - (60) DIN 372-400 kPa (54-58 psi)
 - (50) Pin Indexed 310-341 kPa (45.0-49.5 psi).
 - If the test device pressure does not stabilize within one minute, replace the cylinder supply.
 - If the test device stabilizes within one minute, but the readings are not within specifications, readjust the regulator (Section 6.1.3).
- 5. Slowly open the cylinder valve.
- 6. Set the system switch to On.
- High Flow Test: Set the flow control valve to the maximum indicated flow on the flow tube.
 - The test device reading must be greater than:
 - (60) DIN 221 kPa (32 psi)
 - (50) Pin Indexed 221 kPa (32 psi)
 - If the test device reading under "high flow" conditions is less than specified, readjust the regulator per the procedure in Section 6.1.3; however, set the regulated pressure higher by the difference you noted in this step plus 7 kPa (1 psi). This adjusts the "low flow" regulated output to the high side of the specification so that the "high flow" regulated pressure can fall within the specification.
 - If the regulator subsequently fails the "low flow" specification (step 4) because the reading is too high, replace the cylinder supply.
- 8. Set the system switch to Standby.
- 9. Close the cylinder valve.
- 10. Bleed the system of all pressure.
- 11. Disconnect the test device and plug the test port (pull on the plug to ensure it is locked in the fitting).
- Replace the rear panel(s).
- 13. Perform the checkout procedure (Section 3).

6.1.3 Adjusting Primary Regulators

Important: Cylinder supplies in an Aespire machine must have all primary regulators set to the same pressure range: (50) Pin Indexed or (60) DIN. If a regulator is replaced, the replacement regulator must be set (as required) to the same specification as the one removed.

Important: Install a full cylinder in the cylinder supply to be adjusted. It is essential that the cylinder be within 10% of its full pressure.

If the cylinder supply being adjusted is N_2O , connect a source of O_2 and set the O_2 flow control to the minimum stop (pilot pressure for secondary regulator).

To adjust the primary regulators, follow the procedure in Section 6.1.1 to gain access to the regulators.

Do not attempt to adjust without flow.

- 1. Slowly open the cylinder valve.
- 2. Set the system switch to On.
- 3. Set and maintain the fresh gas flow of the gas being tested to $0.05 \, \text{L/min}$ (or minimum flow for O_2). When adjusting a regulator on systems that have a single flowtube, open the needle valve 1/8 turn from the minimum stop to achieve a flow close to $0.05 \, \text{L/min}$.
- Close the cylinder valve and allow the pressure to decay to 2068 kPa (300 psi) as indicated on the cylinder gauge (upper limit of the red band). The flow may be temporarily increased to facilitate the decay.
- 5. When the cylinder gauge reaches the upper limit of the red band, adjust the regulator output pressure to:
 - (60) DIN 386-400 kPa (56-58 psi)
 - (50) Pin Indexed 327-341 kPa (47.5-49.5 psi).

Note: It may be necessary to open the cylinder valve and repeat steps 4 and 5 a number of times to achieve the above setting.

- 6. Test the regulator settings per the appropriate test in Section 6.1.2:
 - Test A For primary regulators that supply drive gas to the ventilator.
 - Test B For all gases not used to supply drive gas to the ventilator.



Adjust clockwise to increase setting

6-6 07/04 1009-0356-000


6.2 Secondary Regulators

⚠ WARNING

When testing N_2O regulators, nitrous oxide flows through the system. Use a safe and approved procedure to collect and remove it.

6.2.1 Testing/Adjusting Secondary Regulators or Balance Regulators

- 1. Set the system switch to Standby.
- 2. Remove the flowmeter panel shield (Section 4.9.1).
- 3. Remove the plug from the test port and connect a test device capable of measuring 689 kPa (100 psi) using 1/8-inch nylon tubing.



- 4. Set the flow of the tested gas and of O2 as detailed in the chart.
- 5. Verify that the output of the tested regulator is within the range listed in the chart.

| Regulator | Output | Flow Regulated gas | Flow O ₂ |
|------------------|--|-----------------------|------------------------|
| 02 | 207 ±7 kPa (30 ±1 psi) | 2 L/min | |
| Air | 207 ±7 kPa (30 ±1 psi) | 2 L/min | |
| N ₂ 0 | \pm 14 kPa (\pm 2 psi) of O_2 reading | 10 L/min | 4 L/min |

6. If required, adjust the $\rm O_2$ and Air regulators to meet the above specifications. The $\rm N_2O$ regulator is not adjustable; replace if out of range.

Note: The adjustment screws for these regulators are self-locking.

- 7. Disconnect the test device and plug the test port (pull on the plug to ensure it is locked in the fitting).
- 8. Perform the Flow Control Tests (Section 3.3).

6.3 Flowmeter Needle Valve Calibration

You need to calibrate a needle valve:

- if you install a new one,
- . if minimum and maximum flows are not within specifications.

6.3.1 O₂ Needle Valve Calibration (Minimum Flow)


⚠ CAUTION:

Do not force the needle valve against the seat. Overtightening the valve can cause the minimum flow setting to drift out of specifications.


- 1. Set the system switch to Standby.
- 2. Remove the flowmeter panel shield (Section 4.9.1).

Note: If adjusting an existing needle valve,

- remove the N₂O and O₂ knob and sprocket assemblies (on machines with O₂ only flowhead, remove the O₂ knob),
- and loosen the O₂ stop collar setscrews.
- 3. Slide a stop collar onto the valve stem with the stop tab toward the valve. Do not tighten setscrews.




- 4. Set the system switch to On.
- 5. Adjust the O2 needle valve to maintain a flow of
 - 50 ±25 mL/min for dual tube flowmeters.
 - 200 ±25 mL.min for single tube flowmeters.
- 6. Push the stop collar against the valve body.




6-8 07/04 1009-0356-000


7. Turn the collar clockwise until the collar stop tab contacts the minimum stop tab on the valve body. *Do not turn the valve stem*.



8. Carefully pull the collar back so there is a slight gap between collar and the valve body (but still engages the valve stop).



- 9. Tighten the collar setscrews. Start with the one opposite the tab if possible.
- Turn the valve stem counterclockwise at least one revolution to make sure the collar tab clears the valve stop.



If the stop does not clear:

- a. Turn the valve stem back to minimum position.
- b. Loosen the collar setscrews.
- c. Repeat steps 6 through 9.
- 11. Turn the valve stem clockwise to the minimum stop.
- 12. Verify that the flow is within the
 - 50 ±25 mL/min range for dual tube flowmeters.
 - 200 ±25 mL/min range for single tube flowmeters.
- 13. Set the maximum stop collar if necessary (Section 6.3.4).

Note: Maximum stop collars are required in Canada for all gas flow controls.

14. Calibrate the Link proportioning system (Section 6.4) (on machines with O₂ only flowhead, replace and set the knob so that at minimum flow the label text is horizontal; replace the flowmeter shield).

6.3.2 N₂O Needle Valve Calibration (Minimum Flow)

⚠ WARNING

You must be in a well ventilated room or use a gas evacuation device at this time. Anesthetic vapors exhausted into the room air can be harmful to your health.


△ CAUTION:

Do not force the needle valve against the seat. Overtightening the valve can cause the minimum flow setting to drift out of specifications.


- 1. Disconnect all pipeline supplies and close all cylinder valves.
- 2. Remove the upper rear panel.

Note: If adjusting an existing needle valve,

- remove the N₂O and O₂ knob and sprocket assemblies,
- and loosen the O₂ stop collar setscrews.
- 3. Remove the flowmeter panel shield (Section 4.9.1).
- 4. Disconnect the tube from the pilot port on the N₂O regulator.




5. Disconnect the 4-mm outlet tube from the back of the $\rm N_2O$ pipeline manifold.



6. Using a 4-mm tube/tee fixture (see Service Tools - Section 8.1.3), connect a tube (tee end) to the N₂O supply outlet at the back of the pipeline manifold. Connect the outlet tube to the open connection on the tee connector of the fixture. Connect the other end of the fixture to the N₂O pilot port at the front of the flowhead. This setup supplies pilot pressure to the N₂O balancing regulator during the minimum stop calibration.

6-10 07/04 1009-0356-000


7. Slide a stop collar onto the valve stem with the stop tab toward the valve. Do not tighten setscrews.




- 8. Connect either an N₂O pipeline or cylinder supply.
- 9. Slowly open the N2O cylinder valve.

Important: Do not connect the O_2 pipeline or open the O_2 cylinder valve.


10. Adjust the needle valve until the float is nearly touching the filter, but not quite.




- 11. Disconnect the tubing from the inlet of the vaporizer manifold (closest to flowhead).
- 12. If the machine has an Air option, bleed down the air supply. Air can inflate the bubble (next step) if it is not shut off.
- 13. Apply a small amount of leak detection fluid (Snoop) to the end of the tube to form a bubble.
- 14. Turn the valve stem clockwise until the bubble no longer inflates. Do not turn more than 10 degrees clockwise past this point.




15. Push the stop collar against the valve body.




16. Turn the collar clockwise until the collar stop tab contacts the minimum stop tab on the valve body. *Do not turn the valve stem*.



17. Carefully pull the collar back so there is a slight gap between collar and the valve body (but still engages the valve stop).



- Tighten the collar setscrews. Start with the one opposite the tab if possible.
- 19. Turn the valve stem counterclockwise at least one revolution to make sure the collar tab clears the valve stop.



If the stops do not clear:


- a. Turn the valve stem back to minimum position.
- b. Loosen the collar setscrews.
- c. Repeat steps 14 through 17.

6-12

- 20. Turn the valve clockwise to the minimum stop.
- 21. Verify there is no flow at the end of the tube.
- 22. Thoroughly clean the end of the nylon tube and reconnect it to the vaporizer manifold inlet.
- 23. Set the maximum stop collar if necessary (Section 6.3.4).

Note: Maximum stop collars are required in Canada for all gas needle valves.

- 24. After calibrating minimum flow for N₂0:
 - a. Close the cylinder valve and use the needle valve to bleed the remaining gas.
 - Remove the test fixture connecting the N₂O gas supply to pressure balance regulator pilot port.
 - c. Reconnect the pilot tube to the N₂O pilot port. Pull on the tubing to ensure it is locked into the fitting.



- d. Reconnect the outlet tube to the N_2O pipeline supply manifold.
- 25. Calibrate the Link proportioning system (Section 6.4).

6.3.3 Air Needle Valve Calibration (Minimum Flow)


⚠ CAUTION:

Do not force the needle valve against the seat. Overtightening the valve can cause the minimum flow setting to drift out of specifications.


- 1. Set the system switch to Standby.
- 2. Disconnect all pipeline hoses and close all cylinder valves except for air.
- 3. Remove the flowmeter panel shield (Section 4.9.1).
- 4. Remove the upper rear panel.

Note: If adjusting an existing needle valve,

- · remove the Air knob,
- and loosen the Air stop collar setscrews.
- Slide a stop collar onto the valve stem with the stop tab toward the valve.Do not tighten setscrews.




- 6. Set the system switch to On.
- 7. Adjust the needle valve until the float is nearly touching the filter, but not quite.




- Disconnect the tubing from the inlet to the vaporizer manifold (left end of manifold).
- 9. Apply a small amount of leak detection fluid (Snoop) to the end of the tube to form a bubble.

6-14 07/04 1009-0356-000


10. Turn the needle valve clockwise until the bubble no longer inflates. Do not turn more than 10 degrees clockwise past this point.




11. Push the stop collar against the valve body.




12. Turn the collar clockwise until the collar stop tab contacts the minimum stop tab on the valve body. **Do not turn the valve stem**.



13. Carefully pull the collar back so there is a slight gap between collar and the valve body (but still engages the valve stop).



- Tighten the collar setscrews. Start with the one opposite the tab if possible.
- 15. Turn the valve stem counterclockwise at least one revolution to make sure the collar tab clears the valve stop.



If stops do not clear:

- a. Turn the valve stem back to minimum position.
- b. Loosen the collar setscrews.
- c. Repeat steps 11 through 14.
- 16. Turn the valve stem clockwise to the minimum stop.
- 17. Verify there is no flow at the end of the tube.
- 18. Thoroughly clean the end of the nylon tube and reconnect it to the vaporizer manifold inlet.
- 19. Set the maximum stop collar if necessary (Section 6.3.4).

Note: Maximum stop collars are required in Canada for all gas needle valves.


- 20. Set the knob so that at minimum flow the label text is horizontal and the knob is on an even plane (front to back) with the N_2O and O_2 knobs.
- 21. Replace the flowmeter panel shield and the rear panel.

6-16 07/04 1009-0356-000

6.3.4 Needle Valve Calibration (Maximum Flow)

Note: Maximum stop collars are required in Canada for all gas needle valves.

- 1. Calibrate the needle valve for minimum flow:
 - Section 6.3.1 for 0₂
 - Section 6.3.2 for N₂O
 - Section 6.3.3 for Air
- 2. Turn the valve open 1/2 turn beyond the maximum indicated flow.
- Position the maximum stop collar so that its hook contacts the stop collar tab on the counterclockwise side. The hook and tab should have overlapping contact of about 0.75 mm (about half the thickness of the stop collar tab).



Note: This illustration shows the maximum stop collar in a vertical position. The actual position of the maximum stop collar may vary for each needle valve.

- 4. Tighten the locking screw on the maximum stop collar.
- Turn the valve one full turn clockwise to make sure the hook does not contact the stop collar tab. If there is contact, move the maximum stop collar slightly forward.
- 6. Verify that you can turn the valve open 1/2 turn beyond the maximum indicated flow.
- 7. Turn the valve fully clockwise to the minimum stop.

6.4 Link system calibration

Before you start, make sure that:

- All parts are correctly installed.
- Stops on needle valves are set correctly.
- The machine meets leak check requirements.
- Confirm that the O₂ sensor measures 21% in room air and 100% in pure O₂.
 If not, calibrate the O₂ sensor.

Note: All illustrations in this section show ANSI flowmeter module positions. The order is reversed on ISO machines.

▲ WARNING

You must be in a well ventilated room or use a gas evacuation device at this time. Anesthetic vapors exhausted into the room air can be harmful to your health.

- 1. Set the system switch to Standby.
- 2. Remove the flowmeter panel shield (Section 4.9.1).
- 3. Put the plastic spacer on the N₂O needle valve spindle.
- Turn the O₂ and the N₂O needle valves clockwise to their minimum stop position.




Put the chain onto the O₂ knob/sprocket assembly and the N₂O sprocket.

Note: The N₂O sprocket set screws should be away from the valve.



6-18 07/04 1009-0356-000

 Install the chain and sprockets onto the needle valve stems as an assembly. Press the O₂ knob/sprocket against the O₂ minimum stop collar.




7. Tighten the setscrews in the $\rm O_2$ knob. Do not tighten the $\rm N_2O$ sprocket setscrews.

Note: If O_2 label is on the knob, turn the knob so that the identification label is horizontal before tightening the setscrews.


- 8. Turn on the ${\sf O}_2$ and the ${\sf N}_2{\sf O}$ gas supplies (pipeline or cylinder).
- 9. Set the system switch to On.
- 10. Adjust the needle valves:

• $\mathbf{0_2}$ needle valve: 200 ±10 mL/min. • $\mathbf{N_20}$ needle valve: 600 ±25 mL/min.

11. Turn the sprocket on the $\rm O_2$ knob sprocket assembly counterclockwise until it stops against the tab on the $\rm O_2$ knob. Do not allow the $\rm N_2O$ or $\rm O_2$ valve stems to rotate.




12. Push the N₂O sprocket against the plastic spacer.



- 13. Holding the $\rm O_2$ knob, rotate the $\rm N_2O$ sprocket counterclockwise until all slack is removed from the chain.
- 14. Lightly tighten both N₂O sprocket setscrews.
- 15. Turn the $\mathrm{N}_2\mathrm{O}$ needle valve clockwise to the minimum stop position.

6-20 07/04 1009-0356-000

 Install the N₂O knob. Turn the knob so that the identification label is horizontal before tightening the setscrews.



- 17. Turn the N_2O needle valve counterclockwise, and check that the oxygen flow increases as N_2O flow increases.
- 18. Turn the $\rm O_2$ needle valve clockwise, and check that the $\rm N_2O$ flow decreases as $\rm O_2$ decreases.
- 19. Check the proportioning system concentration (increasing N₂O flow). Observe the following precautions:
 - Start with both valves at the minimum setting.
 - Adjust only the N₂O needle valve.
 - Increase the N₂O flow as specified in the table below and make sure the O₂ concentration is in the allowed range.

Note: Allow the O_2 monitor to stabilize. At the lower flows, the O_2 monitor may take up to 90 seconds to stabilize.

 If you overshoot a setting, turn the O₂ needle valve clockwise until the N₂O flow decreases to the previous setting before continuing the test.

| Set the N_2O flow (L/min) | Measured 0 ₂ |
|-----------------------------|-------------------------|
| 0.8 | 22% to 29% |
| 1 | 22% to 29% |
| 2 | 22% to 29% |
| 6 | 22% to 29% |
| 9 | 22% to 29% |

- 20. Check the proportioning system concentration (decreasing O₂ flow). Observe the following precautions:
 - Turn the N₂O needle valve to the maximum setting.
 - Adjust only the O₂ needle valve.
 - Decrease the O₂ flow as specified in the table and make sure the O₂ concentration is in the allowed range.

Note: Allow the O_2 monitor to stabilize. At the lower flows, the O_2 monitor may take up to 90 seconds to stabilize.

 If you overshoot a setting, turn the N₂O needle valve counterclockwise until the O₂ flow increases to the previous setting before continuing the test.

| Set the 0_2 flow (L/min) | Measured 0 ₂ |
|----------------------------|-------------------------|
| 3 | 22% to 29% |
| 2 | 22% to 29% |
| 1 | 22% to 29% |
| 0.3 | 22% to 29% |

If both tests meet the criteria, calibration is complete (go to the next step). If either test fails to meet the criteria, return to step 10 and adjust N_2O to a lower or higher setting.

| If: | Then: | |
|--------------------|---------------------------|--|
| Concentration Low | Decrease N ₂ O | |
| Concentration High | Increase N ₂ O | |


Note: Adjusting the regulator pressure is not recommended. It has little effect on proportioning. If you have difficulty proportioning the system, you may need to replace either or both needle valves.


- 21. Tighten N₂O sprocket setscrews.
- 22. Set the system switch to Standby.
- 23. Turn off the O_2 and the N_2O gas supplies.
- 24. Check that all setscrews are tight.
- 25. Adjust all needle valves to minimum stop position.
- 26. Install flowmeter panel shield.

6-22 07/04 1009-0356-000

6.5 0₂ Flush Regulator

- 1. Bleed all gas pressure for the machine (Section 4.2).
- 2. Ensure that all cylinder and pipeline gauges read zero before proceeding.
- 3. Remove the upper rear panel (Section 4.3).
- 4. Remove the O₂ Flush Regulator output tubing. Attach a 6-mm tee and a test device to the open port.





- 5. Connect an O_2 pipeline supply or slowly open the O_2 cylinder valve.
- Push the flush button just enough to achieve a slight flow or open the auxiliary flowmeter if equipped with this option. Read the pressure shown on the test device.


The pressure should be 241 ± 7 kPa (35 ± 1.0 psi).

- 7. If adjustment is required:
 - a. Loosen the adjustment screw's jam nut.
 - b. Adjust the regulator (in small steps) to the above specification.
 - c. Tighten the jam nut.
 - d. Verify the reading.
- 8. Disconnect the pipeline supply or close the cylinder valve.
- Bleed gas pressure by pushing the flush button; then, disconnect the tee and test device.
- 10. Reattach the output tubing to the regulator.
- 11. Install the rear panel.

6.6 Airway pressure gauge

6.6.1 Zero the pressure gauge

- Attach a patient circuit to the Breathing System. Leave the patient end open.
- 2. Set the Bag/Vent switch to Bag.
- 3. Adjust the APL valve to maximum.
- 4. Remove the lens from the pressure gauge:
 - Apply a slight pressure with your thumb and finger tips around the outer edge of the lens.
 - Turn the lens counterclockwise to remove it.
- 5. Adjust the pressure gauge to zero.
- 6. Plug the patient circuit.
- 7. Press and release the O₂ flush button to sweep the needle across the pressure gauge.
- 8. Remove the plug from the patient circuit to relieve the pressure in the circuit and recheck the zero setting of the pressure gauge.
- 9. If required, repeat zero and span procedure.
- 10. Replace the lens cover.



6-24 07/04 1009-0356-000


6.6.2 Checking the pressure gauge accuracy

The accuracy of the airway pressure gauge can be checked by using the following:

- a low-pressure test device (digital manometer or test gauge) with an accuracy of ±2% of reading,
- · a low-pressure supply source (typically a syringe),
- and an airway pressure gauge test adapter.
- 1. Ensure that the pressure gauge is zeroed (Section 6.6.1).
- 2. Remove the upper rear panel.
- 3. Remove the existing tube from the back of the pressure gauge and connect the test adapter tube directly to the gauge.
- 4. Connect a low-pressure supply source (syringe) to one of the open tubes of the test adapter.
- 5. Connect a low-pressure test device to the remaining open tube of the test adapter.
- Adjust the pressure source to the following pressures as read on the airway pressure gauge. The test device gauge should read within the values indicated.

| Airway Pressure Gauge | Test Device |
|------------------------|---------------------------|
| 0 cm H ₂ 0 | 0 ±1 cm H ₂ 0 |
| 40 cm H ₂ O | 40 ±2 cm H ₂ O |
| -5 cm H ₂ O | -5 ±2 cm H ₂ O |





Notes

6-26 07/04 1009-0356-000

7 Troubleshooting

In this section

This section covers the troubleshooting procedures for the Aespire machine pneumatic systems. For troubleshooting electrical systems, refer to the 7100 Ventilator Technical Reference manual.

| 7.1 General Troubleshooting | -2 |
|---|-----|
| 7.2 Breathing System Leak Test Guide | -4 |
| 7.2.1 Breathing system leak test | '-5 |
| 7.2.2 Breathing System Troubleshooting Flowcharts | -7 |
| 7.2.3 Leak Isolation Tests 7-1 | 12 |

7.1 General Troubleshooting

⚠ WARNING

Objects in the breathing system can stop gas flow to the patient. This can cause injury or death:

- Do not use a test plug that is small enough to fall into the breathing system.
- Make sure that there are no test plugs or other objects caught in the breathing system.

| Problem | Possible Cause | Action | |
|---|---|--|--|
| High Pressure Leak | Pipeline leak | Use a leak detector or Snoop to check for source of lea Repair or replace defective parts. | |
| | 0 ₂ flush valve | Use a leak detector or Snoop to check for source of leak. Make sure tubing connections are tight. Replace valve if defective. | |
| | System switch | Use a leak detector or Snoop to check for source of leak. Make sure tubing connections are tight. Replace switch if defective. | |
| | Cylinder not installed properly | Make sure cylinder is correctly aligned. Verify that tee handles are tight. | |
| | Cylinder gauges | Use a leak detector or Snoop to check for source of leak. Replace gauge if defective. | |
| | Cylinder gaskets | Use a leak detector or Snoop to check for source of leak. Replace gasket if defective. | |
| | Relief valves | Use a leak detector or Snoop to check for source of leak. Replace valve if defective. | |
| Low Pressure Leak (with vaporizer mounted) | Vaporizer not installed properly | Reseat vaporizer if not installed properly. Have vaporizer serviced at vaporizer center if vaporizer leaks. | |
| | Missing or damaged o-ring on vaporizer manifold | Check condition of o-ring. Replace if missing or damaged. | |
| | Loose fill port | Check fill port. Tighten if loose. | |

7-2

| Problem | Possible Cause | Action |
|---|--|---|
| Low Pressure Leak (with or without vaporizer) | Leaking port valve on vaporizer manifold | Use the Vaporizer Manifold Valve Tester to check for leak. See Section 4.10.2 for instructions. If test fails, tighten, repair, or replace as needed. |
| | Leak at flowmeter head | If vaporizer manifold passed previous tests: Remove tubing from input side of head and occlude port. Perform leak test. If test fails: Check for damaged o-rings between flowmeter modules. Replace as needed. Check for damaged o-rings at top and bottom of flow tubes. Replace as needed. Check for cracked flow tube. Replace as needed. |
| | | If secondary regulator leaks, replace the complete module. |
| | | Note: An alternate method is to pressurize the system and use a leak detector or Snoop to check for source of leak. |
| | Leaking relief valve on vaporizer manifold | Remove relief valve. Occlude opening. Perform leak test. If test passes, replace valve. |
| | Leaking flush valve | Attach pressure measuring device on CGO. Replace valve if device shows increased pressure. |
| | Leaking system switch | Attach pressure measuring device on CGO. Replace switch if device shows increased pressure. |
| Bellows leak | Pop-off valve diaphragm not sealing properly | Disassemble pop-off valve; inspect and clean seats; reseat; reassemble. |
| | Bellows mounting rim loose | Remove rim and pop-off valve diaphragm; reseat diaphragm; snap rim (2) into place. |
| | Bellows improperly mounted or has a hole or tear | Check that only the last bellows convolute is mounted to the rim and that the ring roll is in the groove under the rim. Inspect the bellows for damage; replace. |
| Breathing System Leak | Absorber canister not installed properly | Install canister properly. |
| | Soda lime dust on canister seals | Clean seals and mating surfaces. |
| Breathing System Leak (Intermittent) | ACGO O ₂ sense check valve | Replace. |
| $\ensuremath{\mathrm{N_20}}$ flow does not decrease with $\ensuremath{\mathrm{O_2}}$ flow | Defective pilot regulator | Check pilot regulator. Replace if needed. |
| Unit displays low O ₂ pressure with pipeline but not with cylinders | Low O ₂ supply switch | Check switch. Calibrate or replace as appropriate. |
| Unable to begin mechanical | ABS not fully engaged | Remount ABS. |
| ventilation | No O ₂ supply | Check O ₂ supply. |
| | Defective Bag/Vent switch | Check Bag/Vent switch. |

7.2 Breathing System Leak Test Guide

Note

Always perform the low-pressure leak test (Section 3.5) on the machine before proceeding with these breathing system leak tests.

The procedure in Section 7.2.1 helps you isolate the leak: to Bag Mode components, to Vent Mode components, or to components that are common to both modes.

- If you have a similar leak in both the bag mode and the ventilator mode, you
 must consider the Flow Sensor Module, the Circuit Module, the Absorber
 Canister area, and the bulkhead components (including CGO tubing).
 Carefully inspect the circuit module for damaged seals or misassembly, and
 the seating of the O₂ sensor.
- If you have a larger leak in one area than the other (Vent or Bag), the leak is most likely NOT in the Flow Sensor Module, the Circuit Module, the Absorber Canister area, or the bulkhead ports.

Follow the troubleshooting flowcharts in Section 7.2.2 to determine the best sequence of tests for locating a breathing system leak.

The procedures in Section 7.2.3 test specific components of the breathing system for leaks.

⚠ WARNING

Objects in the breathing system can stop gas flow to the patient. This can cause injury or death:

- Do not use a test plug that is small enough to fall into the breathing system.
- Make sure that there are no test plugs or other objects caught in the breathing system.

7-4 07/04 1009-0356-000

7.2.1 Breathing system leak test

This test checks for leaks in Vent Mode and Bag Mode components. It is part of the overall checkout procedure, Section 3.7 "Breathing system tests." It is repeated here for testing convenience.

⚠ WARNING

Objects in the breathing system can stop gas flow to the patient. This can cause injury or death:

- Do not use a test plug that is small enough to fall into the breathing system.
- Verify that AGSS is operating. For systems that have a flow indicator on the side, make sure that the flow indicator shows a flow in the green (normal) region.
- 2. Zero the pressure gauge (Section 6.6.1).

Check Valves

- Make sure that the check valves on the breathing circuit module work correctly:
 - The Inspiratory check valve rises during inspiration and falls at the start of expiration.
 - The Expiratory check valve rises during expiration and falls at the start of inspiration.

Ventilator Bellows

- 4. Ventilator bellows test:
 - a. Set the system switch to Standby.
 - b. Set the Bag/Vent switch to Ventilator.
 - c. Set all flow controls to minimum.
 - d. Close the breathing circuit at the patient connection. Use the test plug located on the side of the ABS.
 - e. Push the O₂ flush button to fill the bellows.
 - f. The pressure must not increase to more than 15 cm H₂O on the pressure gauge.
 - g. If the bellows falls more than 100 mL/min (top of indicator), it has a leak.

Service Mode Tests

- 5. Enter the Service Mode: Push and hold the adjustment knob on the ventilator's display and set the system switch to On.
 - a. Select and confirm "Service Modes."
 - Follow the menu structure outline below to reach the adjustment for the inspiratory flow valve. Select and confirm at each step.
 - "Diagnostics Tests/Tools"
 - "Breathing System Leak Test"
 - c. Follow the instructions on the screen.
 - d. The leak rate should be less than 200 mL/min.
 For machines with a single-tube O₂ flowmeter, the pressure reading should reach 30 cm H₂O at minimum flows greater than 200 mL/min.

Note: If test fails, see Section 7.2, "Breathing System Leak Test Guide".

Bag Circuit 6. Test the Bag circuit for leaks:

- a. Set the system switch to On.
- b. Set the Bag/Ventilator switch to Bag.
- c. Plug the Bag port (use your hand or the approved test plug).
- d. Close the APL valve (70 cm H₂0).
- e. Set the O2 flow to 0.25 L/min.
- f. Close the patient connection (using a hand or test plug on the side of the breathing system) and pressurize the bag circuit with the O₂ flush button to approximately 30 cm H₂O.
- g. Release the flush button. The pressure must not decrease. A pressure decrease large enough to see on the gauge indicates an unacceptable leak.

Note: If test fails, see Section 7.2, "Breathing System Leak Test Guide".

APL Valve


- 7. Test the APL valve:
 - a. Fully close the APL valve (70 cm H₂0).
 - b. Set the total fresh gas flow to approximately 3 L/min and make sure that the value on the inspiratory pressure gauge does not exceed 85 cm H₂O. Some pressure fluctuation is normal.
 - c. Fully open the APL valve (to the MIN position).
 - d. Set O2 flow to 3 L/min. Turn any other gases off.
 - e. Make sure that the value on the inspiratory pressure gauge is less than approximately $5\ cm\ H_2O$.
 - f. Push the O₂ flush button. Make sure that the value on the inspiratory pressure gauge stays less than 10 cm H₂O.
 - g. Set the O₂ flow to minimum and make sure that the value on the inspiratory pressure gauge does not decrease below 0 cm H₂O.
- 8. Remove your hand or the test plug from the patient connection.
- 9. Set the System switch to Standby.

⚠ WARNING


Make sure that there are no test plugs or other objects caught in the breathing system.

7-6 07/04 1009-0356-000


7.2.2 Breathing System Troubleshooting Flowcharts




Flowchart 1




7-8 07/04 1009-0356-000





7-10 07/04 1009-0356-000



7.2.3 Leak Isolation Tests

The previous flowcharts refer you to the following tests.

These tests require the use of the Low Pressure Leak Test Device and the Leak Test Tool Kit (refer to Section 8.1, "Service tools — Anesthesia machine".

The Leak Test Tool Kit includes:

- the Machine Test Tool
- the Circuit Test Tool
- and various Test Plugs

When performing these tests, ensure that the ACGO selector switch is set to the ABS position.

| Test 1 | Verifying the integrity of the test tools |
|---------|--|
| Test 2 | Low-pressure leak testing the machine |
| Test 3 | Testing the airway pressure gauge, and Port 1 and Port 3 u-cup seals 7-15 |
| Test 4 | Testing the bag port cover, the APL valve, the Bag/Vent switch, and the negative pressure relief valve |
| Test 5 | Testing the APL diaphram |
| Test 6 | Testing the bellows module and the Bag/Vent switch 7-18 |
| Test 7 | Testing the bellows, the bellows pop-off valve, the bellows base manifold, and the Bag/Vent switch |
| Test 8 | Testing the bellows assembly |
| Test 9 | Testing the flow sensor module, the circuit module, and the soda lime canister . 7-21 |
| Test 10 | Testing the circuit module and the canister |
| Test 11 | Testing the circuit module |
| Test 12 | Testing the inspiratory side of the circuit module |
| Test 13 | Testing the negative pressure relief valve |
| Test 14 | Testing the flow sensors only |
| Test 15 | Testing a flow sensor including the Ventilator Monitoring Assembly and interfacing components |

⚠ WARNING

Objects in the breathing system can stop gas flow to the patient. This can cause injury or death:

- Do not use a test plug that is small enough to fall into the breathing system.
- Make sure that there are no test plugs or other objects caught in the breathing system.


⚠ CAUTION

Do not use $\rm O_2$ Flush for leak isolation tests. Do not leave pressurized systems unattended. High pressure and equipment damage may result.


7-12 07/04 1009-0356-000

Test 1 Verifying the integrity of the test tools

Machine Test Tool Front View



Back View





- 1. Verify integrity of low-pressure leak test device.
 - Put your hand on the inlet of the leak test device. Push hard for a good seal.
 - Squeeze the bulb to remove all air from the bulb.
 - If the bulb completely inflates in less than 60 seconds, replace the leak test device.



- 2. Attach the low-pressure leak test device to the Machine Test Tool.
- 3. Plug the two pressure orifices.
- 4. Repeatedly squeeze and release the hand bulb until it remains collapsed.
- 5. If the bulb inflates in less than 30 seconds, locate and correct the leak.

Test 2 Low-pressure leak testing the machine





- 1. Remove the breathing system from the machine.
- Attach the Machine Test Tool (using only the Thru Port) and the lowpressure leak test device to **Port 3** of the breathing system interface as shown above.

Note: To prevent damage to the airway pressure gauge, ensure that the gauge port (**Port 1**) is not connected to the Test Tool.

- 3. Set the ACGO selector switch to ABS.
- 4. Set the system switch to Standby.
- 5. Turn off all vaporizers.
- 6. Open the flow controls one and a half turns counterclockwise.
- 7. Compress and release the bulb until it is empty.
- 8. The vacuum causes the floats to move. This is usual. If the bulb completely inflates in 30 seconds or less, there is a leak in the low-pressure circuit.


7-14 07/04 1009-0356-000





- 1. Attach the Machine Test Tool to the breathing system interface ports (using the alignment post) as shown above.
- 2. Turn all of the flow controls fully clockwise (minimum flow).
- 3. Set the system switch to On.
- 4. Occlude the tapered plug.
 - At minimum flow, the airway pressure gauge reading should increase.
 - If not, there is a leak in the tested circuit.

Test 4 Testing the bag port cover, the APL valve, the Bag/Vent switch, and the negative pressure relief valve




- Separate the Bellows Module from the Circuit Module and re-install the Bellows Module.
- 2. Occlude the Bag Port connector.
- 3. Connect the Machine Test Tool to the interface ports as shown above.
- 4. Set the Bag /Vent switch to Bag and close the APL Valve (70 cm H_2O).
- 5. Slowly increase the O_2 flow to achieve 30 cm H_2O .
 - The leak rate is equal to the flow needed to maintain 30 cm H₂O.
 - The leak rate should be less than 200 mL/min.

Note: If the bellows rises, it indicates a leak in the Bag / Vent Switch.

7-16 07/04 1009-0356-000


Test 5 Testing the APL diaphram



Note If required, set up the Machine Test Tool and breathing system as shown in Test 4.

- 1. Slide the Bellows Module away from the machine.
- 2. Remove the APL ramp and diaphragm.
- 3. Insert a Test Plug into the APL scavenging port, as shown above.
- 4. Slide the Bellows Module partially back onto the machine casting.
- 5. Ensure that the Bag Port is plugged and that the Bag/Vent switch is set to Bag.
- 6. Slowly increase the O₂ flow to achieve 30 cm H₂O.
 - Leak rate is equal to the flow needed to maintain 30 cm H₂0.
 - The leak rate should be less than 200 mL/min.

Test 6 Testing the bellows module and the Bag/Vent switch




- Separate the Bellows Module from the Circuit Module and re-install the Bellows Module.
- 2. Enter the Service Mode: Push and hold the adjustment knob on the ventilator's display and set the system switch to On.
 - a. Select and confirm "Service Modes."
 - b. Follow the menu structure outline below to reach the adjustment for the inspiratory flow valve. Select and confirm at each step.
 - "Diagnostics Tests/Tools"
 - · "Breathing System Leak Test"
- 3. Follow the instructions on the screen.
 - At step 5, connect the Machine Test Tool to the interface ports as shown above, instead. Continue with steps 6 and 7 on the screen.
 - The leak rate should be less than 200 mL/min.



7-18 07/04 1009-0356-000

Test 7 Testing the bellows, the bellows pop-off valve, the bellows base manifold, and the Bag/Vent switch





- 1. Separate the Bellows Module from the Circuit Module.
- 2. Insert appropriate test plugs into the bellows base manifold as shown to the left.


Note: Position the bellows assembly so that the bellows remain collapsed as you plug the ports.

- 3. Set Bag/Vent switch to Vent.
- 4. Position the bellows upright with the bellows collapsed.
- 5. Connect the Machine Test Tool to the interface ports as shown above.
- 6. Slowly increase the O2 flow to achieve 30 cm H2O.

Note: The bellows will rise until the pressure equalizes.


- The leak rate is equal to the flow needed to maintain 30 cm H₂O.
- The leak rate should be less than 200 mL/min.

Test 8 Testing the bellows assembly



Note

If required, set up the Machine Test Tool and breathing system as shown in Test 7.




2. Insert appropriate test plugs into the bellows base manifold as shown to the left.

Note: Position the bellows assembly so that the bellows remain collapsed as you plug the ports.

- 3. Connect the tapered plug of the Machine Test Tool to the bellows base inlet as shown to the left.
- 4. Position the bellows upright with the bellows collapsed.
- 5. Slowly increase the O_2 flow to achieve 30 cm H_2O .

Note: The bellows will rise until the pressure equalizes.

- The leak rate is equal to the flow needed to maintain 30 cm H₂O.
- The leak rate should be less than 200 mL/min.



7-20 07/04 1009-0356-000

Test 9 Testing the flow sensor module, the circuit module, and the soda lime canister



- Separate the Bellows Module from the Circuit Module and re-install the Circuit/Flow Sensor Module.
- Connect short tubing between the inhalation and exhalation ports of the breathing system.
- 3. Insert an appropriate test plug in the outlet port of the Circuit Module.
- 4. Enter the Service Mode: Push and hold the adjustment knob on the ventilator's display and set the system switch to On.
 - a. Select and confirm "Service Modes."
 - b. Follow the menu structure outline below to reach the adjustment for the inspiratory flow valve. Select and confirm at each step.
 - · "Diagnostics Tests/Tools"
 - "Display A/D Channels"
 - c. Record the Inspiratory and Expiratory Flow actual values.

Note: The Inspiratory and Expiratory Flow actual values should be near zero.

- Slowly increase the O₂ flow to achieve 30 cm H₂O.
 - The leak rate is equal to the flow needed to maintain 30 cm H₂O.
 - The leak rate should be less than 200 mL/min.
- 6. Observe the Inspiratory and Expiratory Flow actual values. The values should be near zero, as previously recorded in step 4.

Note: If one channel indicates flow, see "Inaccurate Volume Ventilation Troubleshooting" in the 7100 Ventilator Service Manual.


7. Release Pressure.

Test 10 Testing the circuit module and the canister



- 1. Remove the Flow Sensor module.
- 2. Connect the Circuit Test Tool to the Circuit Module as shown above.
- 3. Slowly increase the O_2 flow to achieve $30 \text{cmH}_2 O$. The leak rate is equal to the flow required to maintain $30 \text{ cmH}_2 O$.
 - The leak rate should be less than 200 mL/min.

Test 11 Testing the circuit module




Note: If required, set up the machine as in Test 10.

- 1. Remove the Soda Lime Canister.
- 2. Using appropriate Test Plugs, plug the three canister ports in the Circuit Module as shown above.
- 3. Slowly increase the 0_2 flow to achieve 30 cm H_2O .
 - The leak rate is equal to the flow needed to maintain 30 cm H₂0.
 - The leak rate should be less than 200 mL/min.

7-22 07/04 1009-0356-000

Test 12 Testing the inspiratory side of the circuit module




Note: If required, set up the machine as in Test 10 and 11.

- 1. Connect the Circuit Test Tool to the Circuit Module as shown above.
- 2. Insert an appropriate test plug in the inspiratory outlet to the canister as shown above.
- 3. Slowly increase the O_2 flow to achieve 30 cm H_2O .
 - The leak rate is equal to the flow needed to maintain 30 cm H₂0.
 - The leak rate should be less than 200 mL/min.

Test 13 Testing the negative pressure relief valve

- 1. Separate the Bellows Module from the Circuit Module.
- 2. Remove the Bellows Interface Manifold.
- 3. Insert test plug (recessed end) into the rear Bag/Vent switch port as shown.




- 4. Install the Bellows Module.
- 5. Connect the Machine Test Tool to the interface ports and the Bellows Module as shown above.



- 6. Set the Bag/Vent Switch to Vent.
- 7. Slowly increase the O_2 flow to achieve 30 cm H_2O .
 - The leak rate is equal to the flow needed to maintain 30 cm H₂0.
 - The leak rate should be less than 200 mL/min.


7-24 07/04 1009-0356-000

Test 14 Testing the flow sensors only



- 1. Remove the Flow Sensor Module.
- 2. Plug each Flow Sensor as shown above.
- Connect the low-pressure leak test device to the open end of the Flow Sensor.
- 4. Block the connector end of the Flow Sensor with your hand.
- 5. Compress and release the bulb until it is empty.
- 6. If the bulb inflates in 30 seconds or less, there is a leak in the flow sensor.
- 7. If there are no leaks in the flow sensors, go to Test 15.

Test 15 Testing a flow sensor including the Ventilator Monitoring Assembly and interfacing components



- 1. Remove Flow Sensors from the Flow Sensor Module.
- 2. Attach the Flow Sensor to the bulkhead connector.
- 3. Plug each Flow Sensor as shown.
- Connect the low-pressure leak test device to the open end of the Flow Sensor.
- 5. Compress and release the bulb until it is empty.
- If the bulb inflates in 30 seconds or less, there is a leak. The leak may be through the connector o-rings, in the internal tubbing, or in the Transducer on the VMB.

7-26 07/04 1009-0356-000

8 Illustrated Parts

| In this section | 8.1 Service tools – Anesthesia machine |
|-----------------|---|
| | 8.1.1 Test Devices |
| | 8.1.2 Test Tools |
| | 8.1.3 Secondary regulator pilot pressure tool8-5 |
| | 8.2 External components - front view |
| | 8.3 External components - front view references |
| | 8.4 External Components - rear view |
| | 8.5 Control module mounting for a ProTIVA machine |
| | 8.6 Aespire 100 - exclusive components |
| | 8.6.1 AC Inlet (Aespire 100) |
| | 8.6.2 Display mount (Aespire 100) |
| | 8.7 Front panel, gauges and system switch |
| | 8.8 Rear panel components |
| | 8.9 Tabletop components |
| | 8.10 Right-side Components |
| | 8.11 External components - lower assembly |
| | 8.12 Vent Engine Housing |
| | 8.13 Display cables, serial board, AGSS flowtube, and sample return |
| | 8.14 AC Power cords |
| | 8.15 AC Inlet/Outlet Components |
| | 8.16 Pipeline inlet fittings |
| | 8.17 Cylinder Gas Supplies |
| | 8.17.1 Cylinder inlet fittings |
| | 8.18 Vaporizer manifold |
| | 8.19 Flowmeter components |
| | 8.19.1 Flowtube parts |
| | 8.19.2 Secondary regulator components |
| | 8.20 ABS to machine Interface Components |
| | 8.20.1 Flush Regulator, Flush Valve, and ACGO Selector Switch 8-35 |

| 8.21 Breathing system interface |
|--|
| 8.22 Breathing System 8-37 |
| 8.22.1 APL Valve |
| 8.22.2 Bag/Vent Switch |
| 8.22.3 Absorber canister |
| 8.22.4 Flow Sensor Module |
| 8.22.5 Breathing Circuit Module |
| 8.22.6 Exhalation valve |
| 8.22.7 Bellows 8-43 |
| 8.22.8 Bellow base |
| 8.22.9 Bag Arms 8-45 |
| 8.23 Drawer 8-46 |
| 8.24 Legris quick-release fittings |
| 8.25 Vent Drive and low-pressure tubing |
| 8.26 Tubing for use with Legris fittings |
| 8.27 Cables and harnesses |
| 8.28 Cables and harnesses (Aespire 100) |
| 8.29 Anesthetic Gas Scavenging System – AGSS |
| 8.29.1 Passive AGSS 8-56 |
| 8.29.2 Adjustable AGSS 8-58 |
| 8.29.3 Active AGSS 8-60 |
| 8.30 Integrated Suction Regulator |
| 8.30.1 Major Components (Continuous and Venturi suction) |
| 8.30.2 Suction Control Module |
| 8.30.3 Venturi assembly |
| 8.31 Auxiliary O_2 Flowmeter |
| 8.32 Display mounts |
| 8.33 Cable management arm |
| 8.34 Display arm mounting kits for optional equipment |

8-2 07/04 1009-0356-000

8.1 Service tools — Anesthesia machine

8.1.1 Test Devices

| Item | Tool | | Stock Number |
|---------|---|---------------------------|------------------------|
| 1 | Test flowmeter, 6-50 L/min (Suction Flow Test) | | 1006-8431-000 |
| Not Sho | own | | |
| | Low-pressure Leak Test Device | (negative pressure) | 0309-1319-800 |
| | Low-pressure Leak Test Device | (positive pressure - ISO) | 1001-8976-000 |
| | Low-pressure Leak Test Device | (positive pressure - BSI) | 1001-8975-000 |
| | Flow test device capable of measuring 0–15 L/min with an accuracy of $\pm 2\%$ of reading | | Refer to section 3.8 |
| | Vacuum test gauge capable of measuring 0 to 550 mm Hg with an accuracy of $\pm 1\%$ of reading | | Refer to section 3.9 |
| | Test device capable of measuring 0–30 L/min (see Item 1 above | e) | Refer to section 3.9 |
| | Leakage current test device | | Refer to section 3.11 |
| | Test device capable of measuring 689 kPa (100 psi) | | Refer to section 6.1 |
| | Low-pressure test device (digital manometer or test gauge) with an accuracy of $\pm 2\%$ of reading | | Refer to section 6.6.2 |


Test Devices

1



8.1.2 Test Tools


| Item | Tool | Stock Number |
|------------|---|----------------|
| 1 | Leak Test Tool Kit, ABS breathing system | 1407-7013-000 |
| 1a | Test Tool, bulkhead | 1407-8500-000 |
| 1b | Plug, tapered 27x12 mm | 1407-8505-000 |
| 1c | Plug, tapered 24x18 mm | 1407-8506-000 |
| 1d | Test Tool, circle module (2 each) | 1407-8502-000 |
| 1e | Plug, service B/S 11 mm (2 each) | 1407-8504-000 |
| 1 f | Plug, service BTV 18 mm (2 each) | 1407-8503-000 |
| 2 | Plug, stopper | 2900-0001-000 |
| 3 | Adapter, positive low-pressure leak test | 1009-3119-000 |
| 4 | Vaporizer Manifold Valve Test Tool | 1006-3967-000 |
| Not Sho | own . | |
| | Tool to help disconnect tubing from Legris fittings | 2900-0000-000 |
| | Test Lung | 0219-7210-300 |
| | Leak detection fluid, Snoop | obtain locally |



8-4 07/04 1009-0356-000


8.1.3 Secondary regulator pilot pressure tool

Assemble the secondary regulator pilot pressure tool using a 4-mm tee and tubing as shown. This tool is used with $\rm N_2O$ needle valve calibration.



| Item | Description | Stock Number |
|------|---|---------------|
| 1 | Tee, 4 mm, tube/tube/tube | 1202-3653-000 |
| 2 | Tubing, 4 mm (approximately 450 mm - 18 inches) | 1001-3060-000 |
| 3 | Tubing, 4 mm (approximately 50 mm - 2 inches) | 1001-3060-000 |


8.2 External components - front view



| Item | Description | Stock Number |
|------|---|---------------|
| 1 | Caster, 125-mm with brake (front) | 1006-3070-000 |
| 2 | Caster, 125-mm no brake (rear) | 1006-3071-000 |
| 3 | Cover, cable channel | 1009-3020-000 |
| 4 | Upper shelf | 1009-3022-000 |
| 5 | Bolt, M6x40 | 0144-2131-911 |
| 6 | Lockwasher, M6 internal | 0144-1118-130 |
| 7 | Gauge, airway pressure (includes mounting hardware) | 1009-3034-000 |
| 8 | Handle, side | 1009-3033-000 |
| 9 | Handle, Medirail | 1009-3101-000 |
| 10 | Screw, M6x12 Sems | 0144-2436-106 |
| 11 | Spacer | 1009-3102-000 |
| 12 | Screw. M6x70 | 0144-2131-923 |
| 13 | Lockwasher M6 external | 9213-0560-003 |
| 14 | Screw, M6x20 | 0144-2131-921 |
| 15 | Shim | 1009-3131-000 |
| 16 | Task Light PCB | 1009-5504-000 |
| 17 | Lens, Task Light | 1011-3308-000 |
| | Screw | 0142-4254-106 |
| 18 | Switch Assembly, task light | 1009-5587-000 |
| 19 | Plate, switch mounting retainer | 1009-3143-000 |
| | Screw | 0140-6226-107 |


8-6 07/04 1009-0356-000

8.3 External components - front view references



| Item | Description | Section number |
|------|---|--|
| 1 | "Serial board" | Refer to section 8.13 |
| 2 | "AGSS gauge, and sample return" | Refer to section 8.13 |
| 3 | "Vent Engine Housing" | Refer to section 8.12 |
| 4 | "Anesthetic Gas Scavenging System — AGSS" | Refer to section 8.29 |
| 5 | "Breathing System" | Refer to section 8.22 |
| 6 | "Drawer" | Refer to section 8.23 |
| 7 | "Tabletop components" | Refer to section 8.9 |
| 8 | "Display cables" "Display mounts" | Refer to section 8.13 Refer to section 8.32 |
| 9 | "Flowmeter components" | Refer to section 8.19 |
| 10 | "Auxiliary O ₂ Flowmeter" | Refer to section 8.31 |
| 11 | "Integrated Suction Regulator" | Refer to section 8.30 |
| 12 | "Front panel, gauges and system switch" | Refer to section 8.7 |
| 13 | "ABS to machine Interface Components" | Refer to section 8.20 |
| 14 | "Vaporizer manifold" | Refer to section 8.18 |
| 15 | "Right-side Components" | Refer to section 8.10 |
| 16 | "External components - lower assembly" | Refer to section 8.11 |


8.4 External Components - rear view



| Item | Description | Stock Number |
|------|---|---|
| 1 | AC Inlet | Refer to section 8.14 |
| 2 | Pipeline Inlets Label, pipeline inlet blank | Refer to section 8.16 1009-3197-000 |
| 3 | Cylinder Gas Supplies | Refer to section 8.17 |
| 4 | Electrical Power Outlet Blank panel (no outlets) Screw, M4x12 | Refer to section 8.15 1011-3329-000 1009-3109-000 |
| 5 | Suction items | Refer to section 8.30 |
| 6 | Rear panel items | Refer to section 8.8 |
| 7 | Thumbscrew Ring, retainer | 1406-3304-000 1406-3319-000 |

8-8

8.5 Control module mounting for a ProTIVA machine




| Item | Description | Stock Number |
|------|--|--------------------------------|
| 1 | Thumbscrew, M6x14 | 1505-3005-000 |
| 2 | Bracket, mount 7100 display TIVA | 1009-3290-000 |
| 3 | Screw, M4x12 Pozidriv Lockwasher, M4 external | 0140-6226-111 9213-0540-003 |
| 4 | Block, mount 7100 display TIVA | 1009-3294-000 |
| 5* | Screw, M4x8 Pozidriv Flat HD | 0140-6226-107 |
| 6 | Clamp, cable TIVA | 1009-3291-000 |
| 7 | Screw, M6x14 SKT HD CAP | 0144-2131-922 |
| 8 | Manifold, display mount TIVA | 1009-8232-000 |
| 0.0 | V 1916 191729 | |

^{*} Apply Loctite 242.


8.6 Aespire 100 - exclusive components

| Item | Desc | ription | Stock Number |
|------|-------|--|------------------------|
| 1 | AC In | llet | Refer to section 8.6.1 |
| 2 | Disp | lay mount components | Refer to section 8.6.2 |
| 3 | (the | er supply, 40W universal display/control module in the Aespire 100 a different power supply than in a standard ire machine) | 1609-3040-000 |
| 4 | (with | head components the following exceptions, all components of 2 only flowhead are covered in the wmeter components" section. | Refer to section 8.19 |
| | 4a | Knob, O ₂ (does not include label) | 1006-3634-000 |
| | 4b | Plate, flowtube blank | 1006-1506-000 |
| | 4c | Screw, tread forming 0.88 inch (this screw is slightly longer to compensate for the flowtube blank plate) | 1009-3410-000 |
| 5 | (whe | k Label Set n replacing the flowmeter bezel, use opriate label to cover unused flow control iions) | 1009-3409-000 |
| 6 | Harn | esses | Refer to section 8.28 |
| 7 | Cast | er, non-locking 4-inch | 1009-3038-000 |
| 8 | Cast | er, locking 4-inch | 1009-3039-000 |
| 9 | Pane | l, lower rear cover | 1009-3407-000 |
| 10 | Scre | w, M4x8 | 1006-3178-000 |
| 11 | Thun | nbscrew | 1406-3304-000 |
| 12 | Ring | retaining | 1406-3319-000 |

8-10 07/04 1009-0356-000




8.6.1 AC Inlet (Aespire 100)



| Item | Description | Stock Number |
|------|-----------------------------|-----------------------|
| 1 | AC Inlet, with line filter | 1609-3019-000 |
| 2 | Fuse drawer | 1609-3020-000 |
| 3 | Fuse, 2A | 1009-5778-000 |
| 4 | Screw, M3x8 | 9211-0530-083 |
| 5 | Nut, M3 Keps | 0144-3717-302 |
| 6 | Transient Suppression board | 1609-3101-000 |
| 7 | Standoff, 52-mm Long | 1609-3320-000 |
| 8 | Screw, M4x8 | 0140-6226-113 |
| 9 | Stud, Equal Potential, 6-mm | 0208-0070-300 |
| 10 | Screw, M4x12 | 1009-3109-000 |
| 11 | Power Cord | Refer to section 8.14 |
| 12 | Clamp, power cord retainer | 1009-3405-000 |
| 13 | Screw, M4x8 | 1006-3178-000 |

8-12 07/04 1009-0356-000


8.6.2 Display mount (Aespire 100)



| Item | Description | Stock Number |
|------|--------------------------------------|---------------------|
| 1 | Dovetail insert | 1009-3408-000 |
| 2* | Screw, M8x20 SKT HD | 0144-2440-821 |
| 3* | Screw, M6x12 SKT HD | 0144-2436-101 |
| 4 | Extrusion, Vent bracket mount | 1504-3514-000 |
| 5 | Bolt, carriage | 1006-1433-000 |
| 6 | Casting, Vent Bracket, 20 degree (2) | 1504-3526-000 |
| 7 | Bearing, white plastic (2) | 1006-3228-000 |
| 8 | Handle, T-clamping | 1301-3001-000 |
| 9 | Washer, flat | 9213-0180-006 |
| 10 | Set screw, M6x10 | 0141-4227-111 |
| 11 | Screw, M4x16 | 9211-0440-163 |
| 12 | Lockwasher, M4 | 0144-1118-128 |
| 13 | Clamp | 1504-3527-000 |

^{*} Apply Loctite 242


8.7 Front panel, gauges and system switch



| Item | Description | Stock Number | Stock Number (pipeline) | Stock Number (cylinder) |
|------|--|--|--|---|
| 1 | Panel, gauge front | 1009-3018-000 | | |
| 2 | Switch, D-O system | 1006-8452-000 | | |
| 3 | Gauge, low pressure (includes mounting hardware) | | 1009-3079-000 | |
| 4 | Connector, 1/8 inch Legris to 10-32 | | 1006-3711-000 | |
| 5 | Gauge, high pressure (includes mounting hardware) | | | 1009-3080-000 |
| 6 | Connector, 1/8 inch copper tube to 5/16-24 | | | 1006-3712-000 |
| 7 | Label, gauge | O ₂ ANSI N ₂ O ANSI Air ANSI O ₂ ISO N ₂ O ISO Air ISO O ₂ Neutral N ₂ O Neutral Air Neutral | 1009-3081-000 1009-3082-000 1009-3083-000 1009-3202-000 1009-3203-000 1009-3234-000 1009-3235-000 1009-3236-000 | 1009-3199-000 1009-3201-000 1009-3200-000 1009-3204-000 1009-3201-000 1009-3205-000 1009-3237-000 1009-3238-000 1009-3239-000 |
| 8 | Tube Kit, copper tube and fittings (inboard cylinder) Tube Kit, copper tube and fittings (3rd gas) | 1006-8371-000 1006-8372-000 | 1009-3230-000 | 1003-3233-000 |
| 9 | Plate, gauge blanking | 1009-3045-000 | | |
| 10 | Plate, gauge blank backing | 1009-3147-000 | | |
| 11 | Palnut | 1009-3090-000 | | |

8-14 07/04 1009-0356-000


8.8 Rear panel components



| Description | Stock Number |
|---|--|
| Cover, rear upper | 1009-3073-000 |
| Cap, hose reel | 1009-3075-000 |
| Screw, M5.5x20 | 1009-3384-000 |
| Strap, hook/loop | 1009-3233-000 |
| Screw, M6x1.0 captive | 1009-3114-000 |
| Door, access (not functional for Aespire) | 1009-3074-000 |
| Screw, M4x12 | 1009-3109-000 |
| Spring, cantilever | 1009-3124-000 |
| Screw, M3x8 | 0142-4254-106 |
| Cover, trap bottle (if no internal suction) | 1009-3173-000 |
| Cover, regulator yoke (if no regulator) | 1009-3121-000 |
| Plate, clip cover | 1009-3185-000 |
| Wrench, DIN cylinder (with cable) | 1202-3651-000 |
| Wrench, pin index cylinder (with cable) | 0219-3415-800 |
| Cable | 1010-3049-000 |
| Ferrule, cylinder wrench cable retainer | 1001-3708-000 |
| Handle, P-grip | 1009-3343-000 |
| | Cover, rear upper Cap, hose reel Screw, M5.5x20 Strap, hook/loop Screw, M6x1.0 captive Door, access (not functional for Aespire) Screw, M4x12 Spring, cantilever Screw, M3x8 Cover, trap bottle (if no internal suction) Cover, regulator yoke (if no regulator) Plate, clip cover Wrench, DIN cylinder (with cable) Wrench, pin index cylinder (with cable) Cable Ferrule, cylinder wrench cable retainer |

^{*} Clean mounting surface with isopropyl alcohol.

8.9 Tabletop components




| Item | Description | Stock Number |
|------|---|---------------|
| 1 | Tabletop, work surface | 1009-3029-000 |
| 2 | Screw, relieved | 1504-3001-000 |
| 3 | Washer, retainer | 1009-3178-000 |
| 4 | Window, check-valve | 1009-3088-000 |
| 5 | Palnut | 1009-3090-000 |
| 6 | Hook, breathing circuit | 1009-3086-000 |
| 7 | Bolt, shoulder | 1009-3172-000 |
| 8 | Washer, wave | 1009-3035-000 |
| 9 | Washer, Nylon | 1009-3150-000 |
| 10 | Clip, with tape (used with bag arm) | 1009-8196-000 |
| 11 | Clip, with tape (used with bag on hose) | 1009-8197-000 |

8-16 07/04 1009-0356-000


8.10 Right-side Components





| Item | Description | Stock Number |
|------|---------------------------|---------------|
| 1 | Extrusion cover | 1009-3021-000 |
| 2 | Screw, M6x20 | 0144-2131-921 |
| 3 | Lockwasher, M6 internal | 0144-1118-130 |
| 4 | Dovetail, RH upright | 1009-3129-000 |
| 5 | Screw, M4x10 self-tapping | 1009-5534-000 |
| 6 | Cover, pipeline inlet | 1009-3091-000 |
| 7 | Screw, M4x8 | 1006-3178-000 |


8.11 External components - lower assembly



| Item | Description | Stock Number |
|------|----------------------------|---------------|
| 1 | Panel, access | 1009-3059-000 |
| 2 | Panel, service | 1009-4141-000 |
| 3 | Screw, M4x8 | 1006-3178-000 |
| 4 | Thumbscrew | 1406-3304-000 |
| 5 | Ring, retaining | 1406-3319-000 |
| 6 | Cover, scavenger reservoir | 1009-3027-000 |
| 7 | Bracket, suction reservoir | 1009-3107-000 |
| 8 | Screw, M4x16 | 9211-0440-163 |
| 9 | Lockwasher, M4 external | 9213-0540-003 |
| 10 | Clip, suction bag hose | 1407-3327-000 |
| 11 | Screw, M5x16 PAN HD | 9211-8350-163 |
| 12 | Lockwasher, M5 external; | 0144-1118-220 |
| 13 | Bumper, absorber | 1009-3105-000 |


8-18

8.12 Vent Engine Housing



| Item | Description | Stock Number | Qty |
|------|--|---------------|-----|
| 1 | Vent Engine Cover Plate Assy | 1407-7009-000 | |
| 2 | TAB GUIDE BELLOWS BASE | 1407-3313-000 | (2) |
| 3 | SCR M3X16 POSI DR PAN HD A4 SST | 1504-3003-000 | (2) |
| 4 | Cap, Plug | 1406-3524-000 | |
| 5 | FITTING PNL MOUNT 3.18 HOSE BARB UNION | 1504-3014-000 | (2) |
| 6 | PLUG HOLE 15.9 DIA NYLON MICRO PLASTICS | 1006-1473-000 | |
| 7 | PLATE CONN VENT | 1407-3321-000 | |
| 8 | SCR M4X8 POZI-DR DIN84 PAN SERRATED | 1006-3178-000 | (5) |
| 9 | HARN VENT ENG BRD TO COIN SLOT | 1009-5545-000 | |
| 10 | BLOCK LATCHING DSUB CONN | 1504-3617-000 | (2) |
| 11 | SCR 4-40 X 3/8 SKT BCG HD CAP | 0144-2117-206 | (2) |
| 12 | Bracket, scavenging guard in machines without scavenging | 1407-3922-000 | |
| 13 | CLIP-SUCTION BAG HOSE | 1407-3327-000 | |
| 14 | SCR M5 X 16 PAN PH HD SST | 9211-8350-163 | (2) |
| 15 | Lockwasher | 0144-1118-220 | (2) |


8.13 Display cables, serial board, AGSS flowtube, and sample return




| Item | Description | Stock Number |
|------|--|---------------|
| 1 | CABLE SER ISLN CONN BRD TO CTRL MOD 7100 | 1009-5691-000 |
| 2 | CABLE MONITORING BOARD | 1504-5604-000 |
| 3 | CABLE PNEUMATIC ENGINE | 1504-5605-000 |
| 4 | Cable, power | 1009-5711-000 |
| 5 | "Harness, Serial ISO to 02 flush SW" | 1009-5567-000 |
| 6 | "Harness, Serial ISO to 02 supply SW" | 1009-5568-000 |
| 7 | "Harness, Serial ISO to on/ standby SW" | 1009-5566-000 |
| 8 | Serial Isolation Board | 1009-5500-000 |
| 9 | Standoff | 1202-3092-000 |
| | Lockwasher | 0144-1104-331 |
| 10 | Screw, thread forming | 1009-3400-000 |
| 11 | Label, flow indicator AGSS | 1406-3527-000 |
| | Label, flow indicator AGFS (for German variant) | 1009-3301-000 |
| | Label, blank (for machines without flow indicator) | 1009-3241-000 |
| 12 | Flowtube, AGSS | 1406-3560-000 |
| 13 | Clip, AGSS flowtube | 1009-3181-000 |
| 14 | Coupling, Colder (Kit includes mounting nut) | 1009-8321-000 |
| | Label set, blank (for machines without scavenging) | 1009-3351-000 |
| 15 | Label, AGFS (for German variant) | 1009-3300-000 |

8-20 07/04 1009-0356-000

8.14 AC Power cords






| Item | Description | Stock Number |
|------|--|---------------|
| 1 | Power Cord | |
| | Australia and China, 220-240 VAC AS 3112 outlets | 1006-3888-000 |
| | EURO and France, 220 VAC with CEE 7/7 | 1001-3380-000 |
| | India and South Africa, 220-240 VAC BS546 | 1006-3885-000 |
| | Japan and US, 100-120 VAC NEMA | 1006-3907-000 |
| | Peru, 220-240 VAC NEMA | 1006-3882-000 |
| | Swiss, 220-240 VAC SEV 1011 | 1006-3889-000 |
| | UK, 220-240 VAC BS1363 | 1006-3884-000 |
| 2 | Clamp, power cord retainer | 1009-3103-000 |
| 3 | Screw, M4x8 PAN serrated | 1006-3178-000 |
| 4 | Clamp, power cord retainer | 1009-3405-000 |

8.15 AC Inlet/Outlet Components

| Item | Description | Stock Number |
|------|--|---|
| 1 | Inlet, 100/120 AC, with line filter and 15 A circuit breaker Inlet, 220/240 AC, with line filter and 8 A circuit breaker | 1009-5698-000 1009-5757-000 |
| 2 | Fuse, 2A - 5x20mm Fuse holder | 1009-5778-000 1009-5674-000 |
| 3 | Circuit board, Inrush, 100-120V | 1006-3245-000 |
| | Circuit board, Inrush, 220-240V | 1006-3246-000 |
| | Circuit board, Transient Suppression, 100-120V | 1006-3788-000 |
| | Circuit board, Transient Suppression, 220-240V | 1006-3789-000 |
| 4 | Filter, AC Line, 100-240V | 1009-5675-000 |
| 5 | Cable, line filter to Display Module | 1009-5711-000 |
| 6 | Harness, 100/120 V to Toroid (or to outlets harness) Harness, 220/240 V to Toroid (or to outlets harness) Harness, if no outlets | 1009-5713-000 1009-5714-000 1009-5715-000 |
| 7 | Harness, to 100/120 V outlets Harness, to 220/240 V outlets | 1009-5716-000 1009-5717-000 |
| 8 | Outlet Receptacle, Australia, AS 3112 | 1001-3305-000 |
| | Outlet Receptacle, EURO, CEE 7/7 | 1202-3551-000 |
| | Outlet Receptacle, France, CEE 7/4 Support Frame, snap in | 1006-4421-000 1006-4422-000 |
| | Outlet Receptacle, India and South Africa, BS 546 | 1006-3805-000 |
| | Outlet Receptacle, Japanese | 1006-3578-000 |
| | Outlet Receptacle, NA, Nema 5-15 | 1006-3555-000 |
| | Outlet Receptacle, Swiss, SEV 1011 | 1006-3807-000 |
| | Outlet Receptacle, UK, BS1363 | 1001-3309-000 |
| 9 | Circuit Breaker, 1A, Rocker | 1009-5722-000 |
| | Circuit Breaker, 2A Rocker | 1009-5721-000 |
| | Circuit Breaker, 3A Rocker | 1009-5720-000 |
| | Circuit Breaker, 4A Rocker | 1009-5719-000 |
| 10* | Toroid, 100-240V Heatshrink tubing | 1009-5692-000 1202-3268-000 |
| 11 | Screw, M6x70 | 0144-2131-923 |
| 12 | Lockwasher, M6 | 9213-0560-003 |
| 13 | Washer | 0402-1107-500 |
| 14 | Cover, transformer | 1009-3063-000 |
| 15 | Screw, M4x8 DIN84 (for transformer cover) | 1006-3178-000 |
| 16 | Stud, Equal Potential, 6 mm | 0208-0070-300 |

^{*} Apply heatshrink tubing to terminals of black and white wires if not being used.

8-22







AS 3112 Australia/China





Nema 5-15 Japanese





CEE 7/7 EURO





SEV 1011 Swiss,





CEE 7/4 France





BS1363 UK





BS 546 India and South Africa






Nema 5-15 NA

1009-0356-000 07/04


8.16 Pipeline inlet fittings




| Item | Description | | Stock Number |
|------|--|------------------|------------------------------|
| 1 | Pipeline inlet - 0 ₂ fittings | | |
| | Body, O ₂ DISS | | 1006-5149-000 |
| | Body, O ₂ NIST | | 1006-5158-000 |
| | Body, O ₂ DIN | | 1006-5161-000 |
| | Body, O ₂ G 3/8 BSPP | | 1006-5170-000 |
| | Pipeline inlet assembly O ₂ France | | 1006-8363-00 |
| | Pipeline inlet assembly 02 Canada | | 1006-8360-00 |
| | Pipeline inlet assembly O ₂ Australia | | 1006-8396-00 |
| 1 | Pipeline inlet - N20 fittings | | |
| | Body, N ₂ O DISS | | 1006-5150-000 |
| | Body, N ₂ O NIST | | 1006-5159-000 |
| | Body, N ₂ O DIN | | 1006-5162-000 |
| | Body, N ₂ O G 3/8 BSPP | | 1006-5171-000 |
| | Pipeline inlet assembly N ₂ O France | | 1006-8362-00 |
| | Pipeline inlet assembly N ₂ O Canada | | 1006-8359-00 |
| | Pipeline inlet assembly N ₂ O Australia | | 1006-8397-00 |
| 1 | Pipeline inlet Air fitting | | |
| | Body, Air DISS | | 1006-5151-000 |
| | Body, Air NIST | | 1006-5160-000 |
| | Body, Air DIN | | 1006-5163-00 |
| | Body, Air G 3/8 BSPP | | 1006-5172-000 |
| | Pipeline inlet assembly Air France (service kit) | | 1006-8361-00 |
| | Pipeline inlet assembly Air Canada (service ki | 3 | 1006-8358-00 |
| 0 | Pipeline inlet assembly Air Australia (service k | it) | 1006-8398-00 |
| 2 | O-ring, bore seal | | 0210 0470 200 |
| | O_2 and N_2O | | 0210-0479-300 |
| 3 | Air Sintered metal filter with o-ring | | 0210-0539-30 1006-8351-00 |
| 4 | Pipeline checkvalve with o-ring | | 1006-8351-000 |
| 5 | Gas Inlet Manifold (replacement) | 02 | 1009-8066-00 |
| ă | add mot mamora (representedly | N ₂ O | 1009-8067-00 |
| | | Air | 1009-8068-00 |
| 6 | Relief valve, 758 kPa (110 psi) | | 1011-3049-00 |
| 7 | Screw, M4x20 | | 0144-2124-218 |
| | Lockwasher, M4 external | | 9213-0540-00 |

8-24 07/04 1009-0356-000

8.17 Cylinder Gas Supplies






| Item | Description | Pin Index (Inboard) | DIN (Inboard) | DIN, Large Cylinder (Inboard) |
|------|-----------------------------|----------------------|---------------|----------------------------------|
| 1 | Gas supply 02 | 1006-3201-000 | 1006-3207-000 | 1006-3880-000 |
| 1 | Gas supply N ₂ O | 1006-3202-000 | 1006-3208-000 | 1006-3881-000 |
| 1 | Gas supply Air | 1006-3203-000 | 1006-3209-000 | |
| | | Pin Index (Outboard) | | |
| 2 | Gas supply N ₂ O | 1009-8210-000 | | |


| Item | Description | Stock Number |
|------|--|-------------------------|
| 3 | Standoff (3 per supply) | 1009-3085-000 |
| | Screw, M6x80 socket head cap (3 per supply) | 0144-2131-913 |
| | Lockwasher, M6 external (for above screw) | 9213-0560-003 |
| 4 | Screw, M6x25 socket head cap (3 per supply) | 9211-0660-254 |
| | Lockwasher, M6 external (for above screw) | 9213-0560-003 |
| 5 | Cylinder inlets (Pin Index or DIN for external cylinder) | Refer to section 8.17.1 |
| 6 | Tee handle beige | 0219-3372-600 |
| 7 | Clamp, yoke | 1001-4076-000 |
| 8 | Spacer, gas block (2) | 1001-4077-000 |
| | Screw, M8 x 25 long socket head cap (2) | 9211-0680-253 |
| 9 | Elbow fitting for cylinder pressure gauge (copper tube connection of gas supply) | 1006-3713-000 |
| 10 | Label Set, cylinder supply, O ₂ | 1006-3854-000 |
| | Label Set, cylinder supply, N ₂ O | 1006-3855-000 |
| | Label Set, cylinder supply, Air | 1006-3856-000 |


8.17.1 Cylinder inlet fittings

1 Pin Index



2 DIN (external cylinder)






| Item | | Description | Stock Number |
|------|----|--|---|
| 1 | | Cylinder inlets (Pin Index) | *************************************** |
| 1 | 1a | Gasket | 0210-5022-300 |
| 1t | b* | O-ring | 9221-3013-116 |
| 3 | 1c | Adapter, inlet | 1001-4075-000 |
| 1 | 1d | Filter, sintered bronze | 9914-6380-000 |
| 1 | 1e | Retaining ring, filter | 1001-5954-000 |
| 2 | | Cylinder inlets (DIN) | |
| 2 | 2a | Screw, M8x16 | 0144-2140-242 |
| 2 | 2b | Sealing ring (DIN) | 1009-3356-000 |
| | 2c | DIN Adapter (O ₂) | 1006-4000-000 |
| | | DIN Adapter (N ₂ 0) | 1006-4001-000 |
| | | DIN Adapter (Air) | 1006-4002-000 |
| 2 | 2d | O-ring, 0.687 ID, 0.812 OD | 0210-0544-300 |
| 2 | 2e | Filter, sintered bronze | 9914-6380-000 |
| | 2f | Sealing ring, N ₂ O DIN Conn 11 | 1202-3641-000 |
| | 2g | Adapter, large cylinder N ₂ O | 1006-4028-000 |
| 2 | 2h | Adapter, large cylinder O ₂ | 1006-4027-000 |

Lubricate sparingly with Krytox


8-26 07/04 1009-0356-000

8.18 Vaporizer manifold



| Item | Description | Stock Number |
|------|---|---------------|
| | Manifold assembly, complete, two position | 1006-8355-000 |
| | Manifold assembly, complete, one position | 1009-8065-000 |
| 1 | O-ring, 0.687 inch ID 0.812 inch OD | 0210-0544-300 |
| 2 | Spring, compression | 1006-3736-000 |
| 3 | Valve kit, includes seal | 1006-8373-000 |
| 4 | Seal | 1006-3690-000 |
| 5 | O-ring, 14.3 mm ID | 1102-3043-000 |
| | (Package of 6 o-rings) | 1102-3016-000 |
| 6 | Nipple, vaporizer port (New Style) | 1006-4215-000 |
| 7 | Screw, M2.5 - 0.45x6 PAN, Pozidriv, SST | 1006-3037-000 |
| 8 | Spring, Dzus | 1102-3056-000 |
| 9 | Seat, check valve | 1006-1352-000 |
| 10 | O-ring 27.1 OD 21.89 mm ID | 1006-3866-000 |
| 11 | Flapper | 0211-1451-100 |
| 12 | Housing | 1006-1351-000 |
| 13 | Screw, M4 x 30, cap head | 9211-0640-304 |
| 14 | Valve, relief, 5.5 psi, 7/16-20 THD | 1006-4128-000 |
| 15 | Flexible tubing, 1/4 inch, mixed gas | 1001-3064-000 |
| | | |

8.19 Flowmeter components



| Item | Description | Stock Number |
|------|---|---------------|
| 1 | Flowhead Module: includes regulator, flowtube module, flowtubes, needle valve, intermodule tube and associated o-ring, and label plate; does not include labels, link-25, or knobs (order separately). | |
| | O ₂ flowhead module with dual flowtubes | 1006-8380-000 |
| | O ₂ flowhead module with single flowtube | 1009-8069-000 |
| | N ₂ O flowhead module with dual flowtubes | 1006-8381-000 |
| | N ₂ O flowhead module with single flowtube | 1009-8070-000 |
| | Air flowhead module with dual flowtubes | 1006-8383-000 |
| | Air flowhead module with single flowtube | 1006-8382-000 |
| 2 | Secondary regulators/Balance Regulators | |
| | Regulator Kit, O2 (adjustable), without pressure switch | 1006-8341-000 |
| | Regulator Kit, N ₂ O (pressure balancing) | 1006-8344-000 |
| | Regulator Kit, Air (adjustable) | 1006-8340-000 |

8-28

| 7 Frenche Module: includes housing, o-rings, and plug fewtubes, label or table peanel (order separately). Flowtube module, O ₂ -single Frowtube module, O ₂ -single Frowtube module, N ₂ O-dual Fowtube module, N ₂ O-single Fowtube module connect) 7 Screw, MiSt-26 (O ₂) N ₂ O-3 organ module connect) 8 Screw, MiSt-26 (O ₂) N ₂ module connect) 9 Screw, MiSt-26 (O ₂) N ₂ module connect) 9 Nut, MiS Reps (O ₂ - Mi module connect) 9 Nut, MiS Reps (O ₂ - Mi module connect) 9 Screw, Mishel panel 11 Clly small formater module) 11 Clly small panel 11 Clly small side panel 11 Clly small side panel 11 Clly small side panel 12 Elbow, 1, M inch 13 Beave, (M-M-M-M-S-28) 8 Beave, (M-M-M-M-S-288) 8 Beave, (M-M-M-M-S-2888) 8 Beave, (M-M-M-M-M-M-M-M-M-M-M-M-M-M-M-M-M-M-M- | | Stock Number |
|---|--|----------------------|
| | Flowtube Module: includes hausing, o-rings, and plug ball; does not include, flowtubes, label or label panel (order separately). | |
| | e, 0 ₂ - dual | 1006-8338-000 |
| \$ 10 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 | e, 0 ₂ - single | 1009-8234-000 |
| | e, N ₂ 0 - dual | 1006-8337-000 |
| | e, N ₂ O - single | 1009-8235-000 |
| | e, Air - dual | 1006-8333-000 |
| | e, Air - single | 1006-8334-000 |
| | le connector | 1006-3628-000 |
| | ule connector | 1006-3613-000 |
| 000 | retaining | 1006-4350-000 |
| | module mounting) | 1102-3049-000 |
| (this is a second of the secon | ∂₂/№20/3rd gas module connect) | 1006-3607-000 |
| . 47 | (0 ₅ - Air module connect) | 1006-3080-000 |
| | - Air module connect) | 1006-1459-000 |
| | 2 - Air module connect) | 0144-3717-324 |
| | uri e | refer to chart below |
| | vmeter module | 1006-1290-000 |
| 1000 0000 00000000000000000000000000000 | - Air flowmeter module) | 1009-3186-000 |
| | la. | 1006-3608-000 |
| | ube outlet | 1009-3056-000 |
| Lancova | | 1009-3309-000 |
| | | 1202-3804-000 |
| 14 Bezel flowmeter 2-plas | 3-gas | 1009-3104-000 |
| and a second formation of | 7.2-gas | 1009-3108-000 |

by with Knytox.

The first by with Chydrox and the bards should face to the left.

Bell on the panel is other the right edge of the label is flush with the right edge of the panel. The left has not the panel is other the right of the label is flush with the right edge of the left edge of the panel. When mounted in the flowmeter module, half will extend slightly beyond the left edge of the panel. When mounted in the flowmeter module, bare in the right should slightly overlap the flowmeter label directly to its left.

set contains a label for both the dual-tube (A) and single-tube (B) modules. Make sure to install

| | Stock Number | Stock Number | Stock Number |
|---------|--------------------------------|----------------------------|------------------|
| ANSI | Air (Yellow) | N ₂ 0 (Blue) | 02 (Green) |
| | 1006-0209-000 | 1009-3209-000 | 1009-3210-000 |
| 051 | 0₂ (Neutral) | N ₂ 0 (Blue) | Air (Black/White |
| | 1009-3211-000 | 1009-3209-000 | 1009-3242-000 |
| leutral | 0 ₂ (Neutral) | N ₂ 0 (Neutral) | Air (Neutral) |
| | 1009-3211-000 | 1009-3240-000 | 1009-3240-000 |


1009-0356-000 07/04

8.29

8.19.1 Flowtube parts

| Item | Description | Stock Number |
|----------|---|---------------|
| 1 | Large flowtube kits (includes float, filter, o-rings, tube) | |
| | Flowtube kit, Air, large | 1006-8325-000 |
| | Flowtube kit, N ₂ O, large | 1006-8329-000 |
| | Flowtube kit, O2, large | 1006-8331-000 |
| 2 | Small flowtube kits (includes float, filter, o-rings, tube) | |
| | Flowtube kit, Air, small | 1006-8326-000 |
| | Flowtube kit, N ₂ O, small | 1006-8330-000 |
| | Flowtube kit, 0 ₂ , small | 1006-8332-000 |
| 3 | Single-tube flowtube kits (includes float, filter, o-rings, tube) | |
| | Flowtube kit, single-tube Air flowmeters use the Large flowtube kits | |
| | Flowtube kit, N2O, single-tube, dual-taper with filter | 1009-8199-000 |
| | Flowtube kit, 02, single-tube, dual-taper with filter | 1009-8198-000 |
| 4 | Spring, top of flowtubes | 1006-3624-000 |
| 5 | Float stop, O ₂ large | 1006-1225-000 |
| | Float stop, N ₂ O large | 1006-1226-000 |
| | Float stop, Air large | 1006-1227-000 |
| 6 | O-ring, 17.6 OD, 12.37 ID, large flowtube, top | 1006-3615-000 |
| 7 | Filter, large flowtube | 1006-3584-000 |
| 8 | O-ring, 17.6 OD, 12.37 ID, large flowtube, bottom (red) | 1006-3968-000 |
| 9 | Float stop, O ₂ small | 1006-1233-000 |
| | Float stop, N ₂ O small | 1006-1234-000 |
| | Floet stop, Air small | 1006-1235-000 |
| 10 | O-ring, 11.26 OD, 6.02 ID, small flowtube, top | 1006-3617-000 |
| 11 | Filter, small flowtube | 1006-3583-000 |
| 12 | O-ring, 11.26 OD, 6.02 ID, small flowtube, bottom (red) | 1006-3969-000 |
| 13 | Ball, 6 mm (plug fresh gas end) | 1006-1353-000 |
| Not Show | TIL. | |
| | O-ring Kit, flowtubes (includes 4 each of top and bottom o-rings for large flowtube and 3 each of top and bottom o-rings for small flowtube). | 1006-8393-000 |
| | Silicon tube kit, long, including cable ties | 1006-8378-000 |
| | Silicon tube kit, short, including cable ties | 1006-8379-000 |


8-30 07/04 1009-0356-000




8.19.2 Secondary regulator components

| Item | Description | Stock Number |
|------|--|---------------|
| 1 | Knob (N20) (Air) without label | 1006-3633-000 |
| 2 | O_2 Proportioning assembly (includes knob, sprocket, set screws, without knob label) | 1006-8339-000 |
| 3 | Set screw | 0141-4227-105 |
| 4 | O ₂ Proportioner chain | 1006-3610-000 |
| 5 | Sprocket, N ₂ O | 1006-3625-000 |
| 6 | Plug, 1/8 inch | 1006-3611-000 |
| 7 | Spacer, link system, $N_2\mathrm{O}$ needle valve | 1006-5140-000 |
| 8 | Valve, needle (O ₂) (Air) | 1006-8346-000 |
| | Valve, needle $N_2\mathrm{O}$ (has notch around valve body) | 1006-8345-000 |
| 9 | O-ring, 10.1 ID 13.3 OD | 9221-3010-116 |
| 10 | 0-ring, 0.250 inch ID 0.375 inch 0D | 0210-0687-300 |
| 11 | Stop collar (all gases) | 1006-3632-000 |
| 12 | Maximum stop collar kit (includes item 13) Maximum stop collars are required in Canada for all gas flow controls. | 1006-8055-000 |
| 13 | Mounting screw, M3x8, SKT HD CAP | 1006-3865-000 |
| 14 | Pressure switch, O ₂ supply alarm | 1006-3623-000 |
| 15 | O-ring, 0.250 inch ID 0.375 inch OD | 0210-0687-300 |
| 16 | Screws, M4x12 Pozidriv PAN | 0140-6226-111 |
| 17 | Plug, pressure switch cavity | 1006-3665-000 |
| 18 | Screws, M4x8 Pozidriv PAN | 1006-3178-000 |
| 19 | Fitting, O ₂ pilot, plug-in elbow | 1006-3533-000 |
| 20 | Fitting, 0_2 pilot, thread-in elbow | 1006-3663-000 |
| 21 | Tubing, 4-mm (RH head 144 mm - LH head 164 mm) | 1001-3060-000 |

8-32 07/04 1009-0356-000




8.20 ABS to machine Interface Components



| Item | Description | Stock Number |
|------------------------------|--|----------------------------|
| 1 | PCA VENT MONITORING BRD [TESTED VMB] | 1009-8002-000 |
| 2 | PCA B/S task light interface board (used when no monitoring is ordered) | 1009-5681-000 |
| 3 | SCR M4X12 POZI-DR PAN SSTTYPE 316 | 0140-6226-111 |
| 4 | WASH/LOCK EXT M4 DIN 6797 SSTTYPE 316 | 9213-0540-003 |
| 5 | SPACER.171 IN ID 0.250IN OD 0.250 IN LG ALUMINUM | 1202-3693-000 |
| 6 | Cables and harnesses | Refer to section 8.27 |
| 7 | 0 ₂ Flush Valve | Refer to section 8.20.1 |
| 8 | ACGO Selector Valve | Refer to section 8.20.1 |
| 9 | Port, ACGO body | 1009-3096-000 |
| 10 | Screw, M4x30 | 9211-0640-304 |
| 11 | Lockwasher, M4 | 9213-0540-003 |
| 12 | Cap, ACGO check valve | 1009-3095-000 |
| 13 | Screw, M4x8 | 9211-1040-069 |
| 14 | Disk, ACGO check valve | 1009-3062-000 |
| 15 | Flapper, ACGO check valve | 1009-3097-000 |
| 16* | 0-ring | 0210-0543-300 |
| 17 | Fitting, elbow barbed | 1009-3160-000 |
| 18* | O-ring | 0210-0691-300 |
| 19 | Screw, M3x6 | 9211-1030-055 |
| 20 | Tie wrap | 0203-5915-300 |
| Lubricat | te sparingly with Krytox. | |

3-34 07/04 1009-0356-000


8.20.1 Flush Regulator, Flush Valve, and ACGO Selector Switch



| Item | 1 | Description | Stock Number |
|------|----|---|---------------|
| 1 | | Flush valve, without button | 1006-8357-000 |
| 2 | | Flush Button with rod | 1011-3354-000 |
| 3 | | Spring | 1006-3186-000 |
| 4 | | E-clip | 0203-5225-300 |
| 5 | | Bracket | 1011-3355-000 |
| 6 | | Screw, M4x8 | 1006-3178-000 |
| 7 | | Screw, M4x12 | 0140-6226-111 |
| 8 | | Lockwasher, M4 | 9213-0540-003 |
| 9 | | ACGO Selector Switch, complete (without guard - item 10) | 1009-3099-000 |
| | 9a | Rush pressure switch | 1006-3972-000 |
| | 9b | 0-ring | 1006-3213-000 |
| | 9c | Screws | 0144-2124-201 |
| 10 | | Guard | 1009-3140-000 |
| 11 | | Tubing, silicone | 1009-3164-000 |
| 12 | | Tie wrap | 0203-5915-300 |
| 13 | | Regulator, O ₂ Flush | 1011-3168-000 |

1009-0356-000 07/04

8.21 Breathing system interface




| Item | Description | Stock Number | Qty |
|------|---|---------------|-----|
| 1 | ASSEMBLY MAIN SUPPORT CASTING | 1407-7010-000 | |
| 2 | Bolt, M6x16 flange | 1009-3125-000 | (5) |
| 3 | HANDLE GRIP | 1407-3317-000 | |
| 4 | SCR M6X16 Sems | 0144-2436-109 | (2) |
| 5 | LATCH PUSH TO CLOSE | 1407-3309-000 | |
| 6 | LATCH PUSH TO CLOSE W/MICROSWITCH | 1407-3310-000 | |
| 7 | SCR SKT HD CAP M3-0.5X8 SST | 1006-3865-000 | (2) |
| 8 | WASHER LOCK EXTERNAL M3 | 9213-0530-003 | (2) |
| 9 | Port, plug circuit | 1407-3333-000 | |
| 10 | Port, fresh gas | 1407-3314-000 | |
| 11* | SEAL U-CUP 12.7 ID BCG 19.05 OD EPR | 1407-3320-000 | (4) |
| 12 | RING RET 15.88 SHAFT DIA TYPE E SST | 1406-3446-000 | (4) |
| 13 | Port, sample gas | 1407-3318-000 | (2) |
| 14 | Connector, BULKHEAD 02 CELL, with harness | 1009-5586-000 | |
| 15 | RING RETAINING 9.53 SFT DIAMETER TYPE E SST | 1406-3277-000 | |
| 16 | SW SUBMINITURE W/QDISC TERMINALS | 1406-3296-000 | |
| 17 | SCR M2.5 X 10 | 1009-3153-000 | (2) |
| 18 | BRACKET BTV SWITCH | 1407-3319-000 | |
| 19 | LEVER BTV SWITCH | 1407-3325-000 | |
| 20 | CAP BRACKET BTV | 1407-3324-000 | |
| 21 | SCR SKT HD CAP M3-0.5X8 SST | 1006-3865-000 | (2) |
| 22 | WASHER LOCK EXTERNAL M3 | 9213-0530-003 | (2) |

^{36 07/04 1009-0356-000}


8.22 Breathing System

8.22.1 APL Valve



| Item | Description | Stock Number | QTY |
|------|---|---------------|-----|
| 1 | APL Valve Assy (includes items 2 through 6) | 1009-8200-000 | |
| 2 | SPRING CPRSN 53.14 OD 36.8 L 1.48 N/MM | 1406-3328-000 | |
| 3 | RETAINER SPRING APL | 1407-3404-000 | |
| 4 | DIAPHRAGM APL | 1406-3331-000 | |
| 5 | CAGE APL | 1406-3333-000 | |
| 6 | POPPET APL VALVE | 1406-3332-000 | |
| 7 | RAMP APL | 1407-3400-000 | |
| 8 | COVER APL | 1407-3405-000 | |
| 9 | MANIFOLD APL/BTV | 1407-3401-000 | |
| 10 | Cover, Manifold APL/BTV (with 22-mm male bag port) | 1407-3402-000 | |
| | Cover, Manifold APL/BTV (with Australian bag port - 22 mm female) | 1407-3412-000 | |
| 10 | COVER MANIFOLD APL/BTV | 1407-3402-000 | |
| 11 | WEIGHT DEAD 14CM H20 BCG ABS NEG RELIEF | 1407-3406-000 | |
| 12 | SEAL ABS NEG RELIEF VLV | 1407-3407-000 | |
| 13 | 0-RING 22 ID 30 OD 4 W SI 40 DURO | 1407-3104-000 | |
| 14 | O-RING 88.49 ID 95.55 OD 3.53 W SILICONE 50 DURO | 1407-3403-000 | |
| 15 | O-RING 1.049ID 1.2550D 0.103W EPDM NO 121 | 1407-3408-000 | |
| 16 | SCR M4X16 BT SKT HD SSTTYPE 316 | 0140-6226-115 | (2) |
| 17 | Lockwasher, M4 external | 9213-0540-003 | (2) |
| 18 | SCR THUMB M4 SHLDR 7.5 X 7 | 1407-3410-000 | (3) |
| 19 | RING RETAINING 3.96 SFT DIA CRESCENT SST | 1407-3411-000 | (3) |
| 20 | SCR M4 X 40 FL HD SST PH | 0140-6226-122 | (2) |
| 21 | O-RING 2.9 ID 6.46 OD 1.78 W EP 70 DURO | 1407-3409-000 | (2) |
| 22 | SCR SEMS M4X8 BT SKT HD W/EXT L/W SST 316 | 0144-2436-108 | (3) |

8.22.2 Bag/Vent Switch




| Item | Description | Stock Number | QTY |
|------|---|---------------|-----|
| | BTV Switch Cartridge | 1407-7003-000 | |
| 1 | COVER BTV | 1407-3500-000 | |
| 2 | SCR SEMS M4X8 BT SKT HD W/EXT L/W SST 316 | 0144-2436-108 | (2) |
| 3 | 0-RING 44.02 ID 51.1 OD 3.53 W SI 70 DURO | 1407-3507-000 | |
| 4 | SEAL, BTV | 1407-3506-000 | |

8.38 07/04 1009-0356-000


8.22.3 Absorber canister





| Item | Description | Stock Number | Qty |
|------|---|---------------|-----|
| | Absorber Canister, Reusable | 1407-7004-000 | |
| 1 | CANISTER, CO2 | 1407-3200-000 | |
| 2 | O-RING 110.72 ID 117.78 OD 3.53 W EPR 50 DURO | 1407-3204-000 | |
| 3 | FOAM, CO2 CANISTER (PKG 40) | 1407-3201-000 | |
| 4 | COVER, CO2 CANISTER with LOCKING RING (does not include items 5 and 6) | 1407-3203-000 | |
| 5 | SCREEN, CO2 CANISTER COVER | 1407-3205-000 | |
| 6 | SCREW, M3X8 FL PH HD SST | 9211-0530-083 | (2) |
| | | | |


8.22.4 Flow Sensor Module



| Item | Description | Stock Number | Qty |
|------|--|---|-----|
| | Flow Sensor Module (does not include Item 1) | 1407-7001-000 | |
| 1 | Flow Sensor (plastic) Flow Sensor (metal - autoclavable) Flow Port Adapter | 1503-3856-000 1503-3244-000 1503-3849-000 | |
| 2 | COVER FLOW SNSR | 1407-3000-000 | |
| 3 | HOLDER FLOW SNSR UPPER | 1407-3002-000 | |
| 4 | HOLDER FLOW SNSR LOWER | 1407-3003-000 | |
| 5 | SCR THUMB M6X43 SST | 1406-3304-000 | |
| 6 | SCR M4 X 10 SKT CAP BUTTON HEAD SST | 0144-2117-718 | (2) |
| 7 | CUFF FLOW SNSR | 1407-3004-000 | (2) |
| 8 | LATCH FLOW SNSR | 1407-3001-000 | |
| 9 | SPR TORSION FLOW SNSR LATCH | 1407-3005-000 | |
| 10 | RING TRUARC 0.188 SHAFT E-RING SST | 0203-5225-300 | |

8.40 07/04 1009-0356-000


8.22.5 Breathing Circuit Module



| Item | Description | Stock Number | Qty |
|------|---|---------------|-----|
| | Breathing Circuit Module (*) | 1407-7002-000 | |
| 1 | LENS CIRCUIT CHK VALVES | 1407-3101-000 | |
| 2 | 0-RING 44.02 ID 51.1 OD 3.53 W SI 70 DURO | 1407-3507-000 | (2) |
| 3 | Check Valve Assembly | 1406-8219-000 | (2) |
| 3a | RETAINER DISK 26.97D 12.7H 0.76T | 1400-3017-000 | (2) |
| Зb | DISC CHK V RVSBL 1.025D | 0210-5297-100 | (2) |
| 3c | SEAT UNIDIRECTIONAL V B/S | 1406-3396-000 | (2) |
| 3d | 0-RING 20.35 ID 23.90 OD 1.78W | 1406-3397-000 | (2) |
| 4 | PLATE CIRCUIT FLANGE | 1407-3110-000 | |
| 5 | SCR SEMS M4X8 BT SKT HD W/EXT L/W SST 316 | 0144-2436-108 | (6) |
| 6 | HOOK LATCH | 1407-3604-000 | |
| 7** | 0-RING, 22 ID 30 OD 4 W SI 40 DURO | 1407-3104-000 | (2) |
| 8** | O-RING, 12.37 ID 17.6 OD | 1006-3968-000 | |
| 9 | MANIFOLD CIRCUIT | 1407-3100-000 | |
| 10 | PIN CANISTER PIVOT | 1407-3109-000 | |
| 11 | RING TRUARC 0.188 SHAFT NO 5133-18H E-RING SST | 0203-5225-300 | (5) |
| 12 | LEVER CANISTER LATCH | 1407-3102-000 | |
| 13 | PIN CANISTER LEVER | 1407-3108-000 | |
| 14* | O ₂ Cell | 6050-0004-110 | |
| | O-ring, cell | 1406-3466-000 | |
| 15* | Plug with o-ring (for units without circuit O2 sensing) | 1503-3857-000 | |
| | O-ring, plug | 1407-3112-000 | |
| 16* | Cable, O ₂ Cell | 1009-5570-000 | |

* The $\rm O_2$ cell (or plug) and the cell cable are not included in the breathing circuit module. ** Lubricate sparingly with Krytox.


8.22.6 Exhalation valve



| Item | Description | Stock Number | Qty |
|---------|--|---------------|-----|
| | Exhalation Valve Assy | 1407-7005-000 | |
| 1 | BASE EXHALATION VALVE | 1407-3700-000 | |
| 2 | DIAPHRAGM ASSY EXH VALVE | 1503-8121-000 | |
| 3 | SEAT EXHALATION VLV ABS | 1407-3704-000 | |
| 4 | COVER EXHALATION VALVE | 1407-3701-000 | |
| 5 | SCR M4X16 PH PAN HD SSTTYPE 316 | 9211-0440-163 | (3) |
| 6 | O-RING 2.9 ID 6.46 OD 1.78 W EP 70 DURO | 1407-3409-000 | (3) |
| 7 | SCR THUMB M6X43 10MM HEAD B/S | 1406-3306-000 | (2) |
| 8 | O-RING 4.47 ID X 8.03 OD 1.78 W EPR 70 DURO | 1407-3703-000 | (2) |
| 9 | RETAINER DISK 26.97D 12.7H 0.76T SST FLUTTER V | 1400-3017-000 | |
| 10 | WEIGHT DEAD 10 CMH20 BCG PASSIVE AGSS | 1406-3572-000 | |
| 11 | SEAT POSITIVE PRESS BCG VALVE PASSIVE AGSS | 1406-3571-000 | |
| 12* | O-RING OD 19.16 BCG ID 15.6 EPDM DURO 70 -016 | 1006-3616-000 | |
| 13 | RING RETAINING 19.05 SHAFT DIA SST | 1406-3577-000 | |
| * Lubri | cate sparingly with Krytox. | | |


8-42 07/04 1009-0356-000

8.22.7 Bellows



| Item | Description | Stock Number |
|------|----------------------------|------------------------|
| 1 | Bellows housing | 1500-3117-000 |
| 2 | Bellows | 1500-3378-000 |
| 3 | Rim | 1500-3351-000 |
| 4 | Pressure relief valve assy | 1500-3377-000 |
| 5 | Latch, base | 1500-3352-000 |
| 6 | Seal, base | 1500-3359-000 |
| 7 | Base, bellows | Refer to section 8.22. |
| 8 | Manifold, bellows base | 1407-3702-000 |


8.22.8 Bellow base



| Item | Description | Stock Number |
|------|-------------------|---------------|
| 1 | Bellows Base Assy | 1407-7006-000 |
| 1 | a Latch Assy | 1407-7007-000 |
| 2 | HOOK LATCH | 1407-3604-000 |
| 3 | E-Ring | 0203-5225-300 |
| 3 | | |


8.44 07/04 1009-0356-000

8.22.9 Bag Arms



| Item | 1 | Description | Stock Number | Qty |
|------|-------|-----------------------------------|---------------|-----|
| 1 | | Bag Arm Assembly (complete) | 1009-8159-000 | |
| 2 | | Bag Arm Upper Assembly | 1407-7011-000 | |
| - 3 | 2a | Cover, bag port housing | 1407-3807-000 | |
| | | Screw, M3x20 | 0140-6719-103 | |
| | | Lockwasher, M3 internal | 9213-0430-003 | |
| 3 | 2b | Housing, bag port | 1407-3806-000 | |
| | 2c | Lever, lock release | 1407-3808-000 | |
| 3 | 2d | Ring, retaining | 1406-3577-000 | |
| - 33 | 2e | Nut, M3 Nyloc | 0144-3536-112 | |
| 3 | | Bag Arm Lower Assembly | 1407-7012-000 | |
| | 3a | PAD FRICTION MATERIAL | 1407-3818-000 | |
| - 33 | 3b | SPRING CPRSN | 1406-3270-000 | |
| 33 | 3с | WASHER, shoulder | 1407-3815-000 | |
| - 3 | 3d | WASHER LOCK M3 INT SST | 9213-0430-003 | (2) |
| - 8 | 3e | SCREW M3X16 POSI DR PAN HD A4 SST | 1504-3003-000 | (2) |
| | 3f | Nut, M5 Nyloc | 9212-0350-006 | |
| 4 | | Pin, dowel 3.18 DIA 31.8 L SST | 1407-3804-000 | |
| Item | ıs If | no Bag Arm | | |
| 5 | | Clip, patient tubing | 1407-3810-000 | |
| 6 | | Washer, shoulder | 1407-3814-000 | |
| 7 | | Lockwasher, M4 external | 9213-0540-003 | |
| 8 | | Screw, M4x16 | 9211-0440-163 | |


8.23 Drawer



| Item | Description | Stock Number |
|------|----------------------------------|---------------|
| 1 | Slide, drawer | 1009-3084-000 |
| 2 | Screw, M4x8 Nylac | 1009-3183-000 |
| 3 | Drawer, body | 1009-3078-000 |
| 4 | Drawer Front, lower (down arrow) | 1009-3032-000 |
| 5 | Drawer Front, upper (up arrow) | 1009-3031-000 |
| 6 | Screw, M4x12 | 1009-3109-000 |

8-46 07/04 1009-0356-000

8.24 Legris quick-release fittings




Note: Not every fitting is used in all machines.

8.25 Vent Drive and low-pressure tubing

| 6700-0005-300 | 180 mm - 1/2 inch | Typon | unmarked | 150 |
|---------------|-------------------|----------------|------------------------------------|----------|
| 1001-3063-000 | 40 mm - 8 mm | | unmarked | 149 |
| 6700-0005-300 | 465 mm - 1/2 inch | Tygon | unmarked | 148 |
| 6700-0005-300 | 260 mm - 1/2 inch | Tygon | unmarked | 146 |
| 1009-3296-000 | 300 mm - 8 mm | (black) | Venturi Drive | 145 |
| 1009-3363-000 | 330 mm - 4 mm | (black) | Venturi Pilat | 144 |
| 1605-1001-000 | 300 mm - 1/4 inch | (low-pressure) | RGM to Circuit | 137 |
| 1605-1001-000 | 200 mm - 1/4 inch | (low-pressure) | RGM to Scavenge | 136 |
| 1605-1001-000 | 50 mm - 1/4 inch | (low-pressure) | unmarked | 135 |
| 1605-1001-000 | 25 mm - 1/4 inch | (low-pressure) | unmarked | 134 |
| 1605-1001-000 | 330 mm - 1/4 inch | (low-pressure) | PAW | 133 |
| 1605-1001-000 | 260 mm - 1/4 inch | (low-pressure) | PAW | 132 |
| 1605-1001-000 | 600 mm - 1/4 inch | (low-pressure) | PAW | 131 |
| 1605-1001-000 | 750 mm - 1/4 inch | (low-pressure) | AGSS flowtube | 130 |
| 1605-1001-000 | 151 mm - 1/4 inch | (low-pressure) | unmarked | 129 |
| 1605-1001-000 | 300 mm - 1/4 inch | (low-pressure) | unmarked | 128 |
| 1605-1001-000 | 750 mm - 1/4 inch | (low-pressure) | RGM return | 127 |
| 1605-1001-000 | 250 mm - 1/4 inch | (low-pressure) | Aux 02 0UT | 120 |
| 1009-3296-000 | 900 mm - 8 mm | (black) | VENT DRIVE | 125 |
| | Length - Size | ld only) | Tube Markings (factory build only) | |
| 1009-3164-000 | 42 mm - 3/8 inch | | Tubing (silicone) | 13 |
| 1009-3164-000 | 72 mm - 3/8 inch | | Tubing (silicone) | 12 |
| 1605-1001-000 | 151 mm - 1/4 inch | | Tubing, low-pressure | 11 |
| 1006-3530-000 | | | Plug, 4-mm | 10 |
| 1009-3077-000 | | | Fitting, coupler barb ends | 9 |
| 1406-3524-000 | | | Cap, plug | 00 |
| 1009-3011-000 | | | Tee (male barb) | 7 |
| 1503-3131-000 | | | Coupler, male - yellow | Ø |
| 1503-3132-000 | | | Coupler, female - yellow | ຜ |
| 1503-3236-000 | | | Coupler, male - white | 4 |
| 1503-3119-000 | | | Coupler, female - white | ω |
| 1503-3237-000 | | | Coupler, male - black | 2 |
| 1503-3128-000 | | | Coupler, female - black | <u>ц</u> |
| Stock Number | Length - Size | | Description | Item |

8-48

Refer to section 8.26
 Sample gas return is directed to the scavenging system as a factory default. A qualified service representative can reroute the sample gas back to the breathing system (refer to Section 4.13).



8.26 Tubing for use with Legris fittings


Except for the Tygon tubing (Nems 147 and 150), this tubing is a flexible, Myon-type tubing for use with quick-release fittings.

| *901 | N20 PLINE - N20 FLWMOD | 1200 mm - 4 mm | 1001-3060-000 |
|--------|------------------------|-------------------|---------------|
| 107* | N20 PLINE - N20 CYL | 430 mm - 4 mm | 1001-3060-000 |
| 108** | UNMARKED | 40 mm - 4 mm | 1001-3060-000 |
| *601 | N20 CYL - N20 FLWMOD | 800 mm - 4 mm | 1001-3060-000 |
| 124** | Unmarked | 250 mm - 4 mm | 1001-3060-000 |
| 126 | Unmarked | 50 mm - 4 mm | 1001-3060-000 |
| | SW4-02 FLWMOD | 1350 mm - 6 mm | 1001-3062-000 |
| 114 | REGULATED 02 | 400 mm - 6 mm | 1001-3062-000 |
| | REGULATED 02 | 600 mm - 6 mm | 1001-3062-000 |
| | REGULATED 02 | 250 mm - 6 mm | 1001-3062-000 |
| | 02 PLINE - REG IN | 330 mm - 6 mm | 1001-3062-000 |
| 118 | 02 PLINE - 02 SW3 | 270 mm - 6 mm | 1001-3062-000 |
| 123*** | unmarked | 175 mm - 6 mm | 1001-3062-000 |
| 219 | 02 CYL-02 PLINE | 215 mm - 6 mm | 1001-3062-000 |
| | AIR PLINE - AIR SW3 | 460 mm - 8 mm | 1001-3063-000 |
| 110 | SW4-AIR FLWMOD | 1350 mm - 8 mm | 1001-3063-000 |
| | AIR CYL- AIR PLINE | 270 mm - 8 mm | 1001-3063-000 |
| 100 | AIR PLINE - AIR GAGE | 470 mm - 1/8 inch | 1006-3718-000 |
| | N20 PLINE - N20 GAGE | 470 mm - 1/8 inch | 1006-3718-000 |
| 102 | 02 PLINE - 02 GAGE | 470 mm - 1/8 inch | 1006-3718-000 |
| 104 | VAP OUT- ACGO | 840 mm - 1/4 inch | 1001-3064-000 |
| 105 | FLWMOD-VAP IN | 840 mm - 1/4 inch | 1001-3064-000 |
| | FLUSH VLV-ACGO | 280 mm - 1/4 inch | 1001-3064-000 |
| | unmarked (Tygon) | 290 mm - 1/2 inch | 6700-0005-300 |
| 150 | | | |

07/04 1009-0356-000

9.80


With no N₂ O cylinder supply, items 107 and 109 are replaced with Item 106.
 With an inboard N₂O cylinder supply, Items 124 is replaced with Item 108.
 With two inboard O₂ cylinder supplies, Item 123 connects the second O₂ cylinder supplies, Item 123 connects the second O₂ cylinder supply.



8.27 Cables and harnesses

| Item | Description | Stock Number |
|------|--|--------------------------------|
| 1 | Harness, Vent Monitoring board to ABS flow sensors (includes tubing) | 1009-8165-000 |
| 2 | Cable, Vent Monitoring board | 1504-5604-000 |
| 3 | Cable, Serial Isolation | 1009-5691-000 |
| 4 | Cable, Pneumatic Vent Engine | 1504-5605-000 |
| 5 | Cable, power | 1009-5711-000 |
| 6 | Harness, Serial ISO to O2 flush SW | 1009-5567-000 |
| 7 | Harness, Serial ISO to 02 supply SW | 1009-5568-000 |
| 8 | Harness, Serial ISO to on/ standby SW | 1009-5566-000 |
| 9 | Harness, O2 CELLTO FTLR BRD HARN | 1009-5586-000 |
| 10 | Harness, BAG TO VENT SW TO HARN FILR BRD | 1009-5585-000 |
| 11 | Harness, FLTR BRD TO B/S 02 SNSR AND SW | 1009-5531-000 |
| 12 | Harness, VENT ENG BRD TO CONN PLATE | 1009-5545-000 |
| 13 | Harness, Vent Mon Brd to Task Light | 1009-5533-000 |
| 14 | Harness, to 100/120 V outlets Harness, to 220/240 V outlets | 1009-5716-000 1009-5717-000 |
| 15 | Harness, ACGO Switch | 1009-5762-000 |
| 16 | Harness, Task Light | 1009-5584-000 |
| 17 | Power Cord | Refer to section 8.14 |
| | | |


8-52 07/04 1009-0356-000



8.28 Cables and harnesses (Aespire 100)

| Item | Description | Stock Number |
|------|--|-----------------------|
| 1 | Harness, Vent Monitoring board to ABS flow sensors (includes tubing) | 1009-8165-000 |
| 2 | Cable, Vent Monitoring board | 1504-5604-000 |
| 3 | Cable, 7100 module to machine switches | 1009-6060-000 |
| 4 | Cable, Pneumatic Vent Engine | 1504-5605-000 |
| 5 | Cable, power to Control module | 1009-5711-000 |
| 6 | Power Cord, machine | Refer to section 8.14 |
| 7 | Harness, ACGO Switch | 1009-5762-000 |
| 8 | Harness, 02 CELL TO FTLR BRD HARN | 1009-5586-000 |
| 9 | Harness, BAG TO VENT SW TO HARN FTLR BRD | 1009-5585-000 |
| 10 | Harness, FLTR BRD TO B/S 02 SNSR AND SW | 1009-5531-000 |
| 11 | Harness, VENT ENG BRD TO CONN PLATE | 1009-5545-000 |

8 54 07/04 1009-0356-000




8.29 Anesthetic Gas Scavenging System — AGSS

8.29.1 Passive AGSS Items 1 through 12 are included in all AGSS kits.

| Item | Description, Common Parts | Stock Number | Qty |
|---------|--|---------------|----------|
| 1 | Seal, Receiver Body | 1407-3901-000 | |
| 2 | Reservoir | 1407-3903-000 | |
| 3 | Seal and scavenging down-tube | 1407-3904-000 | |
| 4 | Thumbscrew, M6x28.5 | 1406-3305-000 | |
| 5 | 0-ring, 4.42 ID, 9.65 0D | 1407-3923-000 | (2) |
| 6 | Thumbscrew, M6x43 | 1406-3304-000 | |
| 7 | Valve, unidirectional (negative pressure relief) | 1406-8219-000 | |
| 7a | Seat, Valve, Negative Pressure | 1406-3396-000 | |
| 7b | Retainer, disc | 1400-3017-000 | |
| *7c | 0-ring, 20.35 ID, 23.90 0D | 1406-3397-000 | |
| 7d | Disc, check-valve | 0210-5297-100 | |
| 8* | O-ring, 22 ID, 30 OD silicone | 1407-3104-000 | (2) |
| 9* | 0-ring, 21.95 ID, 25.51 OD | 1406-3558-000 | |
| 10 | Screw, M4x8 | 9211-0640-083 | (2) |
| 11 | Cap, 3.18 Barb, Silicone | 1406-3524-000 | |
| 12 | Adapter, auxiliary inlet, 30-mm male to 30-mm male | M1003134 | |
| 13 | Adapter, auxiliary inlet, 30-mm male to 19-mm male | M1003947 | |
| Passive | AGSS Specific Parts | | |
| 14 | Receiver, Passive/Adjustable | 1407-3908-000 | |
| 15 | Plug Assembly, tethered | 1407-3909-000 | (2) |
| 16 | Screw, shoulder M3 | 1407-3915-000 | |
| 17 | Connector, 30-mm ISO, Male | 1406-3555-000 | |
| 18 | Adapter, scavenging, 30-mm female to 19-mm male | 1500-3376-000 | (5 pack) |

^{*} Lubricate sparingly with Krytox

8-56 07/04 1009-0356-000




8.29.2 Adjustable AGSS Items 1 through 12 are included in all AGSS kits.

| Item | Description, Common Parts | Stock Number | Qty |
|---------|--|---------------|-----|
| 1 | Seal, Receiver Body | 1407-3901-000 | |
| 2 | Reservoir | 1407-3903-000 | |
| 3 | Seal and scavenging down-tube | 1407-3904-000 | |
| 4 | Thumbscrew, M6x28.5 | 1406-3305-000 | |
| 5 | O-ring, 4.42 ID, 9.65 OD | 1407-3923-000 | (2) |
| 6 | Thumbscrew, M6x43 | 1406-3304-000 | |
| 7 | Valve, unidirectional (negative pressure relief) | 1406-8219-000 | |
| 7a | Seat, Valve, Negative Pressure | 1406-3396-000 | |
| 7b | Retainer, disc | 1400-3017-000 | |
| *7c | 0-ring, 20.35 ID, 23.90 OD | 1406-3397-000 | |
| 7d | Disc, check-valve | 0210-5297-100 | |
| 8* | O-ring, 22 ID, 30 OD silicone | 1407-3104-000 | (2) |
| 9* | O-ring, 21.95 ID, 25.51 OD | 1406-3558-000 | |
| 10 | Screw, M4x8 | 9211-0640-083 | (2) |
| 11 | Cap, 3.18 Barb, Silicone | 1406-3524-000 | |
| 12 | Adapter, auxiliary inlet, 30-mm male to 30-mm male | M1003134 | |
| 13 | Adapter, auxiliary inlet, 30-mm male to 19-mm male | M1003947 | |
| Adjusta | ble AGSS Specific Parts | | |
| 14 | Receiver, Passive/Adjustable | 1407-3908-000 | |
| 15 | Plug Assembly, tethered | 1407-3909-000 | |
| 16 | Screw, shoulder M3 | 1407-3915-000 | |
| 17 | Needle Valve Assembly (with DISS EVAC connector) | 1407-3918-000 | |
| 18 | Bag with 30 mm male connector | 8004460 | |
| | | | |


^{*} Lubricate sparingly with Krytox

8-58 07/04 1009-0356-000




8.29.3 Active AGSS tens 1 through 12 are included in all AGSS kits.

| 1 Season of the content of the c | | 1407-3903-000 |
|--|-----------------------------|---------------|
| Research Research Research Research Research Functionary ACLD, 9, 625 00 Thumbscrew, M6028, 5 O-fring, ACLD, 9, 625 00 Thumbscrew, M6028, 5 Seat, Mark, Neighthe Pressure ACLD, 9, 625 00 To Disc, cheek-valve O-fring, 22, 10, 30, 00 silicone O-fring, 21, 10, 20, 00 silicone O-fring, 21, 10, 20, 00 silicone O-fring, 21, 10, 20, 20, 51, 10 Server, MAAAS Cap, 31, 8, 84h, 8, 81cone Adapter, auxiliary inlet, 30-mn Adapter, auxiliary inlet, 30-mn Adapter, auxiliary inlet, 30-mn Adapter, auxiliary inlet, 30-mn Adapter, with air brake Seal, for filter and orifice Filter Commector, lowflow WAC Ommector, lowflow WAC OMMECTOR WAC OM | | 3903-000 |
| Seal and savereging down-tail Thumbscrew, M6028.5 O-ring, 4.42.D, 9.65.00 Thumbscrew, M6028.5 O-ring, 4.42.D, 9.65.00 Thumbscrew, M6028.5 Valve, undirectional (regative Valve) (regative | | 20000 |
| Dumbacrae, M6028.5 Oning AZ ID, 9, 605 00 Thumbacrae, M6043 Valve, unidirectional (regative Pressurary Ta Seat, Valve, Negative Pressurary To Paraliner, disc. Oning 20,35 ID, 23,00 00 To Oning 21,95 ID, 25,51 00 Screw, MAAS Cap, 3,18 Barb, Silcone Adapter, auxiliary inlet, 30-nn Onifice, low flow Onifice, low flow Oninector, low flow 12-n mn Onifice, low flow Oninector, low flow 12-n mn Adapter, 12-nn Adapter, 12-n | 1406 | 1407-3904-000 |
| Punite and A.2 (D. 9.65 0D) Punite and A.6 (D. 9.65 0D) Punite and A.6 (D. 9.65 0D) Punite and orifice and orifice and A.6 (D. 9.65 0D) Punite and orifice and orifice and A.6 (D. 9.65 0D) Punite and Orifice and Orifice and A.6 (D. 9.65 0D) Punite and Orifice and Orifice and Orifice and A.6 (D. 9.65 0D) Punite and Orifice and A.6 (D. 9.65 0D) Punite and Orifice and Orifice and Orifice and A.6 (D. 9.65 0D) Punite and A.6 (D. 9.65 0D) Pu | | 1406-3305-000 |
| Phumbscraw, Midod Walker, Walker Walker, Walker Walker, Walker Wa | 1407- | 1407-3923-000 |
| Asset, windirectional (negative Pressure) 72 Seat, Nake, Negative Pressure) 73 Desiried, disc. 74 O-ring, 20.35 (D. 23.90 00 75 Disc, Observative Control of Server, MAAS 85 Server, MAAS 87 Strong, 315 Berh, Silicone Adapter, auxiliary inlet, 30-mn Adapter, auxiliary inlet, 30-mn Adapter, auxiliary inlet, 30-mn Adapter, with air brake 88 Specific Parts 88 Receiver, with air brake 88 Herri of Therr and orifice Filter Connector, Ingh flow M30 thre Orifice, Ingh flow M30 thre Orifice, Ingh flow W40 thre Orifice, by flow Ne Low Flow with 25 mn connector Confice, to w flow Orifice, low flow 10 Connector, tow flow 25 mn Orifice, low flow 10 Connector, tow flow 45 mn Orifice, low flow 10 Connector, tow flow 45 mn Orifice, low flow 10 Connector, tow flow 25 mn Orifice, low flow 10 Connector, tow flow 25 mn Orifice, low flow 10 Connector, tow flow 25 mn Orifice, low flow 10 Connector, tow flow 25 mn Orifice, low flow 10 Connector, tow flow 25 mn Orifice, low flow 10 Connector, tow flow 25 mn Orifice, low flow 10 Connector, tow flow 25 mn Orifice, low flow 10 Connector, tow flow 25 mn Orifice, low flow 10 Connector, tow flow 25 mn Orifice, low flow 10 Connector, tow flow 27 mn hose ba | 1406 | 1406-3304-000 |
| 7a Seat, Vahve, Negathe Pressure 7b Chrig 20.35 (D. 23.90 00) 7d Disc, check-valve Chrig 22 (D. 23.90 00) Chrig 21.95 (D. 25.51 00) Srew, MAAS Cap, 318 Barb, Silnone Adapter, audiliary inlet, 30-nm Adapter, 30-nm Adapt | | 1406-8219-000 |
| The Restation, cliac The Coning, 22 10, 23 00 os silicone O-ring, 22 10, 23 00 os silicone O-ring, 22 10, 25.51 00 Screw, MAAA Cap, 318 Barb, Silicone Adapter, auxiliary inlet, 30-rm Adapter, inlet parts Receiver, with air brake Seed, for filter and orfice Filter Connector, light flow M30 thre Orfice, inligh flow Ordice, inligh flow Connector, tow flow Ordice, two flow Connector, two flow 25 rm Ordice, two flow Connector, two flow 127 rm hose ba Connector, two flow 127 rm hose ba Connector, two flow 127 rm hose ba | | 1406-3396-000 |
| 7d Desc, Debeck-value 7d Desc, Debeck-value O-ring, 21 D, 30 Do silicone O-ring, 21 D, 30 Do silicone O-ring, 21 B, 8th, 8th, 8th, 8th, 8th, 8th, 8th, 8th | | 1400-3017-000 |
| 7d Disc, check-valve O-ring, 22 Ib, 330 Osilicone O-ring, 21 BB Ib, 25 Ib, 50 OS Serew, MAd8 Cap, 31 B Barb, Silicone Adapter, audiliary inte, 30-rm Connector, Ior Row EMC connector S Connector, Ior Row EMC Connector, Ior Row EMC Office, low flow Office, low flow Office, low flow 12.7 mm hose be Connector, Iow flow 12.7 mm hose be Connector, Iow flow 12.7 mm long- | 1406 | 1406-3397-000 |
| O-ring, 22 ID, 30 OO silicone O-ring, 21 IB, 105 ID, 28-51 OD Screw, MAAA Cap, 3.18 Barb, Silicone Adapter, audiliary inlet, 30-mm Connector, inlet flow M30 three Connector, low flow Onfitice, low flow Connector, low flow to 12.7 mm hose be Connector, low flow to 12.7 mm hose be Connector, low flow to 12.7 mm hose be | 0210- | 0210-5297-100 |
| Sorwi, MAAS Cap, 3.18 Barb, Silcone Adapter, auxiliary inlet, 30-nn AdasS Specific Parts Seal, for filter and onifice Filter Connector, high flow M30 thre Onifice, high flow Connector, for flow EVAC Connector, to vifew EVAC Connector, to vifew EVAC Onifice, but flow Ne Low Flow with 25 mm connector Confice, to vifew Ne Low Flow with 25 mm connector Onifice, but flow Ne Low Flow with 127 mm hose ba Connector, low flow 150 mm Onifice, low flow Connector, low flow 150 mm Onifice, low flow Connector, low flow 150 mm Onifice, low flow Connector, low flow 120 mm Onifice, low flow Connector, low flow 120 mm Onifice, low flow | 1407- | 1407-3104-000 |
| Screw, MAAS Cap, 318 Barb, Sillicone Adapter, auxiliary inite, 30-rm AdaSS Specific Parts Seal, for fitner and orifice Connector, losh flow M30 thre Orifice, light flow AND Connector, low flow EMC Orifice, low flow Connector, low flow flow 127 mm hose be Connector, low flow flow 127 mm hose be Connector, low flow flow 127 mm hose be Connector, low flow 127 mm hose be Connector, low flow 127 mm hose be | 1406 | 1406-3558-000 |
| Adapter, audilary inlet, 30-mn Adapter, audilary inlet, 30-mn Adapter, audilary inlet, 30-mn Adapter, audilary inlet, 30-mn AMS Specific Parts Receiver, with air trake Seed, for filter and orfice Filter Connector, flow M30 thro Outlice, high flow Outlice, in figh flow Connector, flow Flow Office, by flow Office, by flow Office, by flow Office, low flow to the ba | 9211- | 9211-0640-083 |
| Adapter, auxiliary inlet, 30-nn Adapter, auxiliary inlet, 30-nn Adapter, auxiliary inlet, 30-nn Receiver, with air brake Seal, for filter and onfice Filter Connector, high flow M30 thre Onfines, high flow M30 thre Onfines, to flow EVAC Connector, for flow EVAC Connector, to wflow EVAC Onfines, to wflow Onfines, to wflow Onfines, to wflow Onfines, to wflow Onfines, tow flow Onfines, tow flow Onfines, tow flow Onfines, tow flow Connector, tow flow 25 mm Onfines, tow flow Connector, tow flow 25 mm Onfines, tow flow | 1406 | 1406-3524-000 |
| Adapter, auxiliary inter, 30-mm Receiver, with air brake Receiver, with air brake Seal, for filter and orifice Filter Connector, light flow M30 thra Orifice, hight flow Connector, low flow EMAC Orifice, low flow Connector, low flow EMAC Orifice, low flow Connector, low flow Connector, low flow Connector, low flow Connector, low flow 25 mm Orifice, low flow Connector, low flow 12.7 mm hose ba Connector, low flow 12.7 mm hose ba | | M1003134 |
| he AGSS Specific Parts Receiver, with air brake Seal, for filter and orifice Filter Connector, high flow M30 through Connector, high flow M30 through Connector, low flow EMC Connector, low flow EMC Confice, low flow Connector, low flow Connector, low flow Connector, low flow Connector, low flow 12,7 mm hose be | | M1003947 |
| Receiver, with air brake Seal, for filter and orifice Filter For filter and orifice Filter Connector, high flow M30 the Orifice, high flow M30 the Orifice, high flow Connector, low flow W4C Orifice, low flow Connector, low flow 12 m accomedite Connector, low flow 12 m base be Ne Low Flow with 12 m m base be Connector, low flow 12 mm hase be Connector, low flow 12 mm hase be Connector, low flow 12 mm hase be | | |
| Seal, for fitner and orifice Filter Connector, light flow M30 three Orifice, hight flow Townector, low flow BVAC Connector, low flow BVAC Office, but flow Orifice, low flow Connector, low flow BVAC Office, low flow Connector, low flow 25 mm Orifice, low flow Connector, low flow 25 mm Orifice, low flow Connector, low flow 12.7 mm hose be Connector, low flow 12.7 mm long- | 1407 | 1407-3900-000 |
| Filter Connector, high flow M30 three Onflice, high flow Onflice, high flow Connector, flow the W4C Onflice, by flow Connector, tow flow EVAC Onflice, by flow Ne Low Flow with 25 mm connector Connector, low flow 25 mm Onflice, low flow Connector, low flow 12 mm hose ba Connector, low flow 12 mm hose ba Connector, low flow 12.7 mm hose ba | 1407 | 1407-3902-000 |
| No High Flow Spacific Parts Connector, high flow M30 three Onflice, high flow Onflice, high flow Connector, low flow EVAC Onflice, low flow EVAC Onflice, low flow 25 mm Onflice, low flow 22 mm hose ba Connector, low flow 12.7 mm hose ba | 1406 | 1406-3521-000 |
| Connector, high flow M30 three Offices, high flow Ne Lew Flow with EVMC connector S Connector, low flow EVMC Offices, low flow Connector, low flow EVMC Connector, low flow 25 mm Onfices, low flow Ne Lew Flow with 12.7 mm hose ba Connector, low flow 12.7 mm longer Connector, low flow 12.7 mm longer | | |
| Orifice, high flow Connector, low flow EVAC Orifice, low flow EVAC Orifice, low flow 25 mm connector, low flow 25 mm Connector, low flow 25 mm Orifice, low flow 12.7 mm hose ba Connector, low flow 12.7 mm hose ba Connector, low flow 12.7 mm hose ba I Connector, low flow 12.7 mm hose ba | | 1406-3557-000 |
| he Low Flow with EVIAC connectors S Connector, low flow EVIAC Orifice, low flow Connector, low flow 25 mm Orifice, low flow 25 mm Orifice, low flow Connector, low flow Connector, low flow 1 Connector, low flow 12.7 mm hose ba | 1407- | 1407-3920-000 |
| Connector, low flow BVAC Oriflee, low flow Pow Connector, low flow 25 mm Oriflee, low flow 25 mm Oriflee, low flow Connector, low flow Connector, low flow 12.7 mm hose ba Connector, low flow 12.7 mm lone | pecific Parts | |
| Orifice, low flow Connector, low flow 25 mm Orifice, low flow 25 mm Orifice, low flow 12.7 mm hose be Connector, low flow 12.7 mm 1-none- | | 1406-3597-000 |
| he Low Flow with 25 mm connection Connector, low flow 25 mm Orifice, low flow Ne Low Flow with 12.7 mm hose be Connector, low flow 12.7 mm 1-none- | 1407- | 1407-3919-000 |
| Oranector, low flaw 25 mm Orifice, low flow Ne Low Flow with 12.7 mm hose be Connector, low flow 12.7 mm -none- | Specific Parts | |
| Orifice, low flow. Ne Low Flow with 12.7 mm hose be Connector, low flow 12.7 mm -none- | | 1406-3573-000 |
| Ne Low Flow with 12.7 mm hose bar Connector, low flow 12.7 mm inne-none- | 1407- | 1407-3919-000 |
| Connector, low flaw 12.7 mm (| th connector Specific Parts | |
| | 1/2 inch) 1406- | 1406-3574-000 |
| | | |
| Active Low Flow with 30 mm ISO male connector Specific Parts | connector Specific Parts | |
| 17e Connector, low flow 25 mm | 1406 | 1406-3555-000 |
| 18 Orifice, low flow | 1407- | 1407-3919-000 |
| * I wholeste energinely with Krytov | | |
| duricate spallingly with hybrid | | |



8.30 Integrated Suction Regulator


8.30.1 Major Components (Continuous and Venturi suction)



| Item | Description | Stock Number |
|------|----------------------------------|-------------------------|
| 1 | Suction Control Module | Refer to section 8.30.2 |
| 2 | Venturi Assembly | Refer to section 8.30.3 |
| 3 | Cover, blank (if no Suction) | 1009-3271-000 |
| 4 | Bracket, blank cover mounting | 1009-3270-000 |
| 5 | Screw, M4x10 self-tapping | 1009-5534-000 |
| 6 | Manifold | 1009-3123-000 |
| 7 | Screw, M5x20 BHSCS PTTHD FORMING | 1009-3384-000 |
| 8a | Connector, NIST | 1011-3524-000 |
| 8b | Connector, Barb | 0221-0702-300 |
| 8c | Connector, Air Liquide | 1009-8292-000 |
| 8d | Muffler, for Venturi Drive | 1011-3511-000 |
| 9 | Coupling, Colder insert metal | 1009-3135-000 |
| 10 | Coupling, Colder body black | 1009-3373-000 |
| 11 | Coupling, Colder insert black | 1009-3374-000 |
| 12 | Coupling, Colder body white | 1009-3371-000 |
| 13 | Coupling, Colder insert white | 1009-3372-000 |
| 14 | Tubing, Tygon | Refer to section 8.26 |
| 15 | Fitting, barb to 8-mm Lagris. | 1009-3137-000 |
| 16 | Cap, white | 1009-3385-000 |
| 17 | Overflow Safety Trap | 6700-0647-000 |

8-62 07/04 1009-0386-000

8.30.2 Suction Control Module




| Item | Description | Stock Number |
|---------------|--|---|
| 1 1a | Gauge, 760 mmHg Gauge, 1 Bar O-ring, Gauge (included with gauge assy, 2ea. required) | 1009-3227-000 1009-3228-000 6700-0133-500 |
| 2 | Control panel assembly, with suction regulator knob and mode control knob | 1009-3213-000 |
| 3 3a 3b | | 6700-1225-800 6700-0136-500 0210-0527-300 |
| 4 | Manifold Assembly, without Gauge and Regulator Module | 1009-3277-000 |
| 5 | Screw, #6 - 2 inch | 1009-3340-000 |
| 6 | Mounting bracket | 1009-3255-000 |
| 7 | Screw, #6 - 1 inch | 1009-3339-000 |
| 8 | Filter | 0206-5159-300 |
| 9 | Pilot valve adapter assembly (includes plunger, jam nut, and valve assembly) | 1009-3278-000 |
| 10 | Cap, white | 1009-3192-000 |
| | | |

^{**} Lubricate the regulator module o-rings and the mating bore of the manifold sparingly with Dow 111 lubricant.

** Drop the plunger (9a), round end first, into the manifold. Thread the pilot valve into the manifold body. Set the mode switch to raise the plunger. Adjust the pilot valve (9b) so that the plunger actuates the pilot valve approximately half of its travel. Tighten the jam nut (9e).

8.30.3 Venturi assembly




| Item | Description | Stock Number |
|------|----------------------------------|-----------------------|
| 1 | C-clip retainer, Truarc | 1500-3158-000 |
| 2 | Elbow fitting, 4-mm Legris | 1006-3663-000 |
| 3 | Cap | 1011-5002-000 |
| 4 | Spoppet | 1011-5001-000 |
| 5 | Seal, u-cup large | 1503-3090-000 |
| 6 | Orifice | 1011-3508-000 |
| 7 | Screen, 150 mesh monel | 1001-3808-000 |
| 8 | Seal, u-cup small | 1503-3089-000 |
| 9 | Body | 1011-5000-000 |
| 10 | Venturi | 1011-3509-000 |
| 11 | Elbow fitting, 8-mm Legris | 1011-3510-000 |
| 12 | O-ring, large | 9221-3032-116 |
| 13 | O-ring, small | 1503-3108-000 |
| 14 | Check valve | 1011-8002-000 |
| 15 | Bracket, Venturi mounting | 1009-3182-000 |
| 16 | Screw, M5x20 BHSCS PTTHD FORMING | 1009-3384-000 |
| 17 | Cable tie | 0203-5915-300 |
| 18 | Tubing, Tygon | Refer to section 8.25 |

8-64 07/04 1009-0356-000

8-65

8.31 Auxiliary O₂ Flowmeter




| Item | Description | Stock Number |
|------|--|--------------------------------|
| 1 | Flowmeter, 1-10 L/min, Complete with fittings installed Flowmeter, 1-10 L/min, without fittings | 1006-8424-000 1006-3841-000 |
| 2 | Knob, gray | 1011-3471-000 |
| 3 | Set Screw | 9211-0830-053 |
| 4 | Nipple, Panel-Mount, Auxiliary 02 Outlet | 1006-5177-000 |
| 5 | Label, blank (if no Auxiliary O ₂) | 1009-3243-000 |
| 6* | Nut, M12x1.75, SST | 0144-3132-140 |
| 7** | Flowmeter Fitting, 1/8 NPTM straight adapter Flowmeter Fitting, 1/8 NPTM elbow adapter | 0204-8877-300 0204-8788-300 |
| 8** | Flowmeter Fitting Assembly, 6-mm Tubing Adapter | 1006-8423-000 |
| 9 | Plate, Flowmeter Mounting | 1009-3126-000 |
| 10 | Screw, 10-32 x 3/8 | 0140-6631-107 |
| 11 | Screw, M4x10 self-tapping | 1009-5534-000 |
| 12 | Tubing (low-pressure) 250 mm - 1/4 inch | 1605-1001-000 |
| 13 | Plate, blank (if no Auxiliary O ₂) | 1009-3128-000 |

^{*} Apply Loctite 242.** Apply Teflon tape.

1009-0356-000 07/04


8.32 Display mounts



| Item | Description | Stock Number |
|-------|--|---------------|
| 1 | GCX Arm Kit, long | 1009-3262-000 |
| 2 | GCX Bracket Kit, includes display mount (2a), cable guard (2b), and mounting hardware | 1009-3263-000 |
| 3 | GCX Arm Kit, short includes display arm (3a) and fasteners (3b) | 1009-3264-000 |
| 4 | Extrusion, upper dovetail | 1009-3113-000 |
| 5 | Screw, M6x20 | 0144-2131-925 |
| | Cable ties | 0203-5915-300 |
| Refer | to Section 4.21 for mounting the display arms and proper cable dressing. | |


8-66 07/04 1009-0356-000

8.33 Cable management arm



| Item | Description | Stock Number |
|------|---|---------------|
| | Cable management arm, complete assembly | 1009-8181-000 |
| 1 | Bracket, cable management | 1009-3261-000 |
| 2 | Screw, M4x12 | 0140-6226-111 |
| 3 | Lockwasher, M4 external | 9213-0540-003 |
| 4 | Extrusion, cable arm, front loading | 1009-3247-000 |
| 5 | Screw, M4x6 Nyloc | 1009-3283-000 |
| 6 | Retainer, multiple cable | 1009-3252-000 |
| 7 | Retainer, cable small | 1009-3259-000 |
| 8 | Washer M4 flat | 0144-1025-165 |

8.34 Display arm mounting kits for optional equipment



| | 01-1-2-1-00 - 10-10-10-10-10-10-10-10-10-10-10-10-10-1 |
|---|--|
| 1 | Cardiocap 5 mou |
| 2 | C /E Elat Danal m |

3 Spacelabs Flat Panel mount

4 Spacelabs PC Scout mount

Stock Number

1009-3265-000 1009-3266-000 1009-3267-000 1009-3268-000

07/04 1009-0356-000

In this section Schematics are subject to change without notice. Circuit boards are available only as complete assemblies.

| Figure 9-1 | System connection block diagram | 9-2 |
|-------------|--|------|
| Figure 9-2 | Gas scavenging circuits | 9-3 |
| Figure 9-3 | Electrical cabling block diagram | 9-4 |
| Figure 9-4 | Pneumatic circuit diagram | 9-6 |
| Figure 9-5 | Wiring hamesses. | 9-6 |
| Figure 9-6 | Wiring harnesses (Aespire 100). | 9-7 |
| Figure 9-7 | Tubing | 9-8 |
| Figure 9-8 | Schematic, AC Inlet module; 100–120 V (Non-isolated outlets or no outlets) | 9-9 |
| Figure 9-9 | Schematic, AC Inlet module; 100–120 V (Isolated outlets | 9-10 |
| Figure 9-10 | Schematic, AC Inlet module; 220–240 V (Non-isolated outlets or no outlets) | 9-11 |
| Figure 9-11 | Schematic, AC Inlet module; 220–240 V (Isolated outlets) | 9-12 |
| Figure 9-12 | Schematic, AC Power (Aespire 100) | 9-13 |
| | | |






Figure 9-1 • System connection block diagram

9-2 07/04 1009-0356-000





- Key to Symbols

 × = Plugged port (1/8 inch) for sample gas return
- Plugged port (30 mm) for audiliary breathing system scavenging
 Open port (30 mm) for audiliary breathing system scavenging

Note: Active AGSS systems with a 12.7 mm connector do not include the Flow Restrictor and the Flow Indicator.

Figure 9-2 • Gas scavenging circuits




Figure 9-3 • Electrical cabling block diagram 9-4

07/04 1009-0356-000

9-5

Key to Numbered Components 1. Pipeline pressure gauge

- 2. Pipeline inlet
- 3. Cylinder pressure gauge
- 4. Cylinder inlet
- Primary regulator (cylinder pressure)
 High-pressure relief valve (758 kPa / 110 psi)*
- 7. Supply connections for the ventilator
 - a. O₂ drive gas
 - b. Air drive gas
- 8. System switch
- 9. Switch for low 0_2 supply pressure alarm (used with the ventilator)
- 10. O₂ secondary regulator (207 kPa / 30 psi)*
- 11. O₂ flow control valve
- 12. O₂ flow tube(s)
- 13. O₂ flush and auxiliary flowmeter regulator (241 kPa / 35 psi)*
- 14. O₂ Flush
 - a. Flush valve
 - b. Pressure switch (used with the ventilator)
- 15. N₂O balance regulator
- 16. N₂O flow control valve
- 17. N₂O flow tube(s)
- 18. Air secondary regulator (207 kPa / 30 psl)*
- 19. Air flow control valve
- 20. Air flow tube(s)
- 21. Supply connection for Venturi suction
 - a. O₂ drive gas b. Air drive gas
- 22. Vaporizer port valve
- 23. Vaportzer 24. Low-pressure relief valve (38 kPa / 5.5 psi)*
- 25. Auxiliary flowmeter (optional)
- 26. To ABS
- 27. To ACGO
- 28. Test port (primary regulator)
 29. Test port (secondary/balance regulator)
- * Approximate values

Key to Symbols


₩ Pneumatic Connection


0 Filter

D Direction of Flow

0 Check Valve

Figure 9-4 • Preumatic circuit diagram






Key to Symbols

VEB = Vent Engine Board VMB = Ventilator Monitoring Board ACGO = Auxiliary Common Gas Outlet

> Figure 9-5 • Wiring harnesses 9-6

07/04 1009-0356-000




Key to Symbols

VEB = Vent Engine Board
VMB = Ventilator Monitoring Board
ACGO = Auxiliary Common Gas Outlet

Figure 9-6 • Wiring harnesses (Aespire 100) 1009-0395-000 07/04

9-7



The Sample Gas Return is directed to the scavenging system as a factory default. A qualified service representative can reroute the sample gas back to the breathing system. Refer to Section 4.13.

Figure 9-7 • Tubing 9-8

07/04 1009-0356-000




Figure 9-8 • Schematic, AC Inlet module; 100–120 V (Non-isolated outlets or no outlets) 1009-0356-000 07/04

9-9




Figure 9-9 • Schematic, AC Inlet module; 100-120 V (Isolated outlets) 9-10

07/04 1009-0356-000




Figure 9-10 • Schematic, AC Inlet module; 220–240 V (Non-isolated outlets or no outlets) 1009-0356-000 07/04

9-11




Figure 9-11 • Schematic, AC Inlet module; 220-240 V (Isolated outlets)
9-12
07/04
1009-0386-000




Figure 9-12 • Schematic, AC Power (Aespire 100)

9-13

Notes

9-14 07/04 1009-0356-000

S/5 Asspire Machine Technical Reference Manual, English 1009 0386 000 07 04 9 01 01 02 Prioted in USA ©Dates-Ohmeda, Inc. All rights reserved