

## KODAK DryView Laser Imager 8100

The DryView 8100 is an excellent seller worldwide thanks to its remarkable speed and outstanding image quality.

### Includes

- DryView 8100 (14"x17") PACS Link laser imager
- 2 (two) connections: DICOM or Non-DICOM

### Features

- Prints up to 55 sheets of film per hour
- 78-micron laser spot spacing, 325-dpi resolution
- Continuous-tone images at 4,096 gray levels
- Consistent resolution from application to application
- 1:1, 2:1, 4:1, 6:1, 9:1, 12:1, 15:1, 16:1, 20:1 formats
- Accommodates 20:1 format without compression
- Requires no operator intervention, enhancing productivity

Includes installation, applications and 1 year warranty through KODAK.



### Product range

**KODAK DryView Laser Imager 8100  
(2 connections)**

Order No. **K1963024**

**KODAK DryView Laser Imager 8100  
(3 to 5 connections)**

Order No. **K8270811**